


Liberty Tipping Lever



Item Number	Part Number	Description	Quantity
1 - 7	120818	Tipping Lever Assembly Liberty	1
3, 5 - 7	120819	Tipping Lever Mounting Hardware Liberty	1

Item Number	Part Number	Description	Quantity
1	119596	Mount Bracket, P, Tipping Lever	2
4	100804	Tube End Plug 1"	2
5	103385	M8 x 35 HHCS Blk Zn	2
6	112533	HHCS, M8x1.25x50, CLS10.9, BLZ, ISO7380	1
7	101847	M6 x 40 SHCS Blk Zn	2
8	100659	M8 Nylock Nut Blk Zn	2
9	100658	M6 Nylock Nut Blk Zn	2