


Rigid Rear Anti-Tip



If retrofitting, Anti-Tip Receivers will be needed. If retrofitting Anti-Tips, see CS if new camber tube is needed. Used if Alber Emotion/Twion system is present on a Rogue XP

Item Number	Part Number	Description	Quantity
1a, 2, 3-7	119669	Anti Tip Medium 22/24CAT 18/20 Rogue-Clik-LW 22/24/25ETHOS	1
1b, 2, 3 - 7	119670	Anti Tip Long 25/26 CAT 22 ROGUE 22-26 RG ALX/LW XP 26 ETHOS <i>Commonly used with 22 / 24" wheels on Clik, 22" wheels on Rogue, and 22 - 26" wheels on Tsunami and Little Wave</i>	1
1c, 2, 3 - 7	119672	Anti Tip Extra Long 24-26 Rigid Catalyst Spark w/Reclining Back <i>Commonly used with 24 - 26" wheels on Rogue (Replaces X-Large 200524)</i>	1
1d, 2, 3 - 7	119671	Anti Tip Short 18/20 Cat, 18' Clik <i>Commonly used with 18" wheels on Clik</i>	1
1a	000077	Anti Tip Tube Upper Medium <i>Commonly used with 20" wheels on Clik and 18 / 20" wheels on Rogue, Tsunami, and Little Wave</i>	1
1b	000078	Anti Tip Tube Upper Long <i>Commonly used with 22 / 24" wheels on Clik, 22" wheels on Rogue, and 22 - 26" wheels on Tsunami and Little Wave</i>	1
1c	000993	Anti Tip Tube Upper Extra Long <i>Commonly used with 24 - 26" wheels on Rogue</i>	1
1d	102836	Anti-Tip Tube Upper Extra Short <i>Commonly used with 18" wheels on Clik</i>	1
3	100533	Detent Button	1
4	100534	Detent Button Lower	1