


(Discontinued) Spark / Catalyst Aluminum Seat Pan



| Part Number | Tape Size | Seat Pan Depth |
|-----------------|-------------------|----------------|
| 004458 | 3" X 4" | 10" - 11" |
| 004458 & 004459 | 3" X 4" & 3" X 5" | 12" |
| 004459 | 3" X 5" | 13" |
| 004460 | 3" X 6" | 14" - 22" |

This seat pan is used if transit is not present on the chair. The seat pan width and depth measure approximately 0.5" less than the chair size. For a complete seat pan assembly, select the seat pan size needed and hardware kit.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1a, 6a | 506009 | Seat Pan Catalyst / Spark with AB Tape 10 x 10 | 1 |
| 1b, 6a | 506010 | Seat Pan Catalyst / Spark with AB Tape 10 x 11 | 1 |
| 1c, 6a, 6b | 506011 | Seat Pan Catalyst / Spark with AB Tape 10 x 12 | 1 |
| 1d, 6b | 506012 | Seat Pan Catalyst / Spark with AB Tape 10 x 13 | 1 |
| 1e, 6c | 506013 | Seat Pan Catalyst / Spark with AB Tape 10 x 14 | 1 |
| 1f, 6c | 506014 | Seat Pan Catalyst / Spark with AB Tape 10 x 15 | 1 |
| 1g, 6c | 506015 | Seat Pan Catalyst / Spark with AB Tape 10 x 16 | 1 |
| 1h, 6c | 506016 | Seat Pan Catalyst / Spark with AB Tape 10 x 17 | 1 |
| 1i, 6c | 506017 | Seat Pan Catalyst / Spark with AB Tape 10 x 18 | 1 |
| 1j, 6a | 506018 | Seat Pan Catalyst / Spark with AB Tape 11 x 10 | 1 |
| 1k, 6a | 506019 | Seat Pan Catalyst / Spark with AB Tape 11 x 11 | 1 |
| 1l, 6a, 6b | 506020 | Seat Pan Catalyst / Spark with AB Tape 11 x 12 | 1 |
| 1m, 6b | 506021 | Seat Pan Catalyst / Spark with AB Tape 11 x 13 | 1 |
| 1n, 6c | 506022 | Seat Pan Catalyst / Spark with AB Tape 11 x 14 | 1 |
| 1o, 6c | 506023 | Seat Pan Catalyst / Spark with AB Tape 11 x 15 | 1 |
| 1p, 6c | 506024 | Seat Pan Catalyst / Spark with AB Tape 11 x 16 | 1 |
| 1q, 6c | 506025 | Seat Pan Catalyst / Spark with AB Tape 11 x 17 | 1 |
| 1r, 6c | 506026 | Seat Pan Catalyst / Spark with AB Tape 11 x 18 | 1 |
| 1s, 6a | 506027 | Seat Pan Catalyst / Spark with AB Tape 12 x 10 | 1 |
| 1t, 6a | 506028 | Seat Pan Catalyst / Spark with AB Tape 12 x 11 | 1 |
| 1u, 6a, 6b | 506029 | Seat Pan Catalyst / Spark with AB Tape 12 x 12 | 1 |
| 1v, 6b | 506030 | Seat Pan Catalyst / Spark with AB Tape 12 x 13 | 1 |
| 1w, 6c | 506031 | Seat Pan Catalyst / Spark with AB Tape 12 x 14 | 1 |
| 1x, 6c | 506032 | Seat Pan Catalyst / Spark with AB Tape 12 x 15 | 1 |
| 1y, 6c | 506033 | Seat Pan Catalyst / Spark with AB Tape 12 x 16 | 1 |
| 1z, 6c | 506034 | Seat Pan Catalyst / Spark with AB Tape 12 x 17 | 1 |
| 1aa, 6c | 506035 | Seat Pan Catalyst / Spark with AB Tape 12 x 18 | 1 |
| 1ab, 6a | 506036 | Seat Pan Catalyst / Spark with AB Tape 13 x 10 | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1ac, 6a | 506037 | Seat Pan Catalyst / Spark with AB Tape 13 x 11 | 1 |
| 1ad, 6a, 6b | 506038 | Seat Pan Catalyst / Spark with AB Tape 13 x 12 | 1 |
| 1ae, 6c | 506039 | Seat Pan Catalyst / Spark with AB Tape 13 x 13 | 1 |
| 1af, 6c | 506040 | Seat Pan Catalyst / Spark with AB Tape 13 x 14 | 1 |
| 1ag, 6c | 506041 | Seat Pan Catalyst / Spark with AB Tape 13 x 15 | 1 |
| 1ah, 6c | 506042 | Seat Pan Catalyst / Spark with AB Tape 13 x 16 | 1 |
| 1ai, 6c | 506043 | Seat Pan Catalyst / Spark with AB Tape 13 x 17 | 1 |
| 1aj, 6c | 506044 | Seat Pan Catalyst / Spark with AB Tape 13 x 18 | 1 |
| 1ak, 6a | 506045 | Seat Pan Catalyst / Spark with AB Tape 14 x 10 | 1 |
| 1al, 6a | 506046 | Seat Pan Catalyst / Spark with AB Tape 14 x 11 | 1 |
| 1am, 6a, 6b | 506047 | Seat Pan Catalyst / Spark with AB Tape 14 x 12 | 1 |
| 1an, 6b | 506048 | Seat Pan Catalyst / Spark with AB Tape 14 x 13 | 1 |
| 1ao, 6c | 505780 | Seat Pan Catalyst / Spark with AB Tape 14 x 14 | 1 |
| 1ap, 6c | 505781 | Seat Pan Catalyst / Spark with AB Tape 14 x 15 | 1 |
| 1aq, 6c | 505782 | Seat Pan Catalyst / Spark with AB Tape 14 x 16 | 1 |
| 1ar, 6c | 505783 | Seat Pan Catalyst / Spark with AB Tape 14 x 17 | 1 |
| 1as, 6c | 505784 | Seat Pan Catalyst / Spark with AB Tape 14 x 18 | 1 |
| 1at, 6c | 505785 | Seat Pan Catalyst / Spark with AB Tape 14 x 19 | 1 |
| 1au, 6c | 505786 | Seat Pan Catalyst / Spark with AB Tape 14 x 20 | 1 |
| 1av, 6a | 506049 | Seat Pan Catalyst / Spark with AB Tape 15 x 10 | 1 |
| 1aw, 6a | 506050 | Seat Pan Catalyst / Spark with AB Tape 15 x 11 | 1 |
| 1ax, 6a, 6b | 506051 | Seat Pan Catalyst / Spark with AB Tape 15 x 12 | 1 |
| 1ay, 6b | 506052 | Seat Pan Catalyst / Spark with AB Tape 15 x 13 | 1 |
| 1az, 6c | 505787 | Seat Pan Catalyst / Spark with AB Tape 15 x 14 | 1 |
| 1ba, 6c | 505788 | Seat Pan Catalyst / Spark with AB Tape 15 x 15 | 1 |
| 1bb, 6c | 505789 | Seat Pan Catalyst / Spark with AB Tape 15 x 16 | 1 |
| 1bc, 6c | 505790 | Seat Pan Catalyst / Spark with AB Tape 15 x 17 | 1 |
| 1bd, 6c | 505791 | Seat Pan Catalyst / Spark with AB Tape 15 x 18 | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1be, 6c | 505792 | Seat Pan Catalyst / Spark with AB Tape 15 x 19 | 1 |
| 1bf, 6c | 505793 | Seat Pan Catalyst / Spark with AB Tape 15 x 20 | 1 |
| 1bg, 6a | 506053 | Seat Pan Catalyst / Spark with AB Tape 16 x 10 | 1 |
| 1bh, 6a | 506054 | Seat Pan Catalyst / Spark with AB Tape 16 x 11 | 1 |
| 1bi, 6a, 6b | 506055 | Seat Pan Catalyst / Spark with AB Tape 16 x 12 | 1 |
| 1bj, 6b | 506056 | Seat Pan Catalyst / Spark with AB Tape 16 x 13 | 1 |
| 1bk, 6c | 505794 | Seat Pan Catalyst / Spark with AB Tape 16 x 14 | 1 |
| 1bl, 6c | 505795 | Seat Pan Catalyst / Spark with AB Tape 16 x 15 | 1 |
| 1bm, 6c | 505796 | Seat Pan Catalyst / Spark with AB Tape 16 x 16 | 1 |
| 1bn, 6c | 505797 | Seat Pan Catalyst / Spark with AB Tape 16 x 17 | 1 |
| 1bo, 6c | 505798 | Seat Pan Catalyst / Spark with AB Tape 16 x 18 | 1 |
| 1bp, 6c | 505799 | Seat Pan Catalyst / Spark with AB Tape 16 x 19 | 1 |
| 1bq, 6c | 505800 | Seat Pan Catalyst / Spark with AB Tape 16 x 20 | 1 |
| 1br, 6c | 505801 | Seat Pan Catalyst / Spark with AB Tape 16 x 21 | 1 |
| 1bs, 6c | 505802 | Seat Pan Catalyst / Spark with AB Tape 16 x 22 | 1 |
| 1bt, 6c | 505803 | Seat Pan Catalyst / Spark with AB Tape 17 x 14 | 1 |
| 1bu, 6c | 505804 | Seat Pan Catalyst / Spark with AB Tape 17 x 15 | 1 |
| 1bv, 6c | 505805 | Seat Pan Catalyst / Spark with AB Tape 17 x 16 | 1 |
| 1bw, 6c | 505806 | Seat Pan Catalyst / Spark with AB Tape 17 x 17 | 1 |
| 1bx, 6c | 505807 | Seat Pan Catalyst / Spark with AB Tape 17 x 18 | 1 |
| 1by, 6c | 505808 | Seat Pan Catalyst / Spark with AB Tape 17 x 19 | 1 |
| 1bz, 6c | 505809 | Seat Pan Catalyst / Spark with AB Tape 17 x 20 | 1 |
| 1ca, 6c | 505810 | Seat Pan Catalyst / Spark with AB Tape 17 x 21 | 1 |
| 1cb, 6c | 505811 | Seat Pan Catalyst / Spark with AB Tape 17 x 22 | 1 |
| 1cc, 6c | 505812 | Seat Pan Catalyst / Spark with AB Tape 18 x 14 | 1 |
| 1cd, 6c | 505813 | Seat Pan Catalyst / Spark with AB Tape 18 x 15 | 1 |
| 1ce, 6c | 505814 | Seat Pan Catalyst / Spark with AB Tape 18 x 16 | 1 |
| 1cf, 6c | 505815 | Seat Pan Catalyst / Spark with AB Tape 18 x 17 | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1cg, 6c | 505816 | Seat Pan Catalyst / Spark with AB Tape 18 x 18 | 1 |
| 1ch, 6c | 505817 | Seat Pan Catalyst / Spark with AB Tape 18 x 19 | 1 |
| 1ci, 6c | 505818 | Seat Pan Catalyst / Spark with AB Tape 18 x 20 | 1 |
| 1cj, 6c | 505819 | Seat Pan Catalyst / Spark with AB Tape 18 x 21 | 1 |
| 1ck, 6c | 505820 | Seat Pan Catalyst / Spark with AB Tape 18 x 22 | 1 |
| 1cl, 6c | 505821 | Seat Pan Catalyst / Spark with AB Tape 19 x 14 | 1 |
| 1cm, 6c | 505822 | Seat Pan Catalyst / Spark with AB Tape 19 x 15 | 1 |
| 1cn, 6c | 505823 | Seat Pan Catalyst / Spark with AB Tape 19 x 16 | 1 |
| 1co, 6c | 505824 | Seat Pan Catalyst / Spark with AB Tape 19 x 17 | 1 |
| 1cp, 6c | 505825 | Seat Pan Catalyst / Spark with AB Tape 19 x 18 | 1 |
| 1cq, 6c | 505826 | Seat Pan Catalyst / Spark with AB Tape 19 x 19 | 1 |
| 1cr, 6c | 505827 | Seat Pan Catalyst / Spark with AB Tape 19 x 20 | 1 |
| 1cs, 6c | 505828 | Seat Pan Catalyst / Spark with AB Tape 19 x 21 | 1 |
| 1ct, 6c | 505829 | Seat Pan Catalyst / Spark with AB Tape 19 x 22 | 1 |
| 1cu, 6c | 505830 | Seat Pan Catalyst / Spark with AB Tape 20 x 14 | 1 |
| 1cv, 6c | 505831 | Seat Pan Catalyst / Spark with AB Tape 20 x 15 | 1 |
| 1cw, 6c | 505832 | Seat Pan Catalyst / Spark with AB Tape 20 x 16 | 1 |
| 1cx, 6c | 505833 | Seat Pan Catalyst / Spark with AB Tape 20 x 17 | 1 |
| 1cy, 6c | 505834 | Seat Pan Catalyst / Spark with AB Tape 20 x 18 | 1 |
| 1cz, 6c | 505835 | Seat Pan Catalyst / Spark with AB Tape 20 x 19 | 1 |
| 1da, 6c | 505836 | Seat Pan Catalyst / Spark with AB Tape 20 x 20 | 1 |
| 1db, 6c | 505837 | Seat Pan Catalyst / Spark with AB Tape 20 x 21 | 1 |
| 1dc, 6c | 505838 | Seat Pan Catalyst / Spark with AB Tape 20 x 22 | 1 |
| 1dd, 6c | 505839 | Seat Pan Catalyst / Spark with AB Tape 21 x 14 | 1 |
| 1de, 6c | 505840 | Seat Pan Catalyst / Spark with AB Tape 21 x 15 | 1 |
| 1df, 6c | 505841 | Seat Pan Catalyst / Spark with AB Tape 21 x 16 | 1 |
| 1dg, 6c | 505842 | Seat Pan Catalyst / Spark with AB Tape 21 x 17 | 1 |
| 1dh, 6c | 505843 | Seat Pan Catalyst / Spark with AB Tape 21 x 18 | 1 |

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|--|----------|
| 1di, 6c | 505844 | Seat Pan Catalyst / Spark with AB Tape 21 x 19 | 1 |
| 1dj, 6c | 505845 | Seat Pan Catalyst / Spark with AB Tape 21 x 20 | 1 |
| 1dk, 6c | 505846 | Seat Pan Catalyst / Spark with AB Tape 21 x 21 | 1 |
| 1dl, 6c | 505847 | Seat Pan Catalyst / Spark with AB Tape 21 x 22 | 1 |
| 1dm, 6c | 505848 | Seat Pan Catalyst / Spark with AB Tape 22 x 14 | 1 |
| 1dn, 6c | 505849 | Seat Pan Catalyst / Spark with AB Tape 22 x 15 | 1 |
| 1do, 6c | 505850 | Seat Pan Catalyst / Spark with AB Tape 22 x 16 | 1 |
| 1dp, 6c | 505851 | Seat Pan Catalyst / Spark with AB Tape 22 x 17 | 1 |
| 1dq, 6c | 505852 | Seat Pan Catalyst / Spark with AB Tape 22 x 18 | 1 |
| 1dr, 6c | 505853 | Seat Pan Catalyst / Spark with AB Tape 22 x 19 | 1 |
| 1ds, 6c | 505854 | Seat Pan Catalyst / Spark with AB Tape 22 x 20 | 1 |
| 1dt, 6c | 505855 | Seat Pan Catalyst / Spark with AB Tape 22 x 21 | 1 |
| 1du, 6c | 505856 | Seat Pan Catalyst / Spark with AB Tape 22 x 22 | 1 |
| 2 - 5 | 115613 | Spark / Catalyst Aluminum Seat Pan Hardware | 1 |
| 2 | 000990 | Seat Pan Clip CAT | 4 |
| 3 | 100667 | M6 x 12 BTN SHCS BLKZC w/Patch | 4 |
| 4 | 100658 | M6 Nylock Nut Blk Zn | 4 |
| 5 | 102856 | Seat Pan Stop Collar | 2 |
| 6a | 004458 | AB Tape, 3" x 4" | 4 |
| 6b | 004459 | AB Tape, 3" x 5" | 4 |
| 6c | 004460 | AB Tape, 3" x 6" | 4 |