


Rogue ALX Aluminum Seat Pan - Growth (Formerly Tsunami)



Part Number	Tape Size	Seat Pan Depth
004460	3" X 6"	14" - 22"

Screw (Item 2) Qty Per Chair	
Seat Depth	Screw Qty
14" - 15"	8
16" - 18"	10
19" - 21"	12
22"	14

Rigid Seat Pans may be quoted, but must be ordered through Customer Service. If your chair was ordered prior to February 2013 a custom seat pan may be required. The seat pan width measures approximately 0.75" and the seat pan depth measures approximately 0.50" less (however a depth of 21" will be closer to 0.75") than the chair size.

Item Number	Part Number	Description	Quantity
1a, 3	505859	Seat Pan Rigid with AB Tape 14 x 14	1
1b, 3	505860	Seat Pan Rigid with AB Tape 14 x 15	1
1c, 3	505861	Seat Pan Rigid with AB Tape 14 x 16	1
1d, 3	505862	Seat Pan Rigid with AB Tape 14 x 17	1

Item Number	Part Number	Description	Quantity
1e, 3	505863	Seat Pan Rigid with AB Tape 14 x 18	1
1f, 3	505864	Seat Pan Rigid with AB Tape 14 x 19	1
1g, 3	505865	Seat Pan Rigid with AB Tape 14 x 20	1
1h, 3	505866	Seat Pan Rigid with AB Tape 14 x 21	1
1i, 3	505867	Seat Pan Rigid with AB Tape 14 x 22	1
1j, 3	505868	Seat Pan Rigid with AB Tape 15 x 14	1
1k, 3	505869	Seat Pan Rigid with AB Tape 15 x 15	1
1l, 3	505870	Seat Pan Rigid with AB Tape 15 x 16	1
1m, 3	505871	Seat Pan Rigid with AB Tape 15 x 17	1
1n, 3	505872	Seat Pan Rigid with AB Tape 15 x 18	1
1o, 3	505873	Seat Pan Rigid with AB Tape 15 x 19	1
1p, 3	505874	Seat Pan Rigid with AB Tape 15 x 20	1
1q, 3	505875	Seat Pan Rigid with AB Tape 15 x 21	1
1r, 3	505876	Seat Pan Rigid with AB Tape 15 x 22	1
1s, 3	505877	Seat Pan Rigid with AB Tape 16 x 14	1
1t, 3	505878	Seat Pan Rigid with AB Tape 16 x 15	1
1u, 3	505879	Seat Pan Rigid with AB Tape 16 x 16	1
1v, 3	505880	Seat Pan Rigid with AB Tape 16 x 17	1
1w, 3	505881	Seat Pan Rigid with AB Tape 16 x 18	1
1x, 3	505882	Seat Pan Rigid with AB Tape 16 x 19	1
1y, 3	505883	Seat Pan Rigid with AB Tape 16 x 20	1
1z, 3	505884	Seat Pan Rigid with AB Tape 16 x 21	1
1aa, 3	505885	Seat Pan Rigid with AB Tape 16 x 22	1
1ab, 3	505886	Seat Pan Rigid with AB Tape 17 x 14	1
1ac, 3	505887	Seat Pan Rigid with AB Tape 17 x 15	1
1ad, 3	505888	Seat Pan Rigid with AB Tape 17 x 16	1
1ae, 3	505889	Seat Pan Rigid with AB Tape 17 x 17	1
1af, 3	505890	Seat Pan Rigid with AB Tape 17 x 18	1

Item Number	Part Number	Description	Quantity
1ag, 3	505891	Seat Pan Rigid with AB Tape 17 x 19	1
1ah, 3	505892	Seat Pan Rigid with AB Tape 17 x 20	1
1ai, 3	505893	Seat Pan Rigid with AB Tape 17 x 21	1
1aj, 3	505894	Seat Pan Rigid with AB Tape 17 x 22	1
1ak, 3	505895	Seat Pan Rigid with AB Tape 18 x 14	1
1al, 3	505896	Seat Pan Rigid with AB Tape 18 x 15	1
1am, 3	505897	Seat Pan Rigid with AB Tape 18 x 16	1
1an, 3	505898	Seat Pan Rigid with AB Tape 18 x 17	1
1ao, 3	505899	Seat Pan Rigid with AB Tape 18 x 18	1
1ap, 3	505900	Seat Pan Rigid with AB Tape 18 x 19	1
1aq, 3	505901	Seat Pan Rigid with AB Tape 18 x 20	1
1ar, 3	505902	Seat Pan Rigid with AB Tape 18 x 21	1
1as, 3	505903	Seat Pan Rigid with AB Tape 18 x 22	1
1at, 3	505904	Seat Pan Rigid with AB Tape 19 x 14	1
1au, 3	505905	Seat Pan Rigid with AB Tape 19 x 15	1
1av, 3	505906	Seat Pan Rigid with AB Tape 19 x 16	1
1aw, 3	505907	Seat Pan Rigid with AB Tape 19 x 17	1
1ax, 3	505908	Seat Pan Rigid with AB Tape 19 x 18	1
1ay, 3	505909	Seat Pan Rigid with AB Tape 19 x 19	1
1az, 3	505910	Seat Pan Rigid with AB Tape 19 x 20	1
1ba, 3	505911	Seat Pan Rigid with AB Tape 19 x 21	1
1bb, 3	505912	Seat Pan Rigid with AB Tape 19 x 22	1
1bc, 3	505913	Seat Pan Rigid with AB Tape 20 x 14	1
1bd, 3	505914	Seat Pan Rigid with AB Tape 20 x 15	1
1be, 3	505915	Seat Pan Rigid with AB Tape 20 x 16	1
1bf, 3	505916	Seat Pan Rigid with AB Tape 20 x 17	1
1bg, 3	505917	Seat Pan Rigid with AB Tape 20 x 18	1
1bh, 3	505918	Seat Pan Rigid with AB Tape 20 x 19	1

Item Number	Part Number	Description	Quantity
1bi, 3	505919	Seat Pan Rigid with AB Tape 20 x 20	1
1bj, 3	505920	Seat Pan Rigid with AB Tape 20 x 21	1
1bk, 3	505921	Seat Pan Rigid with AB Tape 20 x 22	1
1bl, 3	505922	Seat Pan Rigid with AB Tape 21 x 14	1
1bm, 3	505923	Seat Pan Rigid with AB Tape 21 x 15	1
1bn, 3	505924	Seat Pan Rigid with AB Tape 21 x 16	1
1bo, 3	505925	Seat Pan Rigid with AB Tape 21 x 17	1
1bp, 3	505926	Seat Pan Rigid with AB Tape 21 x 18	1
1bq, 3	505927	Seat Pan Rigid with AB Tape 21 x 19	1
1br, 3	505928	Seat Pan Rigid with AB Tape 21 x 20	1
1bs, 3	505929	Seat Pan Rigid with AB Tape 21 x 21	1
1bt, 3	505930	Seat Pan Rigid with AB Tape 21 x 22	1
1bu, 3	505931	Seat Pan Rigid with AB Tape 22 x 14	1
1bv, 3	505932	Seat Pan Rigid with AB Tape 22 x 15	1
1bw, 3	505933	Seat Pan Rigid with AB Tape 22 x 16	1
1bx, 3	505934	Seat Pan Rigid with AB Tape 22 x 17	1
1by, 3	505935	Seat Pan Rigid with AB Tape 22 x 18	1
1bz, 3	505936	Seat Pan Rigid with AB Tape 22 x 19	1
1ca, 3	505937	Seat Pan Rigid with AB Tape 22 x 20	1
1cb, 3	505938	Seat Pan Rigid with AB Tape 22 x 21	1
1cc, 3	505939	Seat Pan Rigid with AB Tape 22 x 22	1
2	100734	M5 x 16 PHMS ZC	14
3	004460	AB Tape, 3" x 6"	4