


## **CR45 24" MAXX Mag Wheel / Tire / Handrim Kits**

---



*Not all components of the Wheel, Tire and handrim kits will be listed. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size, new axles/receivers and handrims may be needed. If changing wheel size and Attendant Foot Lock or T-Arms are present, spacer 112503 is required which can be found on the CR45 Axle Plate page See CR45 Axles page for axle and receiver utilization.*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
	118359	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118360	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118361	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118362	24" MAXX Mag, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118367	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118368	24" MAXX Mag, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118369	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	118370	24" MAXX Mag, Pneumatic Tire (Black), Projection Handrim Kit, Close Mount	1
	118136	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118137	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount	1
	118138	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118139	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	118144	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118145	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Close Mount	1
	118146	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	118147	24" MAXX Mag, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Close Mount	1
	117696	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	117697	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount	1
	117698	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	117699	24" MAXX Mag, Full Poly Tire, Aluminum Anodized Handrim Kit, Close Mount w/ Micro Grip	1
	117704	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	117705	24" MAXX Mag, Full Poly Tire, Plastic Coated Handrim Kit, Close Mount	1
	117706	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	117707	24" MAXX Mag, Full Poly Tire, Projection Handrim Kit, Close Mount	1