


# Clik Wheels - Spinerger LX



Recommended to order the Spinerger LX Replacement Spoke Tool Kit (110119) when replacing spokes. Replacement Spinerger Spokes and Spinerger Wheels with Flex Rims will have an extended lead time. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

| Item Number | Part Number | Description                         | Quantity |
|-------------|-------------|-------------------------------------|----------|
| 1a          | 100886      | 22" x 1/2" Spinerger LX Black Spoke | 1        |
| 1b          | 100939      | 22" x 1/2" Spinerger LX Red Spoke   | 1        |
| 1c          | 100940      | 22" x 1/2" Spinerger LX White Spoke | 1        |

| Item Number | Part Number | Description                                                                                                                           | Quantity |
|-------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------|----------|
| 1d          | 100938      | 22" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 1j          | 100942      | 24" x 1/2" Spinergy LX Black Spoke                                                                                                    | 1        |
| 1k          | 100941      | 24" x 1/2" Spinergy LX Red Spoke                                                                                                      | 1        |
| 1l          | 100944      | 24" x 1/2" Spinergy LX White Spoke                                                                                                    | 1        |
| 1m          | 100943      | 24" x 1/2" Spinergy LX Blue Spoke                                                                                                     | 1        |
| 1s          | 118348      | 24" x 1/2" Spinergy LX Black Spoke w/<br>Pneumatic w/ Airless Insert Tire (Black)                                                     | 1        |
| 1t          | 118349      | 24" x 1/2" Spinergy LX Red Spoke w/ Pneumatic<br>w/ Airless Insert Tire (Black)                                                       | 1        |
| 1u          | 118350      | 24" x 1/2" Spinergy LX White Spoke w/<br>Pneumatic w/ Airless Insert Tire (Black)                                                     | 1        |
| 1v          | 118351      | 24" x 1/2" Spinergy LX Blue Spoke w/ Pneumatic<br>w/ Airless Insert Tire (Black)                                                      | 1        |
| 2           | 102854      | Spinergy LX Spoke Wrench<br><br><i>No charge when ordered with Spinergy wheel or<br/>replacement spoke</i>                            | 1        |
|             | 110119      | Spinergy LX Replacement Spoke Tool Kit<br><br><i>This tool kit contains: Spinergy LX Spoke Wrench and<br/>Spinergy Spoke T-Wrench</i> | 1        |
| 3a          | 101929      | Replacement Spoke Spinergy LX 22", Black<br>(511291-8.75BLA)                                                                          | 1        |
| 3b          | 101930      | Replacement Spoke Spinergy LX 22", Red<br>(511291-8.75RED)                                                                            | 1        |
| 3c          | 101931      | Replacement Spoke Spinergy LX 22", White<br>(511291-8.75WHT)                                                                          | 1        |
| 3d          | 101932      | Replacement Spoke Spinergy LX 22", Blue<br>(511291-8.75BLU)                                                                           | 1        |
| 3i          | 101937      | Replacement Spoke Spinergy LX 24", Black<br>(511291-9.5BLA)                                                                           | 1        |
| 3j          | 101938      | Replacement Spoke Spinergy LX 24", Red<br>(511291-9.5RED)                                                                             | 1        |
| 3k          | 101939      | Replacement Spoke Spinergy LX 24", White<br>(511291-9.5WHT)                                                                           | 1        |
| 3l          | 101940      | Replacement Spoke Spinergy LX 24", Blue<br>(511291-9.5BLU)                                                                            | 1        |

## Tire Component Chart

| Wheel             |                                 | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|-------------------|---------------------------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel             | Wheel Part #                    | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 22" Spoke         | 116730                          | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spoke         | 116728                          | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 18" MAXX Spoke    | 109285                          | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" MAXX Spoke    | 108244                          | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Spoke    | 105135                          | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Spoke    | 105136                          | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" MAXX Spoke    | 107436                          |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" MAXX Spoke    | 107437                          | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" MAXX Mag      | 111853                          | 116843                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" MAXX Mag      | 111854                          | 116735                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" MAXX Mag      | 111855                          | 116736                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 24" Superlight    | 101159                          | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    |                |            |           |
| 25" Superlight    | 101158                          |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           |                |            |           |
| 20" Spinergy Spox | See<br>Spinergy<br>Spox<br>Page | 116843                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Spinergy Spox |                                 | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spinergy Spox |                                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spinergy Spox |                                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 10145     |
| 26" Spinergy Spox |                                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 22" Spingery LX   | See<br>Spinergy<br>LX<br>Page   | 116735                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Spingery LX   |                                 | 116736                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Spingery LX   |                                 |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Spingery LX   |                                 | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |

See Page 2 for Discontinued Wheels.

### Tire Component Chart (Discontinued Wheels)

| Wheel              |              | 1 3/8"<br>Street Pneumatic |            |           | 1"<br>High Pressure |            |           | 1 3/8"<br>Euro Trek Knobby |            |           | 1"<br>Schwalbe |          |            |           | 1 3/8"<br>Schwalbe HS-127 |            |           | 2"<br>Mountain |            |           |
|--------------------|--------------|----------------------------|------------|-----------|---------------------|------------|-----------|----------------------------|------------|-----------|----------------|----------|------------|-----------|---------------------------|------------|-----------|----------------|------------|-----------|
| Wheel              | Wheel Part # | Tire                       | Inner Tube | Rim Strip | Tire                | Inner Tube | Rim Strip | Tire                       | Inner Tube | Rim Strip | Mara-thon Tire | One Tire | Inner Tube | Rim Strip | Tire                      | Inner Tube | Rim Strip | Tire           | Inner Tube | Rim Strip |
| 18" Ki Spoke       | 200529       | 101587                     | 101588     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 20" Ki Spoke       | 200530       | 100911                     | 100912     | 101137    | 102957              | 102844     | 101137    |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" Ki Spoke       | 200531       | 100634                     | 100588     | 101908    | 102961              | 100609     | 101908    |                            |            |           |                |          |            |           |                           |            |           | 101475         | 100588     | 101908    |
| 24" Ki Spoke       | 200532       | 111818                     | 100589     | 100596    | 102958              | 100580     | 100596    | 102982                     | 100589     | 100596    | 100864         | 103176   | 100580     | 100596    | 102920                    | 100589     | 100596    | 101476         | 101188     | 100596    |
| 25" Ki Spoke       | 200533       |                            |            |           | 102959              | 101457     | 101452    |                            |            |           | 101357         | 103177   | 101457     | 101452    |                           |            |           | 100166         | 100208     | 101452    |
| 26" Ki Spoke       | 200534       | 100636                     | 100590     | 100608    | 102960              | 100582     | 100597    |                            |            |           | 100865         |          | 100582     | 100597    |                           |            |           |                |            |           |
| 20" 5-Spoke X Core | 101961       | 100911                     | 100912     | 101137    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 22" 5-Spoke X Core | 101962       | 100634                     | 100588     | 101908    |                     |            |           |                            |            |           |                |          |            |           |                           |            |           |                |            |           |
| 24" 5-Spoke X Core | 100960       | 111818                     | 100589     | 100596    |                     |            |           | 102982                     | 100589     | 100596    |                |          |            |           | 102920                    | 100589     | 100596    |                |            |           |

| Tire Compatibility Chart |                        |                         |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
|--------------------------|------------------------|-------------------------|--------------------------|---------------------|----------------------------------------------------|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                    |                        | Tire                    |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
| Wheel                    | Wheel Part #           | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 22" Spoke                | 116730                 | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spoke                | 116728                 | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 18" MAXX Spoke           | 109285                 | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Spoke           | 108244                 | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" MAXX Spoke           | 105135                 | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" MAXX Spoke           | 105136                 | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" MAXX Spoke           | 107436                 |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" MAXX Spoke           | 107437                 | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 12" Mag <sup>1</sup>     | N/A                    | 001424                  | 001426                   |                     | 001425                                             |                   |               |                  |                         |                      |                 |                        |             |
| 16" Mag <sup>1</sup>     | N/A                    | 001159                  | 001161                   | 003453              | 001160                                             |                   |               |                  |                         |                      |                 |                        |             |
| 20" MAXX Mag             | 111853                 | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            |               |                  |                         |                      |                 |                        |             |
| 22" MAXX Mag             | 111584                 | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            |               |                  |                         |                      |                 |                        | 101475      |
| 24" MAXX Mag             | 111855                 | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            |               |                  | 102982                  |                      |                 | 102920                 | 101476      |
| 24" Superlight           | 101159                 | 100635                  |                          |                     |                                                    |                   | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 |             |
| 25" Superlight           | 101158                 |                         |                          |                     |                                                    |                   | 101146        | 102959           |                         | 101357               | 103177          |                        |             |
| 20" Spinergy Spox        | See Spinergy Spox Page | 116843                  | 100917                   | 102085              | 116843 / 100913                                    | 102974            |               | 102957           |                         |                      |                 |                        |             |
| 22" Spinergy Spox        |                        | 116735                  | 100621                   | 100632              | 116735 / 100503                                    | 103165            |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy Spox        |                        | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy Spox        |                        |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy Spox        |                        | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |
| 22" Spinergy LX          | See Spinergy LX Page   | 116735                  | 100621                   | 100632              | 116735 / 100503                                    |                   |               | 102961           |                         |                      |                 |                        | 101475      |
| 24" Spinergy LX          |                        | 116736                  | 100622                   | 100633              | 116736 / 100504                                    | 102976            | 100894        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Spinergy LX          |                        |                         |                          |                     |                                                    | 102977            | 101146        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Spinergy LX          |                        | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100895        | 102960           |                         | 100865               |                 |                        |             |

<sup>1</sup> 12" and 16" Mag Wheels are sold as a combination Wheel and Tire Assembly

See Page 2 for Discontinued Wheels

| Tire Compatibility Chart (Discontinued Wheels) |              |                         |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
|------------------------------------------------|--------------|-------------------------|--------------------------|---------------------|----------------------------------------------------|-------------------|---------------|------------------|-------------------------|----------------------|-----------------|------------------------|-------------|
| Wheel                                          |              | Tire                    |                          |                     |                                                    |                   |               |                  |                         |                      |                 |                        |             |
| Wheel                                          | Wheel Part # | 1 3/8" Street Pneumatic | 1 3/8" Full Profile Poly | 1" Low Profile Poly | 1 3/8" Pneumatic w/ Airless Insert (Tire / Insert) | 1" Solid Sentinel | 1" Solid Shox | 1" High Pressure | 1 3/8" Euro Trek Knobby | 1" Schwalbe Marathon | 1" Schwalbe One | 1 3/8" Schwalbe HS-127 | 2" Mountain |
| 18" Ki Spoke                                   | 200529       | 101587                  | 101992                   | 100918              | 101587 / 101092                                    |                   |               |                  |                         |                      |                 |                        |             |
| 20" Ki Spoke                                   | 200530       | 100911                  | 100917                   | 102085              | 100911 / 100913                                    | 102974            | 102717        | 102957           |                         |                      |                 |                        |             |
| 22" Ki Spoke                                   | 200531       | 100634                  | 100621                   | 100632              | 100634 / 100503                                    | 103165            | 102718        | 102961           |                         |                      |                 |                        | 101475      |
| 24" Ki Spoke                                   | 200532       | 111818                  | 100622                   | 100633              | 111818 / 100504                                    | 102976            | 100640        | 102958           | 102982                  | 100864               | 103176          | 102920                 | 101476      |
| 25" Ki Spoke                                   | 200533       |                         |                          |                     |                                                    | 102977            | 101453        | 102959           |                         | 101357               | 103177          |                        | 100166      |
| 26" Ki Spoke                                   | 200534       | 100636                  |                          |                     | 100636 / 100505                                    | 102978            | 100884        | 102960           |                         | 100865               |                 |                        |             |
| 20" 5-Spoke X Core                             | 101961       | 100911                  | 100917                   | 102085              | 100911 / 100913                                    |                   | 102717        |                  |                         |                      |                 |                        |             |
| 22" 5-Spoke X Core                             | 101962       | 100634                  | 100621                   | 100632              | 100634 / 100503                                    |                   | 102718        |                  |                         |                      |                 |                        |             |
| 24" 5-Spoke X Core                             | 100960       | 111818                  | 100622                   | 100633              | 111818 / 100504                                    |                   | 100640        |                  | 102982                  |                      |                 | 102920                 |             |

Effective 2/27/2023

| Handrim Hardware Chart |                |                    |                                                              |                                              |                                                                  |                   |                |                  |                  |                    |                  |                      |                  |          |  |
|------------------------|----------------|--------------------|--------------------------------------------------------------|----------------------------------------------|------------------------------------------------------------------|-------------------|----------------|------------------|------------------|--------------------|------------------|----------------------|------------------|----------|--|
| Wheel                  |                | Handrim Connection |                                                              |                                              |                                                                  | Handrim           |                |                  |                  |                    |                  |                      |                  |          |  |
| Wheel                  | Wheel Part No. | Connection Points  | Rivet/Tab                                                    | Spacer*                                      | Screw*                                                           | Aluminum Anodized | Superlight     | Plastic Coated   | Projection       | Ergonomic Standard | Ergonomic LT     | Natural Fit Standard | Natural Fit LT   | Flex Rim |  |
| 22" Spoke              | 116730         | 5 or 6             | 116732 (Used with 5 Rivet)<br><br>100698 (Used with 6 Rivet) | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection<br>Ergonomic: 100654           | 101898 (5 Rivet)  |                | 101964 (5 Rivet) | 112824 (6 Rivet) | 113082 (6 Rivet)   | 113085 (6 Rivet) | 200538 (6 Rivet)     | 200201 (6 Rivet) |          |  |
| 24" Spoke              | 116728         | 5 or 6             | Not used w/ Superlight Handrims                              |                                              | Natural Fit: 100835<br><br>Superlight: 100669<br><br>Nut: 100657 | 100975 (5 Rivet)  | 101161 (6 Tab) | 100976 (5 Rivet) | 112825 (6 Rivet) | 113083 (6 Rivet)   | 113806 (6 Rivet) | 100793 (6 Rivet)     | 200202 (6 Rivet) |          |  |
| 18" MAXX Spoke         | 109285         | 3                  | 100698 (Not used w/ Superlight Handrim)                      | 100653 (Not used w/ Superlight Handrim)      | Aluminum, Plastic Ctd, Projection, Ergonomic: 100654             | 100206            |                |                  |                  |                    |                  |                      |                  |          |  |
| 20" MAXX Spoke         | 108244         | 6                  |                                                              |                                              | Natural Fit: 100835                                              | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Spoke         | 105135         |                    |                                                              |                                              | Superlight: 100669                                               | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Spoke         | 105136         |                    |                                                              |                                              | Nut: 100657                                                      | 200349            | 101161         | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 25" MAXX Spoke         | 107436         |                    |                                                              |                                              | 200350                                                           |                   | 112822         | 200548           | 113084           | 113087             | 200539           | 200540               |                  |          |  |
| 26" MAXX Spoke         | 107437         |                    |                                                              |                                              | 200351                                                           |                   | 112823         | 200549           |                  |                    | 100907           | 101454               |                  |          |  |
| 20" MAXX Mag           | 111853         | 6                  |                                                              | Aluminum, Plastic Coated, Projection: 100629 | Aluminum, Plastic Coated, Projection, Ergonomic: 103545          | 200536            |                | 112819           |                  |                    |                  |                      |                  |          |  |
| 22" MAXX Mag           | 111854         |                    |                                                              | Natural Fit & Ergonomic: 101756              | Natural Fit: 101893                                              | 100560            |                | 112820           | 112824           | 113082             | 113085           | 200538               | 200201           |          |  |
| 24" MAXX Mag           | 111855         |                    |                                                              |                                              |                                                                  | 105137            |                | 112821           | 112825           | 113083             | 113086           | 100793               | 200202           |          |  |
| 24" Superlight         | 101159         | 6                  |                                                              |                                              | 100536                                                           | 100754            | 101161         | 100836           |                  | 113080             | 113081           | 100830               | 100828           |          |  |
| 25" Superlight         | 101160         |                    |                                                              |                                              |                                                                  |                   |                | 101160           | 101091           |                    |                  |                      | 101464           | 101460   |  |

\* Spacer and and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023

| Handrim Hardware Chart (Continued) |                         |                    |           |         |                    |                   |            |                |            |                    |              |                      |                |                       |
|------------------------------------|-------------------------|--------------------|-----------|---------|--------------------|-------------------|------------|----------------|------------|--------------------|--------------|----------------------|----------------|-----------------------|
| Wheel                              |                         | Handrim Connection |           |         |                    | Handrim           |            |                |            |                    |              |                      |                |                       |
| Wheel                              | Wheel Part No.          | Connection Points  | Rivet/Tab | Spacer* | Screw*             | Aluminum Anodized | Superlight | Plastic Coated | Projection | Ergonomic Standard | Ergonomic LT | Natural Fit Standard | Natural Fit LT | Flex Rim              |
| 20" Spinerger Spox                 | See Spinerger Spox Page | 6                  |           |         | Ergonomic: 100536  | 103125            |            | 103179         |            |                    |              |                      |                |                       |
| 22" Spinerger Spox                 |                         |                    |           |         | All Others: 100669 | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger Spox                 |                         |                    |           |         |                    | 100766            | 101161     | 100615         |            |                    |              | 100830               | 100828         |                       |
| 25" Spinerger Spox                 |                         |                    |           |         | Nut: 100657        | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger Spox                 |                         |                    |           |         |                    | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |
| 22" Spinerger LX                   | See Spinerger LX Page   | 6                  |           |         | Ergonomic: 100536  | 100827            |            | 100808         |            |                    |              | 100889               | 100888         |                       |
| 24" Spinerger LX                   |                         |                    |           |         | All Others: 100669 | 100766            | 101161     | 100615         |            | 113080             | 113081       | 100830               | 100828         | See Spinerger LX Page |
| 25" Spinerger LX                   |                         |                    |           |         |                    | 100767            | 101160     | 100765         |            |                    |              | 101464               | 101460         |                       |
| 26" Spinerger LX                   |                         |                    |           |         | Nut: 100657        | 101477            |            | 101148         |            |                    |              | 200200               | 100950         |                       |

| Handrim Hardware Chart (Discontinued Wheels) |                |                    |           |                                           |                                           |                     |            |                |            |                      |                |          |  |
|----------------------------------------------|----------------|--------------------|-----------|-------------------------------------------|-------------------------------------------|---------------------|------------|----------------|------------|----------------------|----------------|----------|--|
| Wheel                                        |                | Handrim Connection |           |                                           |                                           | Handrim             |            |                |            |                      |                |          |  |
| Wheel                                        | Wheel Part No. | Connection Points  | Rivet/Tab | Spacer*                                   | Screw*                                    | Aluminum Anodized   | Superlight | Plastic Coated | Projection | Natural Fit Standard | Natural Fit LT | Flex Rim |  |
| 18" Ki Spoke                                 | 200529         | 3                  | 100698    | 100653                                    | Aluminum, Plastic Ctd, Projection: 100654 | 100206              |            | 101106         |            |                      |                |          |  |
| 20" Ki Spoke                                 | 200530         | 6                  |           |                                           |                                           | Natural Fit: 100835 | 200536     |                | 200542     |                      |                |          |  |
| 22" Ki Spoke                                 | 200531         |                    |           |                                           |                                           |                     | 100560     |                | 100576     | 100569               | 200538         | 200201   |  |
| 24" Ki Spoke                                 | 200532         |                    |           |                                           |                                           |                     | 200349     |                | 100577     | 200547               | 100793         | 200202   |  |
| 25" Ki Spoke                                 | 200533         |                    |           |                                           | 200350                                    |                     |            | 101870         | 200548     | 200539               | 200540         |          |  |
| 26" Ki Spoke                                 | 200534         |                    |           |                                           | 200351                                    |                     |            | 100578         | 200549     | 100907               | 101454         |          |  |
| 20" 5-Spoke X Core                           | 101961         | 5                  |           | Aluminum, Projection, Natural Fit: 100629 | 100724                                    | 101897              |            | 101963         |            |                      |                |          |  |
| 22" 5-Spoke X Core                           | 101962         |                    |           |                                           |                                           | 101898              |            | 101964         |            |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           |                                           |                                           | 100975              |            | 100976         | 200546     |                      |                |          |  |
| 24" 5-Spoke X Core                           | 100960         |                    |           | Plastic Ctd: 101756                       | 101893                                    |                     |            |                |            | 100768               | 100769         |          |  |

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 100792 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels, part number 100666 for the screw on 5-Spoke X Core wheels, or part number 100654 for the screw on MAXX Mag wheels.

Effective 11/27/2023