



# Zwim Pool Chair

MSRP – USD  
Effective May 1, 2025

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Zwim Pool Chair** – Self-Propel Mobile Pool / Shower Chair includes 24" rear wheels, contoured seat with drain slot, 5" (125mm) Tente casters, swingaway foot supports, molded vinyl back upholstery and locking, flip-up arms. **User weight cap. 150 kg (330 lbs).**

**Need help with this order form? Call Raz customer service at (877) 720-5678**

## 1. Chair Model

ZWIM	Zwim Pool / Shower Chair	\$2,495
------	--------------------------	---------

## 2. Pelvic and Chest Belts

EZPBICL	Infection Control Pelvic Belt – Large	\$95
EZCBICASL	Infection Control Chest Belt – Large	\$95

## 3. Other

EZCSIC19	Infection Control Snap-On Calf Strap	N/C
EZLTS2	Lateral Thoracic Support w/ Large Pad (pr)	\$680
EZHSAS1	Head Support with small pad	\$445
EZATSPR	Anti-Tippers (pr)	\$110
ZQRSP	Quick-Release Axles (pr)	\$45

## Notes:

## Specifications

## Zwim Pool Chair

Weight Capacity	330 lbs (150 kg)
Seat Width	19.5" (50 cm)
Seat Depth (front of back post to front of seat)	17" (43 cm)
Back Width	18.5" (47 cm)
Back Height (from top of seat to top of upholstery)	17" (43 cm)
Arm Support Height Above Seat	9" (23 cm)
Seat-To-Floor Height* (1" increments [2.5cm])	19"-23" (48-58 cm)
Overall Width	26.5" (67 cm)
Product Weight	39.5 lbs (18 kg)
Frame Material	Aluminum