


16" Mag Wheel (Focus Drum Brake)



If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles or axle plate and wheel locks may be needed.

| Item Number | Part Number | Description | Quantity |
|-------------|-------------|---|----------|
| 1a, 2-6 | 002678 | Drum Brake Kit, 16" Mag, Pneumatic Tire PAIR | 1 |
| 1b, 2-6 | 002680 | Drum Brake Kit, 16" Mag, Full Poly Tire PAIR | 1 |
| 1c, 2-6 | 002679 | Drum Brake Kit, 16" Mag, Pneumatic Tire w/ Airless Insert PAIR | 1 |