

# 12" Mag Wheel



*If an Attendant Foot Lock is present on a Flip, please see the 12" & 16" Mag Wheel (Attendant Foot Lock) parts manual page. If Drum Brake is present, see 12" Mag Wheel (Drum Brake) page. If changing wheel size and the Attendant Foot Lock is present on a Focus, a new linkage may be required; reference the corresponding Focus CR Attendant Foot Lock page. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles or axle plate and wheel locks may be needed.*

Item Number	Part Number	Description	Quantity
1a	001424	12" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001425	12" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001426	12" Mag Wheel w/ Full Poly Assembly Tilt	1