





ICON™ TALL BACK



| Icon™ Tall Back Size (in.) | 12 | 14 | 15 | 16 | 17 | 18 | 20 | 22 | 24 |
|--|-------|-------|-------|-------|-------|-------|-------|---------------|---------------|
| A. Back Height (in.) ¹ | 14.6 | 17.0 | 18.2 | 19.4 | 20.5 | 21.75 | 23.6 | Not Available | Not Available |
| B. Bottom of back to beginning of vertical slot (in.) | 2.0 | 2.5 | 2.5 | 3.0 | 3.0 | 3.5 | 3.5 | Not Available | Not Available |
| C. Bottom of back to top of vertical slot (in.) | 9.0 | 9.5 | 10.5 | 11.0 | 11.0 | 11.5 | 11.5 | Not Available | Not Available |
| D. Range of Icon™ hardware adjustment in vertical slot (in.) | 3.5 | 3.5 | 4.5 | 4.5 | 4.5 | 4.5 | 4.5 | Not Available | Not Available |
| E. Outside width of shell (in.) | 9.5 | 11.4 | 12.4 | 13.4 | 14.4 | 15.4 | 17.4 | Not Available | Not Available |
| F. Inside width of support surface (in.) | 8.7 | 10.6 | 11.6 | 12.6 | 13.6 | 14.6 | 16.6 | Not Available | Not Available |
| Will fit wheelchair widths (in.) | 12-13 | 13-15 | 14-16 | 15-17 | 16-18 | 17-19 | 19-21 | Not Available | Not Available |

ICON™ DEEP BACK²



| Icon™ Deep Back Size (in.) | 12 | 14 | 15 | 16 | 17 | 18 | 20 | 22 | 24 |
|--|---------------|------|------|------|------|-------|------|---------------|---------------|
| A. Back Height (in.) ¹ | Not Available | 17.0 | 18.2 | 19.4 | 20.5 | 21.75 | 23.6 | Not Available | Not Available |
| B. Bottom of back to beginning of vertical slot (in.) | Not Available | 2.5 | 2.5 | 3.0 | 3.0 | 3.5 | 3.5 | Not Available | Not Available |
| C. Bottom of back to top of vertical slot (in.) | Not Available | 9.5 | 10.5 | 11.0 | 11.0 | 11.5 | 11.5 | Not Available | Not Available |
| D. Range of Icon™ hardware adjustment in vertical slot (in.) | Not Available | 3.5 | 4.5 | 4.5 | 4.5 | 4.5 | 4.5 | Not Available | Not Available |
| E. Outside width of shell (in.) | Not Available | 12.6 | 14.3 | 15.3 | 16.3 | 17.4 | 19.6 | Not Available | Not Available |
| F. Inside width of support surface (in.) | Not Available | 11.8 | 13.5 | 14.5 | 15.5 | 16.6 | 18.8 | Not Available | Not Available |
| G. Height of fixed lateral support (in.) | Not Available | 11.6 | 12.6 | 13.6 | 14.3 | 15.0 | 16.2 | Not Available | Not Available |

2. Fits the same wheelchair widths as the Icon™ Tall Back