


User Activated Anti-Tips Clik - Rogue XP



Anti-Tip Assemblies					
Wheel Size	Dump				
	1"	2"	3"	4"	5"
20"	114197	116542	116542	114198	114198
22"	114197	116542	116542	116542	114198
24"	114197	114197	116542	116542	116542
25"	114198*	114197	114197	116542	116542
26"	114197*	114197*	116542*	116542*	116542*

* Requires extended lower assembly 116544

Contact Customer Service for chairs with more than 2" of variance between Front and Rear Seat Height. Complete assemblies contain "one" complete user activated anti-tip. Quantity two is required for a pair. A pair of User Activated Anti-Tips is not available if a tipping lever is present

Item Number	Part Number	Description	Quantity
1a, 2-6	122611	UserActivated Anti-Tip Complete Assembly Straight CLIK/Rogue XP	1
1b, 2-6	122612	UserActivated Anti-Tip Complete Assembly Bent CLIK/Rogue XP	1
1c, 2-6	122613	UserActivated Anti-Tip Complete Assembly 15 Degrees CLIK/Rogue XP	1
2	004306	Anti-Tip Anti Rattle Bushing	1
3	119668	Anti-Tip Lower Assembly	1
6	100657	M5 Nylock Nut Blk Zn	1