


Arc Wheels - MAXX Spoke (Drum Brake)



If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. Please note tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity
1	111731	Wheel 20" MAXX Spoke, Drum Brake	1
2	112193	Replacement Spoke, Drum Brake Wheel, 20" Cross	1
3	112192	Replacement Spoke, Drum Brake Wheel, 20" Radial	1