


# Rogue ALX Seat Upholstery - Growth (Formerly Tsunami)



SCREW (ITEM 3) QTY PER CHAIR	
SEAT DEPTH	SCREW QTY
14" - 15"	8
16" - 18"	10
19" - 21"	12
22"	14

*Seat upholstery width measures approximately 1" less than the chair size Seat upholstery depth measures true to the chair size*

Item Number	Part Number	Description	Quantity
1a	101000	Seat Upholstery Rigid/Liberty 1414	1
1b	101005	Seat Upholstery Rigid 1415	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1c	101012	Seat Upholstery Rigid/Liberty 1416	1
1d	101019	Seat Upholstery Rigid 1417	1
1e	101026	Seat Upholstery Rigid/Liberty 1418	1
1f	101033	Seat Upholstery Rigid 1419	1
1g	101074	Seat Upholstery Rigid/Liberty 1420	1
1h	101001	Seat Upholstery Rigid/Liberty 1514	1
1i	101006	Seat Upholstery Rigid 1515	1
1j	101013	Seat Upholstery Rigid/Liberty 1516	1
1k	101020	Seat Upholstery Rigid 1517	1
1l	101027	Seat Upholstery Rigid/Liberty 1518	1
1m	101034	Seat Upholstery Rigid 1519	1
1n	101040	Seat Upholstery Rigid/Liberty 1520	1
1o	101002	Seat Upholstery Rigid/Liberty 1614	1
1p	101007	Seat Upholstery Rigid 1615	1
1q	101014	Seat Upholstery Rigid/Liberty 1616	1
1r	101021	Seat Upholstery Rigid 1617	1
1s	101028	Seat Upholstery Rigid/Liberty 1618	1
1t	101035	Seat Upholstery Rigid 1619	1
1u	101041	Seat Upholstery Rigid/Liberty 1620	1
1v	101542	Seat Upholstery Rigid 1621	1
1w	101543	Seat Upholstery Rigid 1622	1
1x	101003	Seat Upholstery Rigid/Liberty 1714	1
1y	101008	Seat Upholstery Rigid 1715	1
1z	101015	Seat Upholstery Rigid/Liberty 1716	1
1aa	101022	Seat Upholstery Rigid 1717	1
1ab	101029	Seat Upholstery Rigid/Liberty 1718	1
1ac	101036	Seat Upholstery Rigid 1719	1
1ad	101042	Seat Upholstery Rigid/Liberty 1720	1

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
1ae	101544	Seat Upholstery Rigid 1721	1
1af	101545	Seat Upholstery Rigid 1722	1
1ag	101004	Seat Upholstery Rigid/Liberty 1814	1
1ah	101009	Seat Upholstery Rigid 1815	1
1ai	101016	Seat Upholstery Rigid/Liberty 1816	1
1aj	101023	Seat Upholstery Rigid 1817	1
1ak	101030	Seat Upholstery Rigid/Liberty 1818	1
1al	101037	Seat Upholstery Rigid 1819	1
1am	101043	Seat Upholstery Rigid/Liberty 1820	1
1an	101546	Seat Upholstery Rigid 1821	1
1ao	101547	Seat Upholstery Rigid 1822	1
1ap	101010	Seat Upholstery Rigid 1915	1
1aq	101017	Seat Upholstery Rigid/Liberty 1916	1
1ar	101024	Seat Upholstery Rigid 1917	1
1as	101031	Seat Upholstery Rigid/Liberty 1918	1
1at	101038	Seat Upholstery Rigid 1919	1
1au	101044	Seat Upholstery Rigid/Liberty 1920	1
1av	101548	Seat Upholstery Rigid 1921	1
1aw	101549	Seat Upholstery Rigid 1922	1
1ax	101011	Seat Upholstery Rigid 2015	1
1ay	101018	Seat Upholstery Rigid/Liberty 2016	1
1az	101025	Seat Upholstery Rigid 2017	1
1ba	101032	Seat Upholstery Rigid/Liberty 2018	1
1bb	101039	Seat Upholstery Rigid 2019	1
1bc	101045	Seat Upholstery Rigid/Liberty 2020	1
1bd	101550	Seat Upholstery Rigid 2021	1
1be	101551	Seat Upholstery Rigid 2022	1
1bf	101528	Seat Upholstery Rigid/Liberty 2116	1

Item Number	Part Number	Description	Quantity
1bg	101529	Seat Upholstery Rigid 2117	1
1bh	101530	Seat Upholstery Rigid/Liberty 2118	1
1bi	101531	Seat Upholstery Rigid 2119	1
1bj	101532	Seat Upholstery Rigid/Liberty 2120	1
1bk	101533	Seat Upholstery Rigid 2121	1
1bl	101534	Seat Upholstery Rigid 2122	1
1bm	101535	Seat Upholstery Rigid/Liberty 2216	1
1bn	101536	Seat Upholstery Rigid 2217	1
1bo	101537	Seat Upholstery Rigid/Liberty 2218	1
1bp	101538	Seat Upholstery Rigid 2219	1
1bq	101539	Seat Upholstery Rigid/Liberty 2220	1
1br	101540	Seat Upholstery Rigid 2221	1
1bs	101541	Seat Upholstery Rigid 2222	1
2a	000241	Seat Rail 14" Deep Rigid/Liberty	2
2b	000242	Seat Rail 15" Deep RIGID	2
2c	000243	Seat Rail 16" Deep Rigid/Liberty	2
2d	000244	Seat Rail 17" Deep RIGID	2
2e	000245	Seat Rail 18" Deep Rigid/Liberty	2
2f	000246	Seat Rail 19" Deep RIGID	2
2g	000247	Seat Rail 20" Deep Rigid/Liberty	2
2h	000552	Seat Rail 21" Deep RIGID	2
2i	000553	Seat Rail 22" Deep RIGID	2
3	100734	M5 x 16 PHMS ZC	14