

# CR45 Wheels - 12" Mag



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. if changing wheel size and Attendant Foot Lock, spacer 112503 is required which can be found on the CR45 Axle Plate page*

Item Number	Part Number	Description	Quantity
1a	001424	12" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001425	12" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001426	12" Mag Wheel w/ Full Poly Assembly Tilt	1