


Flip Fixed Height with Adjustable Handle Backrest - Growth



Has your width changed? Will the push handle (#5) accommodate your new width? If no, please select the appropriate back assembly. Callout #5: 10"-13" widths receive the curved handle and 14"-18" widths receive the "L" shaped handle. If chair width is 14"-18" wide and has a serial number of FLP010784 through FLP011079, see the Discontinued Flip Fixed Height with Adjustable Handle Backrest page.

| Item Number | Part Number | Description | Quantity |
|-----------------------|-------------|--|----------|
| 3a, 4, 5a, 6, 7, 8, 9 | 112172 | Back Assembly, Fixed Height w/ Adjustable Handle, Short, 10-13" Wide, Left Flip | 1 |
| 3a, 4, 5b, 6, 7, 8, 9 | 112173 | Back Assembly, Fixed Height w/ Adjustable Handle, Short, 10-13" Wide, Right Flip | 1 |
| 3a, 4, 5a, 6, 7, 8, 9 | 112174 | Back Assembly, Fixed Height w/ Adjustable Handle, Tall, 10-13" Wide, Left Flip | 1 |
| 3a, 4, 5b, 6, 7, 8, 9 | 112175 | Back Assembly, Fixed Height w/ Adjustable Handle, Tall, 10-13" Wide, Right Flip | 1 |
| 3b, 4, 5c, 6, 7, 8, 9 | 112176 | Back Assembly, Fixed Height w/ Adjustable Handle, Short, 14-18" Wide, Left Flip <i>Used with serial numbers prior to FLP010784 and serial numbers starting with FLP011080</i> | 1 |

| Item Number | Part Number | Description | Quantity |
|-----------------------|-------------|---|----------|
| 3b, 4, 5d, 6, 7, 8, 9 | 112177 | Back Assembly, Fixed Height w/ Adjustable Handle, Short, 14-18" Wide, Right Flip <i>Used with serial numbers prior to FLP010784 and serial numbers starting with FLP011080</i> | 1 |
| 3b, 4, 5c, 6, 7, 8, 9 | 112178 | Back Assembly, Fixed Height w/ Adjustable Handle, Tall, 14-18" Wide, Left Flip <i>Used with serial numbers prior to FLP010784 and serial numbers starting with FLP011080</i> | 1 |
| 3b, 4, 5d, 6, 7, 8, 9 | 112179 | Back Assembly, Fixed Height w/ Adjustable Handle, Tall, 14-18" Wide, Right Flip <i>Used with serial numbers prior to FLP010784 and serial numbers starting with FLP011080</i> | 1 |
| 1a | 002578 | Growth Strut, Stroller, Short | 1 |
| 1a, 2 | 503107 | Growth Strut Assembly Short, Focus/Flip | 1 |
| 1b | 002579 | Growth Strut, Stroller, Long <i>Used only if growing chair from 13" wide or less to 14" wide or greater</i> | 1 |
| 2 | 101843 | M5 x 20 SHCS Blk Zn | 1 |