


## Arc Wheels - MAXX Spoke (Drum Brake)



*If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. Please note tab mounted handrims must be installed prior to the tire being installed.*

Item Number	Part Number	Description	Quantity
1	111731	Wheel 20" MAXX Spoke, Drum Brake	1
2	112193	Replacement Spoke, Drum Brake Wheel, 20" Cross	1
3	112192	Replacement Spoke, Drum Brake Wheel, 20" Radial	1