



Raz-ART

Attendant Rotational Tilt

Mobile Shower Commode Chair

The **Raz-ART** utilizes a patented, dual-arc, orbital, tilt system that allows for the smallest footprint and lowest tilt effort of any chair on the market. A unique frame design provides two ranges of tool-less height adjustment and a 40° tilt range.



FEATURES

Patented Tilt Geometry:

- Dual arc, orbital tilt, minimizes center-of-mass movement throughout range
- Extremely low tilt effort
- 40° tilt range

Small Footprint:

Provides unparalleled manoeuvrability

Patent-Pending Tilt Lock:

Eliminates cable looping

Lightweight, invertible, aluminum frame:

Provides two tool-less height ranges totaling 7" (18 cm)



Flip-up arm supports with optional locks

Adjustable back angle

Compatible with all standard-width Raz seats, which include the exclusive IPAS

Backrest with 3" (7.5 cm) fore/aft adjustment to vary seat depth

Built-in tilt angle indicator

Dual-locking (swivel and roll) 5" (125 mm) casters

Extensive range of options and accessories is available

SPECIFICATIONS	Raz-ART (Attendant propel Rotational Tilt)	
WEIGHT CAPACITY	350 lbs (160 kg)	
OVERALL WIDTH	22.5" (57 cm)	
OVERALL LENGTH – with head support and standard foot supports in middle position	Upright: 35" (89 cm) / Tilted to 40°: 50.5" (128 cm)	
SEAT WIDTH	18" (46 cm)	
SEAT DEPTH – front of back post to front of seat	16" – 18" (40.5 – 46 cm) or 19" – 21" (48 – 53 cm) ¹	
SEAT-TO-FLOOR HEIGHT – 1" (2.5 cm) increments – above standard molded seat	Lower Range: 20" – 23" (51 – 58.5 cm)	Higher Range: 24" – 27" (61 – 68.5 cm)
CLEARANCE HEIGHT UNDER SEAT – 1" (2.5 cm) increments		
– with pan brackets	Lower Range: 15" – 18" (38 – 46 cm)	Higher Range: 19" – 22" (48 – 56 cm)
– without pan brackets	Lower Range: 15.5" – 18.5" (41 – 49 cm)	Higher Range: 19.5" – 22.5" (51 – 59 cm)
ARMREST HEIGHT – above standard molded seat		
– 18"H Fixed Back	9" (23 cm)	
– 15"H Adjusta-Back	5.5" – 8.5" (14 – 21.5 cm) in 0.5" (1.25 cm) increments	
– 17"H Adjusta-Back	6.5" – 9.5" (16.5 – 24 cm) in 0.5" (1.25 cm) increments	
– 20"H Adjusta-Back	7.5" – 10.5" (19 – 27 cm) in 0.5" (1.25 cm) increments	
– 24"H Adjusta-Back	8.5" – 11.5" (21.5 – 29 cm) in 0.5" (1.25 cm) increments	
PRODUCT WEIGHT (with foot supports)	69 lbs (31 kg)	
FRAME MATERIAL – base frame and seat frame	Powder-Coated Aluminum ²	
FRAME MATERIAL – back frame and footrests	Stainless Steel ²	
HARDWARE MATERIAL	Stainless Steel ²	

¹ Ischial Pelvic Alignment System (IPAS) allows for the following seat adjustments: 2" (5 cm) of fore/aft, 1" (2.5cm) lateral, or 4.5° rotational

² Lifetime warranty against corrosion

A 2-year limited warranty covers the frame against manufacturing defects. A 1-year limited warranty covers all other components, subject to normal wear and tear exclusions. *Specifications and options are subject to change without notice.*