


Arc Wheels - 16" Mag (Non-Drum Brake)



If a drum brake is present, see 16" Mag Wheel (Arc Drum Brake) page for correct wheel assembly. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed.

Item Number	Part Number	Description	Quantity
1a	001159	16" Mag Wheel w/ Pneumatic Assembly Tilt	1
1b	001160	16" Mag Wheel w/ Pneumatic w/ Airless Insert Assembly Tilt	1
2	001161	16" Mag Wheel w/ Full Poly Assembly Tilt	1
3	003453	16" Mag Wheel w/ Low Poly Assembly Tilt	1