


# MAXX Spoke Wheel / Tire / Handrim Kits

---



*Not all components of the Wheel, Tire and Handrim kits will be listed. These kits cannot be used if Drum Brake is present. If changing wheel size and anti-tips are present, new anti-tips may be required. If changing wheel size or style, new axles may be needed. For the Catalyst 4, new axle plates may be needed.*

<b>Item Number</b>	<b>Part Number</b>	<b>Description</b>	<b>Quantity</b>
	110979	18" MAXX Spoke, Pneumatic Tire (Gray), Aluminum Anodized Handrim Kit	1
	110980	18" MAXX Spoke, Pneumatic Tire (Gray), Aluminum Anodized Handrim Kit w/ Micro Grip	1
	110981	18" MAXX Spoke, Pneumatic Tire (Gray), Plastic Coated Handrim Kit	1
	110982	18" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit	1
	110983	18" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	110984	18" MAXX Spoke, Low Poly Tire, Plastic Coated Handrim Kit	1
	110985	18" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit	1
	110986	18" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	110987	18" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit	1
	110988	18" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Gray), Aluminum Anodized Handrim Kit	1
	110989	18" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Gray), Aluminum Anodized Handrim Kit w/ Micro Grip	1
	110990	18" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Gray), Plastic Coated Handrim Kit	1
	118156	20" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118157	20" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118158	20" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	110994	20" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit	1
	110995	20" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	110996	20" MAXX Spoke, Low Poly Tire, Plastic Coated Handrim Kit	1
	110997	20" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	110998	20" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	110999	20" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1
	111000	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Gray), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118159	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118160	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118161	20" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	111003	20" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	111004	20" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	111005	20" MAXX Spoke, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	111006	20" MAXX Spoke, High Pressure, Aluminum Anodized Handrim Kit	1
	111007	20" MAXX Spoke, High Pressure, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111008	20" MAXX Spoke, High Pressure, Plastic Coated Handrim Kit	1
	111009	20" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit	1
	111010	20" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111011	20" MAXX Spoke, Shox Tire, Plastic Coated Handrim Kit	1
	109790	22" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	109908	22" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	119329	22" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120089	22" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120090	22" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119332	22" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120092	22" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119333	22" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120093	22" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109796	22" MAXX Spoke, Full Poly Tire, Natural Fit LT Handrim Omit Thumb Grip Kit	1
	109793	22" MAXX Spoke, Full Poly Tire, Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	109791	22" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	109792	22" MAXX Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	109815	22" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit	1
	109911	22" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	119347	22" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120101	22" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120102	22" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119350	22" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120104	22" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120105	22" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109820	22" MAXX Spoke, High Pressure Tire, Natural Fit LT Handrim Omit Thumb Grip Kit	1
	109817	22" MAXX Spoke, High Pressure Tire, Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	109816	22" MAXX Spoke, High Pressure Tire, Plastic Coated Handrim Kit	1
	109782	22" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit	1
	109907	22" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	119323	22" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120083	22" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120084	22" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119326	22" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120086	22" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120087	22" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109787	22" MAXX Spoke, Low Poly Tire, Natural Fit LT Handrim Omit Thumb Grip Kit	1
	109784	22" MAXX Spoke, Low Poly Tire, Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	109783	22" MAXX Spoke, Low Poly Tire, Plastic Coated Handrim Kit	1
	118189	22" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118206	22" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	118190	22" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118191	22" MAXX Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	119474	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120216	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120217	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119477	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120219	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120220	22" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	118195	22" MAXX Spoke, Pneumatic Tire (Black), Natural Fit LT Handrim Omit Thumb Grip Kit	1
	118192	22" MAXX Spoke, Pneumatic Tire (Black), Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	118102	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount	1
	118103	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	119465	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120207	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119466	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120208	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119462	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120204	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119463	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	120205	22" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	118104	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Natural Fit LT Handrim Omit Thumb Grip Kit	1
	118107	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	118110	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118111	22" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	109890	22" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount	1
	109910	22" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip	1
	119395	22" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120137	22" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120138	22" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119398	22" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120140	22" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120141	22" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109895	22" MAXX Spoke, Sentinel Tire, Natural Fit LT Handrim Omit Thumb Grip Kit	1
	109892	22" MAXX Spoke, Sentinel Tire, Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	109891	22" MAXX Spoke, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	109807	22" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit	1
	109912	22" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	119341	22" MAXX Spoke Shox Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120095	22" MAXX Spoke Shox Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	120096	22" MAXX Spoke Shox Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119344	22" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120098	22" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1

Item Number	Part Number	Description	Quantity
	120099	22" MAXX Spoke Shox Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109812	22" MAXX Spoke, Shox Tire, Natural Fit LT Handrim Omit Thumb Grip Kit	1
	109809	22" MAXX Spoke, Shox Tire, Natural Fit Standard Handrim Omit Thumb Grip Kit	1
	109808	22" MAXX Spoke, Shox Tire, Plastic Coated Handrim Kit	1
	120829	24" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120835	24" MAXX Spoke, Full Poly Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109840	24" MAXX Spoke, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	109915	24" MAXX Spoke, Full Poly Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	119364	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120112	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119365	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120113	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119366	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120114	24" MAXX Spoke Full Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119367	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120115	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119368	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120116	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119369	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120117	24" MAXX Spoke Full Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109841	24" MAXX Spoke, Full Poly Tire, Plastic Coated Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	109842	24" MAXX Spoke, Full Poly Tire, Projection Handrim Kit, Standard Mount	1
	120831	24" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120837	24" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109874	24" MAXX Spoke, High Pressure Tire, Hard Aluminum Anodized Handrim Kit <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1
	109919	24" MAXX Spoke, High Pressure Tire, Hard Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1
	119382	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120124	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119383	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120125	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119384	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120126	24" MAXX Spoke High Pressure Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119385	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120127	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119386	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120128	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119387	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120129	24" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109875	24" MAXX Spoke, High Pressure Tire, Plastic Coated Handrim Kit	1
	120828	24" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120834	24" MAXX Spoke, Low Poly Tire, Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1



Item Number	Part Number	Description	Quantity
	109832	24" MAXX Spoke, Low Poly Tire, Hard Aluminum Anodized Handrim Kit <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	109914	24" MAXX Spoke, Low Poly Tire, Hard Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	119358	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120106	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119359	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120107	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119360	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120108	24" MAXX Spoke Low Poly Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119361	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120109	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119362	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120110	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119363	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120111	24" MAXX Spoke Low Poly Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109833	24" MAXX Spoke, Low Poly Tire, Plastic Coated Handrim Kit	1
	120846	24" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120847	24" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	118199	24" MAXX Spoke, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	118215	24" MAXX Spoke, Pneumatic Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	118200	24" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1

Item Number	Part Number	Description	Quantity
	118201	24" MAXX Spoke, Pneumatic Tire (Black), Projection Handrim Kit, Standard Mount	1
	119479	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120221	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119480	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120222	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119481	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120223	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119482	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120224	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119483	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120225	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119484	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120226	24" MAXX Spoke Pneumatic Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	120844	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120845	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	118112	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	118113	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	119470	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	120214	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119471	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120213	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119472	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120212	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119467	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120211	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119468	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120210	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119469	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120209	24" MAXX Spoke Pneumatic w/ Airless Insert Tire (Black) Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	118120	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Plastic Coated Handrim Kit, Standard Mount	1
	118121	24" MAXX Spoke, Pneumatic w/ Airless Insert Tire (Black), Projection Handrim Kit, Standard Mount	1
	120832	24" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120838	24" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109882	24" MAXX Spoke, Schwalbe Marathon Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1
	109920	24" MAXX Spoke, Schwalbe Marathon Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1

Item Number	Part Number	Description	Quantity
	119388	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120130	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119389	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120131	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119390	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120132	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119391	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120133	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119392	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120134	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119393	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120135	24" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109883	24" MAXX Spoke, Schwalbe Marathon Tire, Plastic Coated Handrim Kit, Standard Mount	1
	120840	24" MAXX Spoke, Schwalbe One Tire, Aluminum Anodized Handrim Kit <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120841	24" MAXX Spoke, Schwalbe One Tire, Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109922	24" MAXX Spoke, Schwalbe One Tire, Hard Aluminum Anodized Handrim Kit <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	109930	24" MAXX Spoke, Schwalbe One Tire, Hard Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	119406	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120148	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1

Item Number	Part Number	Description	Quantity
	119407	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120149	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119408	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120150	24" MAXX Spoke Schwalbe One Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119409	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120151	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119410	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120152	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119411	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120153	24" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109923	24" MAXX Spoke, Schwalbe One Tire, Plastic Coated Handrim Kit	1
	120833	24" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120836	24" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109898	24" MAXX Spoke, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	109917	24" MAXX Spoke, Sentinel Tire, Hard Aluminum Anodized Handrim Kit, Standard Mount w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Klik, and Flip</i>	1
	119400	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120142	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119401	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120143	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119402	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	120144	24" MAXX Spoke Sentinel Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119403	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120145	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119404	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120146	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119405	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120147	24" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109899	24" MAXX Spoke, Sentinel Tire, Plastic Coated Handrim Kit, Standard Mount	1
	120830	24" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	120839	24" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 4, Cat5Vx, Cat E, Rogue ALX, Little Wave, Focus, and Liberty</i>	1
	109857	24" MAXX Spoke, Shox Tire, Hard Aluminum Anodized Handrim Kit <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1
	109921	24" MAXX Spoke, Shox Tire, Hard Aluminum Anodized Handrim Kit w/ Micro Grip <i>Used on Cat 5, Spark, Ethos, Rogue2, Rogue XP, Clik, and Flip</i>	1
	119376	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Omit Thumb Grip Kit	1
	120118	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119377	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Regular Thumb Grip Kit	1
	120119	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119378	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Extra Thumb Grip Kit	1
	120120	24" MAXX Spoke Shox Tire Ergonomic Standard Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	119379	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120121	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119380	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120122	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119381	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1

Item Number	Part Number	Description	Quantity
	120123	24" MAXX Spoke Shox Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	109858	24" MAXX Spoke, Shox Tire, Plastic Coated Handrim Kit	1
	111467	24" MAXX Spoke, Mountain Tire, Plastic Coated Handrim Kit	1
	111022	25" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit	1
	111023	25" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111024	25" MAXX Spoke, Sentinel Tire, Plastic Coated Handrim Kit	1
	119422	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120164	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119423	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120165	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119424	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120166	25" MAXX Spoke Sentinel Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	111031	25" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit	1
	111032	25" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111033	25" MAXX Spoke, High Pressure Tire, Plastic Coated Handrim Kit	1
	119428	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120170	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119429	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120171	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119430	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120172	25" MAXX Spoke High Pressure Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	111040	25" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit	1
	111041	25" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111042	25" MAXX Spoke, Schwalbe Marathon Tire, Plastic Coated Handrim Kit	1
	119434	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120176	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1

Item Number	Part Number	Description	Quantity
	119435	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120177	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119436	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120178	25" MAXX Spoke Schwalbe Marathon Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	111049	25" MAXX Spoke, Schwalbe One Tire, Aluminum Anodized Handrim Kit	1
	111050	25" MAXX Spoke, Schwalbe One Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111051	25" MAXX Spoke, Schwalbe One Tire, Plastic Coated Handrim Kit	1
	119440	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120182	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119441	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120183	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119442	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120184	25" MAXX Spoke Schwalbe One Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	111058	25" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit	1
	111059	25" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111060	25" MAXX Spoke, Shox Tire, Plastic Coated Handrim Kit	1
	119446	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Omit Thumb Grip Kit	1
	120188	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Omit Thumb Grip Kit w/ Micro Grip	1
	119447	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit	1
	120189	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Regular Thumb Grip Kit w/ Micro Grip	1
	119448	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Extra Thumb Grip Kit	1
	120190	25" MAXX Spoke Shox Tire Ergonomic Lite Handrim Extra Thumb Grip Kit w/ Micro Grip	1
	120629	26" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit	1
	120630	26" MAXX Spoke, Pneumatic Tire (Black), Aluminum Anodized Handrim Kit w/ Micro Grip	1
	120631	26" MAXX Spoke, Pneumatic Tire (Black), Plastic Coated Handrim Kit	1



Item Number	Part Number	Description	Quantity
	120632	26" MAXX Spoke, Pneumatic Tire (Black), Projection Handrim Kit	1
	120633	26" MAXX Spoke, Pneumatic (Black) w/ Airless Insert Tire, Aluminum Anodized Handrim Kit	1
	120634	26" MAXX Spoke, Pneumatic (Black) w/ Airless Insert Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	120635	26" MAXX Spoke, Pneumatic (Black) w/ Airless Insert Tire, Plastic Coated Handrim Kit	1
	120636	26" MAXX Spoke, Pneumatic (Black) w/ Airless Insert Tire, Projection Handrim Kit	1
	111087	26" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit	1
	111088	26" MAXX Spoke, Sentinel Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111089	26" MAXX Spoke, Sentinel Tire, Plastic Coated Handrim Kit	1
	111096	26" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit	1
	111097	26" MAXX Spoke, High Pressure Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111098	26" MAXX Spoke, High Pressure Tire, Plastic Coated Handrim Kit	1
	111105	26" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit	1
	111106	26" MAXX Spoke, Schwalbe Marathon Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111107	26" MAXX Spoke, Schwalbe Marathon Tire, Plastic Coated Handrim Kit	1
	111114	26" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit	1
	111115	26" MAXX Spoke, Shox Tire, Aluminum Anodized Handrim Kit w/ Micro Grip	1
	111116	26" MAXX Spoke, Shox Tire, Plastic Coated Handrim Kit	1