

## Catalyst 5 Casters - 1.25" Soft Roll (Poly) Casters



If changing size of casters, new forks may be needed.

| Item Number | Part Number | Description                                                      | Quantity |
|-------------|-------------|------------------------------------------------------------------|----------|
| 1a          | 100809      | Caster 4 x 1.25" Soft Roll<br><i>Replaced part number 100517</i> | 1        |
| 1b          | 100807      | Caster 5 x 1.25" Soft Roll<br><i>Replaced part number 100518</i> | 1        |
| 1c          | 100786      | Caster 6 x 1.25" Soft Roll<br><i>Replaced part number 100521</i> | 1        |
| 1d          | 111862      | Caster 7 x 1.25" Soft Roll<br><i>Replaced part number 100556</i> | 1        |
| 1e          | 111863      | Caster 8 x 1.25" Soft Roll<br><i>Replaced part number 100525</i> | 1        |