


# Liberty 1.25" Soft Roll (Poly) Casters



These casters effective February 8th, 2021 If replacing discontinued casters, a quantity of 2 will be needed. New casters do not match the sizes of the discontinued casters. If changing size of casters, new forks may be needed.

Item Number	Part Number	Description	Quantity
1a	100809	Caster 4 x 1.25" Soft Roll <i>Replaced part number 100517</i>	1
1b	100807	Caster 5 x 1.25" Soft Roll <i>Replaced part number 100518</i>	1
1c	100786	Caster 6 x 1.25" Soft Roll <i>Replaced part number 100521</i>	1
1d	111863	Caster 8 x 1.25" Soft Roll <i>Replaced part number 100525</i>	1