


## Catalyst 5Ti Footplates - Composite



Sizes for footplates are based on seat width. Only compatible with extension tube hangers. If changing footplate style, new extension tubes may be needed.

| Item Number | Part Number | Description                               | Quantity |
|-------------|-------------|-------------------------------------------|----------|
| 2-4         | 112895      | Hardware Assembly for Composite Footplate | 1        |
| 1a          | 100537      | Footplate Composite 14-15" Wide CAT/TILT  | 1        |
| 1b          | 100540      | Footplate Composite 16-17" Wide CAT/TILT  | 1        |
| 1c          | 100542      | Footplate Composite 18-24" Wide CAT/TILT  | 1        |
| 2           | 000053      | Saddle 3/4"                               | 2        |
| 3           | 101849      | M6 x 65 SHCS Blk Zn                       | 1        |
| 4           | 100658      | M6 Nylock Nut Blk Zn                      | 1        |