


Arc Stroller Handle Backrest - Growth



Handle Strut (Item 1)			
Part No.	Growth Range	Used on Widths	Length
002032	10"-12"	10"	4.6"
002033	11"-13"	11"-13"	5.6"
002578	14"-20"	14"-18"	8.8"

Is your width changing? Do you have stroller handle extensions for your stroller handle backrest? Will your new width accommodate the stroller handle growth strut? If no, please select the appropriate kit. When changing backrest types, check to see if new cables and tilt triggers are needed for both foot and hand tilts.

Item Number	Part Number	Description	Quantity
1a, 2	113922	Stroller Handle Backrest Growth Strut, 10-12" Wide, Arc	1
1b, 2	113923	Stroller Handle Backrest Growth Strut, 11-13" Wide, Arc	1

Item Number	Part Number	Description	Quantity
1c, 2	113924	Stroller Handle Backrest Growth Strut, 14-18" Wide, Arc	1