


Rogue Caster Mount



Item Number	Part Number	Description	Quantity	Retail
1-7	501225	Caster Pinion Hardware Assm, Rogue	1	\$103.97
1-7	501226	Caster Pinion Hardware Assm PAIR, Rogue	1	\$207.93
1	000964	Caster Mount Pinion Steel RG	1	\$61.34
2	000933	Caster Adjust CG RG	1	\$25.76
3	000934	Pinion Housing Bolt Plate RG	1	\$6.69

Item Number	Part Number	Description	Quantity	Retail
4	000932	Caster Pinion Cap RG	2	\$2.52
5	105037	M8x16 CLS10.9 ISO7380 BTN SHCS ENVRBLK W/PRE5	2	\$1.58
6	102304	M5 x 6mm FHCS Blk Zc w/Patch	2	\$3.15
7	102303	M5 x 30mm FHCS Blk Zc w/Patch	2	\$2.63