

SEDEO® ERGO™	SEDEO® PRO™	SEDEO® LITE™	CAPTAIN SEAT
			
A highly specialized and technical seating system, the SEDEO ERGO allows for precision fitting for complex clinical needs. Biometric Repositioning and a full range of programmable seating positions make SEDEO ERGO the best ergonomic seating from QUICKIE.	Truly modular in design, the SEDEO PRO seating system provides a new level of function, comfort and design. For moderate to complex clinical needs, the SEDEO PRO combines a highly adjustable seat frame with an excellent selection of postural support options.	The SEDEO LITE is an adjustable seat frame with a sleek look. Its minimal design alleviates weight capacity concerns when using aftermarket equipment. Traditional back canes are width-adjustable from 15" to 20" and easily fit JAY® backs. An array of options cover mild to moderate clinical needs, and intuitive adjustability meets individual requirements.	Designed for comfort and convenience, the Captain Seat adjusts in seat depth, legrest placement, and armrest width. The convertible seat pan allows alternate cushion selection, and the backrest has both semi-recline and fold-down capability.

CLINICAL NEED	HIGH	MED TO HIGH	MED TO HIGH	LOW
Q700 M COMPATIBILITY	✓	✓	✗	✓
Q500 M COMPATIBILITY	✗	✓	✓	✓
Q500 H COMPATIBILITY	✗	✓	✗	✗
JAY® SYNCRA™ COMPATIBILITY	✗	✓	✗	✗
BIOMETRIC REPOSITIONING	✓	✗	✗	✗
REMOTE SEATING APP	✓	✗	✗	✗
MEMORY SEATING	✓	✗	✗	✗
PRESSURE RELIEF ALERT REMINDER (PRAR)	✓	✗	✗	✗
POWER RECLINE w/7" of SHEAR (165°)	✓	✗	✗	✗
POWER RECLINE w/4" of SHEAR (150°)	✗	✓	✗	✗
SEMI POWER RECLINE OPTION (130°)	✗	✓	✗	✗
MANUAL RECLINE OPTION	✗	✓	✗	✓
POWER 50° CENTER GRAVITY TILT	✓	✓	✓	✓
POWER LEGS	✓	✓	✓	✓
12" LIFT	✓	✓	✓	✓
SEAT WIDTH	16" to 22"	16" to 22"	15" to 20"	16", 18", 20"
SEAT DEPTH	15" to 20"	14" to 22"	16" to 22"	16" to 22"
BACKREST HEIGHT	21" to 28"	21" to 30"	18" to 30"	25"

ELEVATE AND DRIVE WITH C-ME® AVAILABLE ON SEDEO ERGO, SEDEO PRO, AND SEDEO LITE










No more waiting around or driving at reduced speed when elevated. C-Me enables you to elevate up to 12" while traveling up to 3 mph! (Weight limitations may apply.)* You can lift or descend at the touch of a button with no need to stop while transitioning. Just rise-up and go at the same time!



A COMPREHENSIVE
GUIDE TO SELECTING
THE RIGHT QUICKIE
POWER WHEELCHAIR

**POWER WHEELCHAIR
REFERENCE GUIDE**



MODEL CODES	Q700 M						Q500 M						Q500 H			Q400 M						P-222 SE																															
1st Digit	MID WHEEL DRIVE						MID WHEEL DRIVE						HYBRID (MID/REAR) DRIVE			MID WHEEL DRIVE						REAR WHEEL DRIVE																															
N- No Power Seating																																																					
2nd Digit	K0848		K0850		K0868		K0861		K0862		K0884		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0868																				
3rd Digit	K0856		K0858		K0877		K0861		K0862		K0884		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0868																				
C- Group 3	K0861		K0862		K0884		K0861		K0862		K0884		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0868																				
D- Group 4	K0861		K0862		K0884		K0861		K0862		K0884		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0868																				
HPCPS CODES	K0848		K0850		K0868		K0861		K0862		K0884		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0849		K0868		K0848		K0868																				
MODELS	GRP 3		GRP 3 HD		GRP 4		GRP 3		GRP 3 HD		GRP 4		GRP 3		GRP 3		GRP 3		GRP 3		GRP 3 HD		GRP 4		GRP 3		GRP 3		GRP 3		GRP 4																						
SEAT TYPE	NPC		NPC HD		NPD		MPC		MPC HD		MPD		NCC		NCD		NPC		NCC		NPC		NPD		NPC		NCC		NPC		BD (Basic Rehab Seat)																						
SEAT TYPE	SPC		SPC HD		SPD		SCC		SCD		MPC		SCC		SCD		SPC		SCC		SPC		SPD		SPC		SCC		SD (Single Power)																								
SEAT TYPE	MPC		MPC HD		MPD		MPC		MPC HD		MPD		MPC		MPC HD		MPD		MPC		MPC HD		MPD		MPC		MPC HD		MPD																								
SEAT TYPE	PRO		ERGO		CAPTAIN / SOLID SEAT		PRO		CAPTAIN / SOLID SEAT		LITE		PRO		PRO		CAPTAIN / SOLID SEAT		LITE		REHAB		REHAB		REHAB		REHAB		REHAB																								
STD MAX SPEED	5 mph, 6 mph, 8 mph		5 mph, 6 mph, 8 mph		6 mph, 8 mph		6 mph		6 mph		6 mph		5 mph, 6 mph, 8 mph		4.5 mph		4.5 mph		4.5 mph		8.5 mph		8.5 mph		8.5 mph		8.5 mph		8.5 mph																								
USER WEIGHT CAPACITY	250 lbs*, 300 lbs, 350 lbs HD, 450 lbs HD		250 lbs*, 300 lbs, 350 lbs HD		300 lbs		300 lbs		300 lbs		300 lbs		300 lbs, 350 lbs HD		300 lbs		300 lbs		300 lbs		250 lbs, 350 lbs		250 lbs, 350 lbs		250 lbs, 350 lbs		250 lbs, 350 lbs		250 lbs, 350 lbs																								
ELECTRONICS	R-NET		R-NET		R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		VR-2 AND R-NET		R-NET		R-NET		R-NET		R-NET		R-NET																								
SUSPENSION	SPIDERTRAC 2.0		SPIDERTRAC 2.0		SPIDERTRAC 2.0		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL		SIX-WHEEL																								
TURNING RADIUS ¹	19.8"		19.8"		19.8"		18.5"		18.5"		18.5"		25.5"		18.5"		18.5"		18.5"		25"		25"		25"		25"		25"																								
SEAT WIDTH	16" to 22"		16" to 22"		16" to 20"		16" to 22"		16" to 20"		15" to 20"		16" to 22"		16" to 22"		16" to 20"		15" to 20"		14" to 24"		14" to 24"		14" to 24"		14" to 24"		14" to 24"																								
SEAT DEPTH	14" to 22"		15" to 20"		16" to 22"		14" to 22"		16" to 22"		16" to 22"		14" to 22"		14" to 22"		16" to 22"		16" to 22"		14" to 20"		14" to 20"		14" to 20"		14" to 20"		14" to 20"																								
SEAT-TO-FLOOR HEIGHT	16.5" to 19"		17" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		16.5" to 19"		18" to 20"		18" to 20"		18" to 20"		18" to 20"		18" to 20"																								
POWERED SEATING	NPC / NPD: NO POWER OPT. SPC / SPD: CHOOSE ONE MPC / MPD: TILT & 150° RECLINE		MPC / MPD: TILT & RECLINE		NCC / NCD: NO POWER OPT. SCC / SCD: TILT OR ELEVATE		NPC: NO POWER OPT. SPC: CHOOSE ONE MPC: TILT & 150° RECLINE		NCC: NO POWER OPT. SCC: TILT OR ELEVATE		NPC: NO POWER OPT. SPC: TILT OR ELEVATE		NPC / NPD: NO POWER OPT. SPC / SPD: CHOOSE ONE MPC / MPD: TILT & 150° RECLINE		NPC: NO POWER OPT. SPC: CHOOSE ONE MPC: TILT & 150° RECLINE		NCC: NO POWER OPT. SCC: TILT OR ELEVATE		NPC: NO POWER OPT. SPC: TILT OR ELEVATE		BD: NO POWER OPT. SD: TILT ONLY		BD: NO POWER OPT. SD: TILT ONLY		BD: NO POWER OPT. SD: TILT ONLY		BD: NO POWER OPT. SD: TILT ONLY		BD: NO POWER OPT. SD: TILT ONLY																								
TILT	POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 50°		POSTERIOR 30°, 50°		POSTERIOR 30°, 50°		POSTERIOR 30°, 50°		POSTERIOR 53°		POSTERIOR 53°		POSTERIOR 53°		POSTERIOR 53°		POSTERIOR 53°																								
RECLINE	150° POWERED RECLINE w/4" SHEAR REDUCTION OR 130° POWERED SEMI-RECLINE ²		POWERED RECLINE w/ 7" SHEAR REDUCTION		25° MANUAL ADJUSTABLE RECLINE		150° POWERED RECLINE w/4" SHEAR REDUCTION OR 130° POWERED SEMI-RECLINE ²		25° MANUAL ADJUSTABLE RECLINE		X		150° POWERED RECLINE w/4" SHEAR REDUCTION OR 130° POWERED SEMI-RECLINE ²		150° POWERED RECLINE w/4" SHEAR REDUCTION OR 130° POWERED SEMI-RECLINE ²		25° MANUAL ADJUSTABLE RECLINE		X		X		X		X		X																										
SEAT ELEVATOR	12", C-Me		12", C-Me		12"		12", C-Me		12"		12", C-Me		12"		12", C-Me		12"		12", C-Me		X		X		X		X		X																								
LEGREST	FIXED CENTER MOUNT POWER CENTER MOUNT SWING-AWAY		FIXED CENTER MOUNT POWER CENTER MOUNT		FIXED CENTER MOUNT POWER CENTER MOUNT		FIXED CENTER MOUNT POWER CENTER MOUNT SWING-AWAY		FIXED CENTER MOUNT POWER CENTER MOUNT		FIXED CENTER MOUNT POWER CENTER MOUNT ANGLE-ADJ. CENTER MOUNT SWING-AWAY		SWING-AWAY		FIXED CENTER MOUNT POWER CENTER MOUNT SWING-AWAY		FIXED CENTER MOUNT POWER CENTER MOUNT		FIXED CENTER MOUNT POWER CENTER MOUNT ANGLE-ADJ. CENTER MOUNT SWING-AWAY		SWING-AWAY FIXED		SWING-AWAY FIXED		SWING-AWAY FIXED		SWING-AWAY FIXED																										
VENT TRAY	✓		✓		X		✓		X		✓		✓		✓		X		✓		X		X		X		X		X																								
BATTERIES	GROUP 24 GROUP 34		GROUP 24 GROUP 34		GROUP 24 GROUP 34		GROUP 34 GROUP 22		GROUP 34 GROUP 22		GROUP 22 GROUP 34		GROUP 24 GROUP 34		GROUP 34 GROUP 22		GROUP 34 GROUP 22		GROUP 22 GROUP 34		GROUP 24		GROUP 24		GROUP 24		GROUP 24		GROUP 24																								
DRIVE WHEELS	14"		14"		14"		13"		13"		13"		13", 14"		13"		13"		13"		13"		13"		12"		12"		12"																								
OVERALL BASE WIDTH	24.5"		24.5"		25"		24"		24"		24"		24.5"		24"		24"		24"		24.5"		24.5"		24.5"		24.5"		24.5"																								
OVERALL LENGTH	36.5"		36.5"		36.5"		33"		33"		33"		36"		33"		33"		33"		31.5"		31.5"		31.5"		31.5"		31.5"																								
OCCUPIED TRANSIT TESTED	✓		✓		✓		✓		✓		✓		✓		✓		✓		✓		✓		✓		✓		✓		✓																								

1. This measurement excludes footplates.

2. The 130° Powered Semi-Recline does not meet the criteria for HCPCS code E1004 or E1007.