PRODUCT CATALOG





CONTENTS

How To Order	3	Li'l Kiddos Pelvic/Thigh Pads	
SEATING SYSTEMS		Modifications	
Seating System Packages	/	Bracket Covers	
Li'l Kiddos Seating Systems		Li'l Kiddos Lateral Bracket Covers	
Li i kiddos deating dysteins	0	Fixed Brackets	
SEATS		Adjustable Bracket And Modular Brackets	
Insert Packages	7	Modular Brackets	
Li'l Kiddos Seat Packages	8	Modular Bracket Components	
Solid Inserts	9	Modular Bracket Components	
Universal Amputee Inserts	10	Li'l Kiddos Lateral Brackets	
Wood Bases	10	Li'l Kiddos Modular Hardware Components	
Modifications & Covers	11	Medial Support Pelvic/Thigh	
		Support Packages	40
BACKS		MEDIAL SUPPORT	
Insert Packages		Pelvic/Thigh	
Li'l Kiddos Back Packages		Support Packages	Д 1
Solid Inserts		-	
Modifications & Covers	15	Li'l Kiddos Medial Thigh Support	
FOAM		Pads	
	10	Modifications	
TypesCustom Molded – Foam Form		Brackets	
	17	Bracket Components	44
HARDWARE		HEAD SUPPORT	
Seat To Back Brackets		Packages	
Universal Seating Mounting Brackets	19	Pads	45
Universal Back Mounting Brackets	19	Modifications	45
Universal Components	20	Hardware	46
E–Z Mount Brackets	21	Hardware Components	46
Bi-Angular Back Brackets	21	Trunk & Pelvic Support Straps	47
J & L Brackets	22	Soft Support Components & Hardware	47
J & L Components		LOWER CURRORT	
Li'l Kiddos Mobility Base Attaching Hardware	24	LOWER SUPPORT	
LATERAL CURRORT		Lower Extremity Support	48
LATERAL SUPPORT		UPPER SUPPORT	
Lateral Thoracic Packages		Arm & Wrist Support	50
Lateral Thoracic Packages		Multi-Adjustable Arm	
Lateral Pelvic/Armrest Packages		Platforms (Trays)	
Lateral Pelvic/Thigh Packages		Platform Mounting Hardware	
Li'l Kiddos Lateral		Platform Accessories	
Thoracic Support Packages		Platform Modifications	
Li'l Kiddos Lateral		Miscellaneous	
Pelvic/Thigh Support Packages	29	IVIISCEIIdIIEUUS	33
Thoracic Pads	30		
Pelvic/Thigh Pads	30		

Li'l Kiddos Thoracic Pads31

ABOUT AES


AES. AES custom seating solutions and hardware address basic to complex client needs. Adjustable and expandable, the products are fabricated to specific client measurements and are available in a variety of options and colors.

AES Lite. AES Lite offers an easy-to-order, simplified mix of AES options.

AES Li'l Kiddos. AES Li'l Kiddos pediatric custom seating solutions and hardware are designed specifically for children who are using wheelchairs with frames 14" wide and under.

HOW TO ORDER

- 1. Complete AES Custom Seating Order Form
 - a. Include wheelchair information. including model, tubing size, frame width/depth, etc.
 - b. Note the client's measurements
 - c. Indicate desired seat dimensions



- d. Select your desired Seating System Package or Individual Seat Package and/or Back Package
- e. Select from all of our unique modifications to customize the size. shape, foam, and countless others to fit your client's unique needs.
 - i. Select additional options, including lateral, medial, anterior trunk, pelvic, and/or knee supports.
- 2. Call, email, or fax in your order to Sunrise Medical Customer Service.

AES Lite orders may be completed conveniently online!

- 1. Login to Sunrise Medical's order entry system
- 2. Start a new quote/order
- 3. Enter part number "EICV" and complete custom seating selections

AES Upholstery Options

AES upholstery options can be found on the AES and AES Lite product pages on www.SunriseMedical.com.



Whitmyer Headrests

Please refer to the Whitmver Product Catalog (MKT3000-3) for detailed information on head support options.



Looking for foot boxes, arm troughs, anterior trunk supports, or pelvic supports? Please refer to the JAY

JAY Positioning Supports

Positioning Supports Catalog (MK-100093) for detailed information.






SEATING SYSTEMS | PACKAGES


When ordering

- Specify seat insert and back insert dimensions
- **№** Specify upholstery fabric and color
- **Standard hardware diameter is 1":** specify if other size is needed
- Dimensions over 20" in width, height, or depth require additional modifications
- Curved Back modification requires change of hardware
- Standard seat foam is 1 ½" HR, standard back foam is 1" HR — see Modifications table for upgrades
- ¥ ½" multi-ply base material standard
- ¥ 1/4-20 T-nuts included



Universal Seating system package


- Includes solid seat insert, solid back insert, standard seat to back bracket, universal seat and back clamps with standard L brackets, and rear seat tabs.
- For each piece of hardware included in this package, please specify exact type/size needed.
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.



PKG10XX

Universal/E-Z Mount package


- Includes solid seat insert, solid back insert, standard seat to back bracket, 1" EZ Mount Seat Brackets, and universal back clamps with standard L brackets.
- For each piece of hardware included in this package, please specify exact type/size needed.
- Adjustment range: Up to 3" drop and up to 7" forward adjustment. NOTE: Seat to back bracket limits forward adjustment.



Space Saver package

- Includes solid seat insert, solid back insert, standard seat to back bracket, universal space saver seat and back clamps with standard L brackets, and rear seat tabs.
- For each piece of hardware included in this package, please specify exact type/size needed.
- Seat mounting hardware is 11/4" wide in this configuration.
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.


SEATING SYSTEMS | PACKAGES



PKG20XX

Adjustable J-Hook package


- Includes solid seat insert, solid back insert, standard seat to back bracket. Seat hardware consist of 2 hole J hooks and 5 hole L brackets. Back hardware consist of 6 hole I hooks and 2 hole L brackets. Also includes standard seat stops and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed and to what size tubing this system will be mounted ($\frac{3}{4}$ ", $\frac{7}{8}$ ", $\frac{1}{9}$ ").
- Between rail or full width back options available.



PKG23XX

Adjustable Flush Hook /J-Hook Seating package

- Includes solid seat insert, solid back insert, standard seat to back bracket, flush hooks seat hardware, 2 hole J hooks and 7 hole L bracket back hardware, seat stops and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed.
- Between rail or full width back options available.



E-Z Mount Seating package

- Includes solid seat insert, solid back insert, standard seat to back bracket, EZ mount brackets, and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed and to what size tubing this system will be mounted $(\frac{3}{4})''$, $\frac{7}{8}$ '', $\frac{1}{9}$ '').
- Seat hardware accommodates up to 1 1/4" drop adjustment.
- Back hardware accommodates up to 1 3/8" forward adjustment.

PKG90XX

Custom Seating Package

Consult with Sunrise Medical Customer Service for details on ordering this package

ı'L KIDDOS

SEATING SYSTEMS





Designed exclusively for wheelchair frames 14" and under, the Li'l Kiddos product line includes many innovations in pediatric and early intervention seating. Li'l Kiddos gives therapists and clinicians a pro-active approach to pediatric postural support. Maximum user weight 65 lbs.

When ordering:

- **अ** Specify seat insert and back insert dimensions up to 14" W
- **¥** Standard materials include: ¼" plywood base with T-nuts 1" pressure-relieving material in back and 1" pressure-relieving material in seat
- Adjustment rails are configured for 10-32 attaching hardware
- **№** Specify upholstery fabric and color
- **Y** Curved Back modification can not be applied to Li'l Kiddos systems: please note some package options include MOD5530 Curve Foam
- Additional Seat Modifications may be found on page 6
- Additional Back Modifications may be found on page 10



PKG40XX

Li'l Kiddos Seating System Package


- Includes solid seat insert, solid back insert, adjustable pediatric seat to back bracket, Li'l Kiddos universal fixed seat brackets, seat tabs, Li'l Kiddos universal adjustable back brackets, seat and back adjustment rails, and 1" push button pelvic positioning belt.
- For each piece of hardware included in this package, please specify exact type/size needed
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.



PKG41XX

Li'l Kiddos Contoured Seating System Package


- Includes contoured solid seat insert with 3" deep notches for back posts, solid I-back insert with curved foam, adjustable pediatric seat to back bracket, Li'l Kiddos universal fixed seat brackets, seat tabs, Li'l Kiddos universal adjustable back brackets, seat and back adjustment rails, and 1" push button pelvic positioning belt.
- For each piece of hardware included in this package, please specify exact type/size needed
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.



PKG42XX

Li'l Kiddos Independent Seat and Back Package

- Includes solid seat insert, solid back insert, Li'l Kiddos universal fixed seat brackets, Li'l Kiddos universal adjustable back brackets, seat and back adjustment rails, and 1" push button pelvic positioning belt.
- For each piece of hardware included in this package, please specify exact type/size
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.



PKG43XX

Li'l Kiddos Contoured Independent Seat and Back Package

- Includes contoured solid seat insert with 3" deep notches for back posts, solid I-back insert with curved foam, Li'l Kiddos universal fixed seat brackets, Li'l Kiddos Universdal adjustable back brackets, seat and back adjustment rails, 1" push button pelvic positioning belt.
- For each piece of hardware included in this package, please specify exact type/size needed
- Optional Seat/Back Modifications
- This seating system is transit approved. Please add part number MODWC20 for WC20 approval.

SEATS | INSERT PACKAGES

When ordering

- Specify seat insert dimensions; see measurement forms online or contact your Rehab Product Specialist
- Specify upholstery fabric and color
- Dimensions over 20" in width or depth require additional modifications
- Standard hardware diameter is 1", specify if other needed
- 1½" HR foam is standard. See Modifications table for upgrades
- ⅓" multi-ply base material standard
- ¥ ¼-20 T-nuts included



SSK10XX

Seat package with Universal/E-Z Mount hardware


- Includes solid seat insert, universal seat clamps with standard L brackets, and 1" EZ Mount brackets.
- For each piece of hardware included in this package, please specify exact type/size needed.



SSK17XX

Universal Seat Hardware package


- Includes solid seat insert, standard seat tabs, universal seat clamps with standard L brackets.
- For each piece of hardware included in this package, please specify exact type/size needed.
- Configured to be ordered with SBK17XX.



SSK20XX

Seat package with Adjustable I-hook and Rotating Hold Downs


- Includes solid seat insert, 6 hole J hooks, 2 hole L brackets, stops with fixed hold, and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed.
- Optional full width option available via 2 hole J hook and 5 hole L bracket.



SSK30XX

Seat package with E-Z Mount Hardware


- Includes solid seat insert and 1" EZ mount brackets.
- For each piece of hardware included in this package, please specify exact type/size needed.



SSK50XX

Seat package with Drop Hooks and Rotating Hold Downs

- Includes solid seat insert, 1 1/2" drop hooks, stops with fixed hold, and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed.
- 1 ½" drop is standard, specify if other is required.




SSK55XX

Seat package with Flush hooks and Rotating Hold Downs

- Includes solid seat insert, Flush hook hardware, stops with fixed hold, and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size needed.

LI'L KIDDOS SEAT PACKAGES






SSK70XX

Li'l Kiddos Solid Seat Insert Package


- Includes solid seat insert, Li'l Kiddos Universal fixed seat brackets, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size needed



SSK71XX

Li'l Kiddos Contoured Solid Seat Insert Package


- Includes contoured solid seat insert with 3" deep notches for back posts, Li'l Kiddos Universal fixed seat brackets, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size needed



SSK77XX

Li'l Kiddos Seat & Hardware Package

- Includes solid seat insert, Li'l Kiddos universal fixed seat brackets, seat tabs, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size
- Designed to be ordered with SBK77XX or SBK78XX back insert.



SSK78XX

Li'l Kiddos Contoured Seat & Hardware Package

- Includes contoured solid seat insert with 3" deep notches for back posts, Li'l Kiddos universal fixed seat brackets, seat tabs, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size needed
- Designed to be ordered with SBK77XX or SBK78XX back insert.

SEATS | SOLID INSERTS

When ordering

- Specify seat insert dimensions; see measurement forms at the back of this catalog
- Specify upholstery fabric and color
- Dimensions over 20" in width or depth require additional modifications. Upgrade to ¾" plywood base is recommended for oversize seats
- 1½" HR foam is standard see Seat Modifications table for premium foam upgrades
- **№** "multi-ply base material standard
- ¥ 1/4-20 T-nuts included
- When ordering as part of a package, please use MOD part number




- T-style seat insert
- Extender bars available, must be purchased separately.
- May also be ordered as MOD4300.



SSP60XX

Contour seat insert

- Four-way stretch fabric recommended.
- Outside lateral height over 4" requires MOD4640 Lateral Reinforcement.
- Specify dimensions on measurement form.
- May also be ordered as MOD4600.



Solid Seat Interface

SSI10XX

- Includes ¼-20 T-Nuts, but does NOT include Foam or upholstery.
- 1/4" Volara Foam.
- Widths greater than 20" require 3/4" plywood.



SSP10XX

Solid seat Insert



SSP20XX

Offset seat insert


- Specify long left side or long right side
- May also be ordered as MOD4200L (Left) or MOD4200R (Right)



SSP40XX

Anti-Thrust seat insert

- Specify height and depth of anti-thrust.
- May also be ordered as MOD4400.



SSP50XX

I-style seat insert

May also be ordered as MOD5000.



SSP70XX

Built-in Medial Thigh Support seat insert

- Four-way Stretch Fabric Face is recommended.
- Specify dimensions on measurement form.
- May also be ordered as MOD4700.



SSP80XX

Adjustable Split-seat insert

- Between Rail width measurement is required at time of ordering.
- Specify Left or Right adjustable side.
- Specify Split dimensions.
- Bilateral adjustment requires additional hinge hardware.
- May also be ordered as MOD4900.




Seat Cushion, Custom

SCI19XX

- Specify width and depth of cushion.
- Specify foam type and thickness.
- Specify upholstery fabric and color.
- Does not include plywood seat base or T-nuts
- Does not include hook & loop attaching strips.
- 2" HR Foam is standard.

SEATS | UNIVERSAL AMPUTEE INSERTS



AMPUTEE SEAT INSERTS

AMP10XX	Solid Insert, 16" W (actual width)
AMP11XX	Solid Insert, 18" W (actual width)
AMP16XX	Drop Base Style, for 16" wide W/C
AMP18XX	Drop Base Style, for 18" wide W/C
AMP19XX	Custom Insert (max: 22" W x 22" D)

- 1" HR foams is standard on seat insert and extension pad
- Extension pads are 8"W x 10"D
- Specify RIGHT or LEFT extension
- 16"D is standard, maximum insert size is 22"W x 22"D
- Adjustable angle brackets require notching the amputee insert

AMPUTEE INSERT COMPONENTS

AMP20XX	Additional Extension – Solid (removable)
AMP30XX	Additional Extension – Folding (non-removable)
AMP40XX	Additional Extension – Adjustable Angle (removable)
AMP70XX	Extension Pad Only
AMP80XX	Amputee Seat Insert Only

MODIFICATIONS

MOD6828 Add Folding Hinge to Extension Pad

MOD6826 Add Adjustable Angle Bracket to Extension Pad

MOD6827 No Seat Cushion


SEATS | WOOD BASES



WOOD SEAT BASES

	W	D	Н
WJB3100	14"	16"	1/4"
WJB3000	14"	16"	1/2"
WJB2100	16"	16"	1/4"
WJB2000	16"	16"	1/2"
WJB1100	18"	16"	1/4"
WJB1000	18"	16"	1/2"
WJB4000	20″	18"	1/4"
WJB4050	20″	18"	1/2"
WJB9000	Custom		

- Dimensions are FULL SIZE. Provide cushion specifications to ensure accurate sizing
- Does NOT include: T-nuts, foam, cover or hook & loop attaching strips



DROP-IN WOOD BASES

	W	D
WJB3500	14"	17"
WJB2506	16"	16"
WJB2500	16"	17"
WJB2508	16"	18"
WJB1506	18"	16"
WJB1500	18"	17"
WJB1508	18"	18"
WJB1206	20"	16"
WJB1200	20"	17"
WJB1208	20"	18"

- Bases are finished in Black
- Hardware not included
- MOD4125 Hook & Loop attaching strips included upon request
- 1/4-20 T-nuts are standard



DROP-IN WOOD BASE KITS

	W	Tube diameter	
WJB3551	14"	1"	
WJB3557	14"	7/8"	
WJB2551	16"	1"	
WJB2557	16"	7/8"	
WJB1551	18"	1"	
WJB1557	18"	7/8"	
WJB1251	20"	1"	
WJB1257	20"	7/8"	
WJB9100	Custom		

- Bases are finished in Black
- Includes mounting hardware: HDW0100 or HDW0078
- 1/4-20 T-nuts are standard
- All kit bases are 17" deep

SEATS | MODIFICATIONS & COVERS


FOAM MOD	DIFICATIONS	STRUCTUR	AL MODIFICATIONS
MOD0050	Omit Foam in Seat	MOD4150	Notches Back Post, 2" x 2"
MOD1180	1" Stimulite Honeycomb, Seat	MOD4153	Notches Back Post, 3" Deep
MOD1185	½" Stimulite Honeycomb, Seat	MOD4155	Notches For Front Hangers
MOD1195	Partial Stimulite Insert, Seat (max 100 sq in)	MOD4175	Mid Seat Cutout For X-Braces
MOD1200	2" HR Foam In Seat	MOD4180	Foam Form Seat
MOD1300	3" HR Foam In Seat	MOD4200L	Offset Seat, Long On Left Side
MOD1400	1" Pudgee Foam In Seat	MOD4200R	Offset Seat, Long On Right Side
MOD1405	½" Pudgee In Seat	MOD4300	T-style Seat
MOD1495	1" Pudgee Partial Insert, Seat (max 100 sq in)	MOD4400	Anti-Thrust Seat
MOD1500	5/8" Gel over 1" HR Foam In Seat	MOD4450	Quick Fit Anti-thrust Seat
MOD1525	½" Cubed Gel over 1" HR Foam In Seat	MOD4460	Wedge Seat
MOD1595	5/8" Gel Partial Insert, Seat (max 100 sq in)	MOD4500	Cubed Seat w/Cover (2x2x2 HR foam)
MOD1600	1" T-Foam In Seat	MOD4525	Custom Cubes in Seat
MOD1605	½" T-Foam In Seat	MOD4600	Contour Seat
MOD1615	1½" T-Foam In Seat	MOD4620	Sculptured Seat
MOD1700	2" T-Foam In Seat	MOD4630	Quick-Fit Contoured Seat
MOD1800	1" Sunmate Foam In Seat	MOD4635	Hip Relief, Sculptured Seat only
MOD1805	½" Sunmate Foam In Seat	MOD4640L	Lateral Reinforcement, Seat, Left
MOD1815	1½" Sunmate Foam In Seat	MOD4640R	Lateral Reinforcement, Seat, Right
MOD1900	2" Sunmate Foam In Seat	MOD4650	Replace Multi-Ply w/ABS, Seat
		MOD4655L	Belt Strap Notch, Seat, Left
UPHOLSTER	RY MODIFICATIONS	MOD4655R	Belt Strap Notch, Seat, Right
MOD4125	Hook & Loop strips, Loop sewn on bottom	MOD4700	Built-In Medial Thigh Support
MOD4125H	Hook & Loop strips, Hook sewn on bottom	MOD4725	Built-In Contour Lateral Adductor in Seat
MOD4126	Hook & Loop strips, Hook sewn on top	MOD4800	Seat Rail Foam Overlay
MOD5120	Add Hook/Loop Closure, Seat	MOD4825	Seat Front-edge Bevel Undercut
MOD5130	Add Faux Sheepskin to contact surface of seat	MOD4850L	Pelvic Obliquity Addition (up to 1"), Left
MOD5150	Incontinent Film Covering, Seat	MOD4850R	Pelvic Obliquity Addition (up to 1"), Right
MOD5160	Removable Surface Area Cover, Zipper at bottom	MOD4900	Adjustable Split-seat
		MOD4950	Ischial Shelf Addition
		MOD5000	I-style Seat
		MOD5700	Replace ½" Multi-ply with ¾", Seat
		MOD8000	Omit T Nuts, Seat
		MOD9005	Pre-Fit Modifications
		MOD9050	Oversize Seat, over 20", up to 24"
		MOD9075	Oversize Seat, over 24" up to 28
			·
		COVERS	
		COV4100	Additional Zippered Seat Cover, Standard, Flat
		COV4109	Additional Zippered Seat Cover, Custom
		COV4110	Additional Slip Over Seat Cover, Standard, Flat
		COV4119	Additional Slip Over Seat Cover, Custom
		COVOOVY	Additional Theorem I Come France

COV99XX Additional Zippered Cover, Foam Form

BACKS | INSERT PACKAGES

When ordering


- **अ** Specify seat insert dimensions; see measurement forms at the back of this
- Specify upholstery fabric and color
- 1" HR foam is standard, see Modifications table for upgrades
- Dimensions over 20" in width, height, or depth require additional modifications
- Standard hardware diameter is 1", specify if other is required
- 1/4"-20 T-nuts included
- ½" Multi-ply base material is standard



SBK20XX

Solid Back with Adjustable | Hook and Fixed/Rotating Hold Downs


- Includes solid back insert, 2 hole J hooks, 7 hole L brackets, stops with fixed hold, and stops with rotating holds.
- For each piece of hardware included in this package, please specify exact type/size needed
- Standard 7 hole L bracket allows for up to 1" drop and up to 1 1/2" forward adjustment and mounts forward of back uprights.
- Optional between uprights available (6 hole J hooks allow for drop range %" up to 3")



SBK10XX

Solid Back w/ Universal Hardware


- Includes solid back insert, universal back clamps with standard L brackets.
- For each piece of hardware included in this package, please specify exact type/size needed
- Up to 3" rear adjustment
- Up to 7" forward adjustment can be attained by reversing orientation of clamp brackets



SBK24XX

Solid Back with Mounting Hooks and Fixed Hold Downs


- Includes solid back insert, mounting hook with 1" lift, and stops with fixed hold.
- For each piece of hardware included in this package, please specify exact type/size needed
- Mounts forward of back uprights
- Standard 1" of lift



SBK12XX

Solid Back with Drop Hooks and Fixed Hold Downs

- Includes solid back insert, 1 1/2" drop hooks, and stops with fixed hold.
- For each piece of hardware included in this package, please specify exact type/size needed
- Mounts between back uprights
- Standard hardware has 1 1/2" drop. Specify if other is required.



SBK25XX

Solid Back with E-Z Mount Hdw and Rotating Hold Downs

- Includes solid back insert, EZ Mount with brackets, and stops with rotating hold.
- For each piece of hardware included in this package, please specify exact type/size
- 10 to 15 degree angle adjustment
- Up to 1 3/8" drop and forward adjustment

Need image

SBK90XX

Consult a Sunrise Customer Service Representative for guidance on ordering a custom AES back package.

LI'L KIDDOS

BACK PACKAGES






SBK70XX

Li'l Kiddos Back Insert Package

- Includes Li'l Kiddos solid back insert, Li'l Kiddos adjustable universal back brackets, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size



SBK71XX

Li'l Kiddos Contoured Back Insert Package

- Includes Li'l Kiddos solid I-back insert with curved foam, Li'l Kiddos adjustable universal back brackets, and adjustment
- For each piece of hardware included in this package, please specify exact type/size needed
- Optional back modifications available



SBK77XX

Li'l Kiddos Back & Hardware Package

- Includes Li'l Kiddos solid back insert, adjustable pediatric seat to back bracket, Li'l Kiddos adjustable universal back brackets, seat tabs, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size
- Designed to be ordered with SSK77XX or SSK78XX seat insert.



SBK78XX

Li'l Kiddos Contoured Back & Hardware Package

- Includes Li'l Kiddos solid I-back insert with curved foam, adjustable pediatric seat to back bracket, Li'l Kiddos fixed universal back brackets, sat tabs, and adjustment rails.
- For each piece of hardware included in this package, please specify exact type/size
- Designed to be ordered with SSK77XX or SSK78XX seat insert.

BACKS | SOLID INSERTS

When ordering


- **■** Specify back insert dimensions; refer to measurement forms online or contact
- Specify upholstery fabric and color
- Dimensions over 20" in width, or height require additional modifications
- 1" HR foam is standard, see Modifications table for foam upgrades
- ¥ 1/4-20 T-nuts included
- ½" multi-ply base material is standard
- Mounting hardware not included
- When ordering as part of a package, please use MOD part number



SBP30XX

T-style Back insert

Allows for enhanced placement of lateral thoracic and lateral pelvic supports



SBP10XX

Solid Back insert

Standard rectangular back made with 1/2" multi-ply base material and covered in 1" HR foam (standard)



SBP45XX

SBP46XX

Bi-angular Back insert


Adjusts at a pivot point to accommodate lower hip tissue



SBP20XX

I-style Back insert


Allows for enhanced placement of lateral thoracic supports



SBP50XX

Curved Back insert


Provides mild lateral pelvic and thoracic support, and offers improved comfort



SBP55XX

Contour Back insert

Provides moderate lateral pelvic support



SBP60XX


Adjustable Lumbar Back insert

Provides a customized placement and level of support of the lumbar spine

BACKS | MODIFICATIONS & COVERS


FOAM MOI	DIFICATIONS	STRUCTUR	AL MODIFICATIONS
MOD0055	Omit Foam In Back	MOD4645L	Lateral Reinforcement, Contour/Sculptured, Left
MOD2150	1 ½" HR foam in Back	MOD4645R	Lateral Reinforcement, Contour/Sculptured, Right
MOD2180	1" Stimulite Honeycomb, Back		E-Z Access Bracket Modification,
MOD2185	½" Stimulite Honeycomb, Back		Contour/Sculptured Back, Left
MOD2195	Stimulite Honeycomb Partial Insert, Back	MOD4646R	E-Z Access Bracket Modification,
MOD2200	2" HR Foam In Back		Contour/Sculptured Back, Right
MOD2300	3" HR Foam In Back	MOD5180	Foam Form Back
MOD2400	1" Pudgee Foam In Back	MOD5200	I Style Back
MOD2405	1/2" Pudgee Foam In Back	MOD5300	T Style Back
MOD2495	1" Pudgee Partial Insert (max 100 sq in)	MOD5425	Half Round Lumbar Roll
MOD2500	5/8" Gel Overlay in Back, (5/8" Gel over 1" HR)	MOD5426	Half Round Lumbar Roll w/Strap
MOD2595	%" Gel Partial Insert (max 100 sq in)	MOD5450	Adjustable Bi-Angular Back, w/2 Adj Brackets
MOD2600	1" T-Foam In Back	MOD5451	Adjustable Bi-Angular Back, w/1 Adj Bracket
MOD2605	1/2" T-Foam In Back	MOD5460	Adjustable Lumbar Back w/ 5.5" Height Adjustment Range
MOD2615	1 ½" T-Foam In Back	MOD5461	Adjustable Lumbar Back w/ 3" Height Adjustment Range
MOD2650	1/2" T-Foam In Thoracic Support	MOD5462	Adjustable Lumbar Back w/ 7.5" Height Adjustment Range
MOD2700	2" T-Foam In Back	MOD5500	Curved Back
MOD2800	1" SunMate Foam In Back	MOD5501	Curved Back, w/2 ea Hinged L Brackets
MOD2805	1/2" SunMate Foam In Back	MOD5502	Curved Back, w/4 ea Hinged L Brackets
MOD2815	1 ½" SunMate Foam In Back	MOD5550	Cubed Back w/Cover (2x2x2 HR cubes)
MOD2900	2" SunMate Foam In Back	MOD5575	Custom Cubes in Back
MOD5525	Back Rail Foam Overlay	MOD5650	Substitute ABS for Multi-Ply In Back
MOD5530	Curve Foam Addition To Flat Back	MOD5705	Replace 1/2" Multi-Ply with 3/4", Back
		MOD5800	Contour Back Modification
UPHOLSTE	RY MODIFICATIONS	MOD5801	Quick Fit Contour/Bi-Angular Back Modification
MOD5125	Add Hook & Loop closure to Back upholstery	MOD5802	Quick Fit Contour Back Only
MOD5135	Add Faux Sheepskin to contact surface of back	MOD5830	Arm Relief-Back, Contour/Sculptured
MOD5165	Removable Top Cover, zipper at center	MOD5835	Harness Strap Notches, Back
MOD5900	Monogram Back w/Name, Block Letter (single line)	MOD5850	Scapular Contour Cutout, Back
MOD5901	Monogram Back w/Name, Script Letter (single line)	MOD5855	Top Edge Angle Cut, back
MOD5902 [†]	Monogram Back, Stock Design	MOD8100	Omit T Nuts, Back
		MOD9105	Pre-Fit Modifications, Back
		MOD9150	Oversize Back; over 20" and up to 24"
†MOD5902 Sto	ck Designs may be viewed online at www.aesys.com	MOD9175	Oversize Back; over 24" and up to 28"
		COVERS	
		COV5100	Additional Zippered Back Cover, Standard, Flat
		COV5109	Additional Zippered Back Cover, Custom
		COV5110	Additional Slip Over Back Cover, Standard, Flat
		COV5119	Additional Slip Over Back Cover, Custom
		COV99XX	Additional Zippered Cover, Foam Form

FOAM | TYPES




PUDGEE FOAM

Extra soft foam with a gel-like consistency conforms to the client and wicks moisture away from the body.




STIMULITE HONEYCOMB

Honeycomb material provides pressure relief, shear reduction, and ventilation, and is machine-washable.




T-FOAM

Memory foam technology that was originally developed by NASA. Designed to evenly distribute pressure and minimize heat build-up.



SUNMATE (VISCO) FOAM

Viscoelastic "memory" foam conforms to the client to distribute weight evenly for pressure reduction and maximum comfort.




GEL

Available in either cubed or solid form, gel inserts provide even pressure relief for the client across the entire surface.



JAY FLUID

Tripad conforms to each individual's shape, adjusts to sitting positions, and ensures proper fluid placement beneath bony prominences to help protect the skin from breakdown.



HIGH-RESILIENCY (HR) FOAM

Durable, high quality foam encourages proper positioning.

CUSTOM MOLDED | FOAM FORM





SSP85XX

Foam Form Solid Seat Insert

- Four-way stretch fabric required on face
- 1 unit liquid Soft SunMate in bottles (FFK1010) is standard, specify if medium density is desired, or if unit is preferred in bag form.
- Base foam layer standard is 1" Soft SunMate
- Standard cover is 4" high and includes dualpull zipper along back edge and both sides; 1/4" soft SunMate sewn into face
- Plastic bag included for molding

SBP80XX

Foam Form Solid Back Insert

- Four-way stretch fabric required on face
- 1 unit liquid Soft SunMate in bottles (FFK1010) is standard, specify if medium density is desired, or if unit is preferred in bag form
- Base foam layer standard is 1" soft SunMate
- Standard cover is 4" deep and includes dualpull zipper on top and both sides, and 1/4" soft SunMate sewn into face
- Lateral supports are additional; refer to FTK part numbers for options
- Plastic bag included for molding

FOAM FORM THORACIC SUPPORT KIT

- Each Kit contains two (2) thoracic supports and four (4) 7-hole L brackets
- Soft density liquid SunMate is standard
- Helps direct foam during pour
- Reinforces lateral section



5H X 5D	7H X 5D	9H X 6D	11H X7D	CUSTOM
FTK10XX	FTK11XX	FTK12XX	FTK13XX	FTK19XX




FOAM FORM SUPPLIES


UNIT	BOTTLE	BAG
1/2	FFK1005	FFK1005B
1	FFK1010	FFK1010B
1½	FFK1015	FFK1015B

- Soft Density is Standard
- Medium Density is Standard

HARDWARE | SEAT TO BACK BRACKETS




- Flush mount hardware for tall back heights.
- Standard height flush or offset mount bracket.
 Wide is double the width for extra stability if back, low profile can be ordered in flush or bracket will be bearing significant weight or repetitive pressure (client tone).
 - back, low profile can be ordered in flush or offset options. Wide is double the width for extra stability if bracket with be bearing significant weight or repetitive pressure (client tone).





Allows for up to 2 3/4" of forward adjustment depending on option ordered. Ideal for clients who may have changing seat depth needs over time.



HARDWARE | UNIVERSAL SEATING MOUNTING BRACKETS



- Adjusts 3" above rail
- Adjusts 1" below rail
- Saddle cutout to mount on double rails
- Up to 3" height adjustment
- No height adjustment
- Limited height adjustment

HARDWARE | UNIVERSAL BACK MOUNTING **BRACKETS**

	UNIVERSAL	UNIVERSAL W/ HINGED L	UNIVERSAL, FIXED
LEFT	HDW6820L	HDW6820HL	HDW6819BL
RIGHT	HDW6820R	HDW6820HR	HDW6819BR

- Adjusts up to 71/4" forward of rails
- Adjusts up to 3" behind rails
- For use on Curved Wood backs
- No depth adjustment

HARDWARE | UNIVERSAL COMPONENTS

	STANDARD	TRANSPORT OPTION	SHORT
SEAT L BRACKETS	21/4" 11/4	2.1/2	
LEFT	HDW6813	HDW6863	HDW6815
RIGHT			

	STANDARD	TRANSPORT OPTION	SHORT
BACK L BRACKET	15%"	2-1/2" 2-1/2"	
LEFT	HDW6823L	HDW6868L	HDW6868L
RIGHT	HDW6823R	HDW6868R	HDW6868R
PAIR	HDW6823	HDW6868	HDW6868

TS	HINGED
BACK LBRACKETS	
LEFT	HDW6867L
RIGHT	HDW6824R
PAIR	HDW6824


Z	100°	100° TRANSIT	107°	110°
BACK LBRACKETS				
LEFT	HDW6803L	HDW6867L	HDW6804L	HDW6805L
RIGHT	HDW6803R	HDW6867R	HDW6804R	HDW6805R
PAIR	HDW6803	HDW6867	HDW6804	HDW6805

For use with Curved Wood Backs only.



HARDWARE | E-Z MOUNT BRACKETS

S	E-Z MOUNT	E-Z MOUNT, SHORT	E-Z MOUNT FIXED
E-Z MOUNT BRACKETS	00		© C
1" TUBE	HDW6850	HDW6850P	HDW6853
7/8" TUBE	HDW6851	HDW6851P	_
3/4" TUBE	HDW6852	HDW6852P	_
CLAMP ONLY	HDW6850C	HDW6851C	HDW6852C



HARDWARE | BI-ANGULAR BACK BRACKETS


NO TEETH	WITH TEETH	LONG
23/4" 2"	23/4" 2"	41/4"
HDW6826T	HDW6826	HDW6826L

HARDWARE | J & L BRACKETS

	FLUSH HOOK SEAT HDW KIT	DROP HARDWARE KIT	FULL WIDTH HDW KIT
1" TUBE	HDW0550	HDW0100	HDW0110
%" TUBE	HDW0560	HDW0078	HDW0178

- All kits listed include fasteners
- All kits listed include fasteners
- All kits listed include fasteners

S	TWO HOLE	THREE HOLE	FIVE HOLE	SEVEN HOLE	TEN HOLE
ADJUSTABLE L BRACKETS	11/2"		3"	4"	51/2"
STANDARD	HDW6120	HDW6130	HDW6150	HDW6170	HDW6110
110 DEGREE	HDW6125	HDW6135	HDW6155	HDW6175	HDW6115
ADJUSTMENT	_	_	½" forward	up to 1.5" forward	up to 3" forward







SEAT TABS		
31/2"	24"	
SHORT	HDW6301	
LONG	HDW6303	

HARDWARE | J & L COMPONENTS







	DROP STOP	DROP STOP W/ HOLD DOWN	DROP STOP W/ ROTATING HD
RAIL MOUNT	HDW6406	HDW6405	_
¾" TUBE	HDW6410	HDW6440	HDW6445
%" TUBE	HDW6420	HDW6450	HDW6455
1" TUBE	HDW6430	HDW6460	HDW6465
11/8" TUBE		HDW6461	



and 1" round tube diameters

LI'L KIDDOS **HARDWARE**

ATTACHING



	ADJUSTABLE	FIXED	SEAT L BRACKET
LI'L KIDDOS SEAT BRACKETS			
LEFT	HDW7810L	HDW7819SL	HDW7861L
RIGHT	HDW7810R	HDW7819SR	HDW7861R

- Hardware is slotted for 10-32 and requires adjustment rails on seats and backs
- Up to 2" height adjustment range and up to 2" drop adjustment range
- Hardware is slotted for 10-32 and requires adjustment rails on seats and backs
- Hardware is slotted for 10-32 and requires adjustment rails on seats and backs

	ADJUSTABLE	FIXED	BACK L BRACKET
LI'L KIDDOS BACK BRACKETS			
LEFT	HDW7820L	HDW7819BL	HDW7866L
RIGHT	HDW7820R	HDW7819BR	HDW7866R

- Hardware is slotted for 10-32 and requires adjustment rails on seats and backs
- Up to 2" height adjustment range and up to 2" drop adjustment range
- adjustment rails on seats and backs
- Hardware is slotted for 10-32 and requires Hardware is slotted for 10-32 and requires adjustment rails on seats and backs

SEAT TRANSPORT OPTION	BACK TRANSPORT OPTION
1-1/4*	1.1/4
HDW7862	HDW7867L
	HDW7867R

LI'L KIDDOS E-Z MOUNT	3½" 156"
<i>7</i> %"	HDW7851
1"	HDW7850
	- 11

- Hardware is slotted for 1/4-20
- Li'l Kiddos clamp is 1" wide for smaller, space-challenged frames

LATERAL SUPPORT | LATERAL THORACIC **PACKAGES**

When ordering

- **№** Specify pad size up to 5H x 6D, larger sizes may incur additional charges
- **¥** Specify if offset is needed
- Adjustment rail included; 6" is standard
- **№** 1/2" Medium pressure-relieving foam is standard, please specify if other is
- **№** Specify upholstery fabric and color
- Specify base material; upcharges may apply
- ¥ 1/4-20 T-nuts included
- Sold as EACH; order two (2) for bilateral application





LTK24XX

Lateral Thoracic Package w/ Medium or High Profile Bracket, Flat Pad

- Standard package includes flush, medium or high profile bracket and flat pad.
- Offset brackets in 1" or 2" available.
- For curved backs, AES recommends package LTK34XX instead.



LTK30XX

Lateral Thoracic Package w/ Front Quick Adjust and E-Z Access Bracket, Flat Pad

- Includes Front Quick Adjust and E-Z Access bracket and flat pad.
- E-Z Access swing-away hinge pivot point can be rotated up to 130°.
- For backs over 2" thick, height adapter is required.
- Adjusts from ½" offset to 1" or 1 ½" offset.
- Package includes modular hardware components are interchangeable



LTK22XX

Lateral Thoracic Package w/ Standard Bracket, Flat Pad

- Standard package includes flush, standard profile bracket and flat pad.
- Offset brackets in 1" or 2" available.
- For curved backs, AES recommends package LTK34XX instead.



LTK25XX

Lateral Thoracic Package w/ E-Z Access Bracket, Flat Pad

- Includes E-Z Access bracket and flat pad.
- E-Z Access swing-away hinge pivot point can be rotated up to 130°.
- For backs over 2" thick, height adapter is reauired.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



LTK34XX

Lateral Thoracic Package w/ Fixed Modular Bracket, Flat Pad

- Includes Fixed Modular bracket and flat pad.
- For backs over 2" thick, height adapter is
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



LTK23XX

Lateral Thoracic Package w/ Standard Bracket, Curved Pad

- Standard package includes flush, standard profile bracket and curved pad.
- Offset brackets in 1" or 2" available.
- For curved backs, AES recommends package LTK34XX instead.



m

LTK27XX

mb

Lateral Thoracic Package w/ E-Z Access Bracket, Curved Pad

- Includes E-Z Access bracket and curved pad.
- E-Z Access swing-away hinge pivot point can be rotated up to 130°.
- For backs over 2" thick, height adapter is reauired.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable



LTK35XX

Lateral Thoracic Package w/ Rear Quick Adjust Bracket, Flat Pad

- Includes Rear Quick Adjust bracket and flat pad.
- Rear Quick Adjust has 2" adjustment range: overall length expands from 3 $\frac{3}{4}$ " to 5 $\frac{3}{4}$ ".
- For backs over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.

LATERAL SUPPORT | LATERAL THORACIC **PACKAGES**





LTK37XX

Lateral Thoracic Package w/ Rear Quick Adjust and E-Z Access Bracket, Flat Pad

- Includes Rear Quick Adjust and E-Z Access bracket and flat pad.
- Rear Quick Adjust has 2" width adjustment range: overall length expands from 33/4" to
- E-Z Access swing-away hinge pivot point can be rotated up to 130°.
- For backs over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



LTK40XX

Lateral Thoracic Package w/ 1" Reverse Offset Bracket, Flat Pad

- Includes 1" Reverse offset bracket and flat pad.
- Reverse offset Profile is 11/2".
- For curved backs or backs over 1½" thick, AES recommends package LTK34XX instead".





LTK45XX

Lateral Thoracic Package w/ Removable Bracket, Flat Pad

- Includes Removable bracket and flat pad.
- For backs over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



Lateral Thoracic Package w/ Rear Quick Adjust and Removable Bracket, Flat Pad

- Includes Rear quick adjust and removable bracket and flat pad.
- Rear Quick Adjust has 2" width adjustment range: overall length expands from 33/4" to 53/4".
- For backs over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.


LATERAL SUPPORT | LATERAL PELVIC/ ARMREST PACKAGES



LPK40XX

Lateral Pelvic/Thigh/Armrest Combo w/ Standard Fixed Bracket


- Includes pelvic pad, armrest pad and fixed brackets (armrest and seat).
- Specify both pelvic pad and armrest pad dimensions.
- Minimum pelvic pad height is 3 /4".



LPK41XX

Lateral Pelvic/Thigh/Armrest Combo w/ Removable Bracket

- Includes pelvic pad, armrest pad, fixed armrest bracket and removable lateral bracket.
- Specify both pelvic pad and armrest pad dimensions.
- Minimum pelvic pad height is 33/4".
- For seats over 2" thick, height adapter is reauired.



LPK42XX

Lateral Pelvic/Thigh/Armrest Combo w/ Removable and Height Adjustable Bracket

- Includes pelvic pad, armrest pad, height adjustable armrest bracket and removable lateral bracket.
- Specify both pelvic pad and armrest pad dimensions.
- Minimum pelvic pad height is 3¾".
- For seats over 2" thick, height adapter is required.

LATERAL SUPPORT | LATERAL PELVIC/ THIGH PACKAGES

When ordering

- Specify pad size up to 5H x 18L, larger sizes may incur additional charges
- Specify upholstery fabric and color
- Standard Foam is ½" premium pressurerelieving material, specify if other is desired
- Specify base material; upcharges may apply
- 14-20 T-nuts included
- Sold as EACH; order two (2) for bilateral application



LPK22XX

Lateral Pelvic/Thigh Package w/ Standard Bracket, Pad

- Includes standard bracket and pad.
- Standard fixed bracket available in 1", 2" or 3" offset.
- Offset profile is 1 5/8" high.



LPK27XX

Lateral Pelvic/Thigh Package w/ Removable Bracket, Pad

- Includes removable bracket and pad.
- Offset adapters available in 1" or 2".
- Height adapters available.
- Package includes modular hardware components are interchangeable.



LPK32XX

Lateral Pelvic/Thigh Package w/ Med/ High Profile Bracket, Pad

- Includes Medium or High Profile Bracket and pad.
- Medium profile bracket available in flush, 1, 2, or 3" offset.
- Medium profile bracket offset is 2 1/4" high.
- High profile bracket available in flush, 1" or 2" offset.
- High profile bracket offset is 2 3/4" high.



LPK33XX

Lateral Pelvic/Thigh Package w/ 20° Adjustable and Removable Bracket, Pad

- Includes 20° Adjustable and Removable Bracket and pad.
- Up to 20° pivot adjustment.
- For seats over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



LPK34XX

Lateral Pelvic/Thigh Package w/ Modular Bracket

- Includes Modular Bracket and pad.
- For seats over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.



LPK55XX

Lateral Pelvic/Thigh Package w/ 1" Reverse Offset Bracket, Pad

- Includes 1" Reverse offset bracket and pad.
- Offset profile is 1 5/8" high.



LPK57XX

20° Adjustable, Fixed Bracket

- and pad.
- Up to 20° pivot adjustment.
- For seats over 2" thick, height adapter is required.
- Offset adapters available in 1" or 2".
- Package includes modular hardware components are interchangeable.





LI'L KIDDOS | LATERAL THORACIC SUPPORT PACKAGES



When ordering

- NOTE: ADJUSTMENT RAILS NOT **INCLUDED** with Lateral Thoracic packages. Some seat packages and some back packages come standard with Adjustment rails, please check before ordering rails separately
- Specify pad size up to 5H x 6D, larger sizes may incur additional charges
- **≥** 1/4" base material and 10-32 T-nuts
- Specify if offset is needed
- 1/2" Medium Sunmate foam is standard, please specify if other is desired
- **■** Specify upholstery fabric and color
- Sold as EACH; order two (2) for bilateral application



LTK72XX

Li'l Kiddos Lateral Thoracic Support Package

• Specify Flush, 1/2" offset or 1" offset bracket



LTK73XX

Li'l Kiddos Lateral Thoracic Support Package w/ curved pad

• Specify Flush, 1/2" offset or 1" offset bracket



LTK74XX

Li'l Kiddos Lateral Thoracic Support Package w/ Fixed Modular Bracket, flat pad

- Specify flush, 1/2" offset, 1" offset, 1 1/2" offset, or 2" offset
- 1" Height adapter included; 2" height adapter available



LTK75XX

Li'l Kiddos Lateral Thoracic Support Package w/E-Z Access Bracket

- Specify flush, ½" offset, 1" offset, 1½" offset, or 2" offset
- 1" Height adapter included; 2" height adapter available



LTK76XX

Li'l Kiddos Lateral Thoracic Support Package, w/Rear Quick Adjust Bracket

- Specify flush, 1/2" offset, 1" offset, 1 $\frac{1}{2}$ " offset, or 2" offset
- 1" Height adapter included; 2" height adapter available



LTK77XX

Li'l Kiddos Lateral Thoracic Support Package w/Rear Quick Adjust and E-Z Access Bracket

- Specify flush, 1/2" offset, 1" offset, 1½" offset, or 2" offset
- 1" Height adapter included; 2" height adapter available



LTK78XX

Li'l Kiddos Lateral Thoracic Support Package w/E-Z Access Bracket, Curved

- Specify flush, 1/2" offset, 1" offset, 11/2" offset, or 2" offset
- 1" Height adapter included; 2" height adapter available

LI'L KIDDOS | LATERAL PELVIC/THIGH SUPPORT PACKAGES



When ordering

- NOTE: ADJUSTMENT RAILS NOT **INCLUDED** with Lateral Pelvic/Thigh packages. Some seat packages and some back packages come standard with Adjustment rails, please check before ordering rails separately
- Specify pad size up to 5H x 14L, larger sizes may incur additional charges
- ¥ 1/4" base material and 10-32 T-nuts
- **अ** Specify if offset is needed
- [№] ½" Medium Sunmate foam is standard, please specify if other is desired
- **■** Specify upholstery fabric and color
- Sold as EACH; order two (2) for bilateral application



LPK72XX

Li'l Kiddos Lateral Pelvic/Thigh Support Package, w/Standard Bracket

- Specify Flush, 1/2" offset or 1" offset
- For offsets greater than 1", switch out to Modular bracket is required at an additional



LPK74XX

Li'l Kiddos Lateral Pelvic/Thigh Support Package, w/Modular Fixed Bracket



LPK75XX

Li'l Kiddos Lateral Pelvic/Thigh Support Package, w/ 20° Adjustable Bracket

- Specify 1" offset or 2" offset
- 1" Height adapter included; 2" height adapter is available



LPK76XX

Li'l Kiddos Lateral Pelvic/Thigh Support Package w/High Profile Bracket

- Specify Flush, 1/2" offset or 1" offset
- For offsets greater than 1", switch out to Modular bracket is required at an additional charge



LPK77XX

Li'l Kiddos Lateral Pelvic/Thigh Support Package, w/Removable Bracket

- Specify flush, ½" offset, 1" offset, 1½" offset, or 2" offset
- 1" Height adapter included; 2" height adapter is available



Package, w/Removable, 20^o Adjustable Bracket

- Specify 1" offset or 2" offset
- 1" Height adapter included; 2" height adapter is available



LATERAL SUPPORT | THORACIC PADS

	D-FORWARD	TAPERED	CURVED
STANDARD	LTS25XX	LTS15XX	LTS35XX
CUSTOM	LTS29XX	LTS19XX	LTS39XX

RECOMMENDED THORACIC PAD SIZES

Н	D
3	4
4	5
4	6
5	6

- The standard and most commonly used thoracic lateral support pad that is in the shape of a "D."
- The tapered pad shape helps prevent irritation from client arm movement across the midline.
- The curved pad provides mild anterior support in addition to lateral support as it was designed to replicate the shape of a hand wrapping around and supporting a client's trunk.
- Specify pad size; see standard Specify pad size; see size chart.
 - standard size chart.
- Specify pad size; see standard size chart.

CONTOUR PAD PROFILE (TOP VIEW)		
$ \begin{array}{c} $	\overrightarrow{MD} $AD \longrightarrow$	
FLAT	CURVED	


RECOMMENDED THORACIC **CONTOURED PAD SIZES**

Н	AD	MD
5	5	2 ½
7	5	3
9	6	3
11	7	3 ½

AD=Anterior Depth MD=Medial Depth

	CONTOURED
STANDARD	LTS45XX
CUSTOM	LTS49XX


Specify both pad profile and pad shape (at right) for contoured pads.



STANDARD

CUSTOM

LATERAL SUPPOR VIC/THIGH PADS



- The D-Up lateral pelvic/thigh pad is in the shape of a "D" pointed up. The D-Up pad provides a soft contact surface for the client while facilitating proper lateral pelvic positioning and helping minimize unwanted weight shifts.
- Specify pad size; see recommended size chart.
- For full length, specify height and length.
- D-up is standard—specify if D-forward is preferred.



4 6

5

6

5 12

5 14

- The Adductor lateral pelvic/thigh pad mounts to the back and is in the shape of a "D." The adductor pad provides a soft contact surface for the client's hips while facilitating proper lateral pelvic positioning and helping minimize unwanted weight shifts.
- Specify pad size; see recommended size chart.
- D-forward is standard—specify if D-up is preferred.

ADD24XX

ADD29XX

LI'L KIDDOS I THORACIC PADS



	D-FORWARD	TAPERED	CURVED
STANDARD	LTS25XXLK	LTS15XXLK	LTS35XXLK
CUSTOM	LTS29XX	LTS19XX	LTS39XX


- Specify pad size; see standard Specify pad size; standard size chart at right
 - sizes are listed in chart at right
- Specify pad size; standard sizes are listed in chart at right


RECOMMENDED THORACIC PAD SIZES

Н	D
3	4
4	5
4	6
5	6

CONTOURED STANDARD LTS45XXLK LTS49XXLK **CUSTOM**

• Specify both pad profile and pad shape (at right) for contoured pads





RECOMMENDED THORACIC **CONTOURED PAD SIZES**

Н	AD	MD
5	5	2 1/2
7	5	3
9	6	3
11	7	3 ½

AD=Anterior Depth MD=Medial Depth

LI'L KIDDOS

PELVIC/THIGH PADS



8 10

12


14

	D-UP	
STANDARD	LPS24XXLK	
CUSTOM	LPS29XXLK	

- Standard is D-up—specify if D-forward is preferred
- For full-length pads, specify height and length

STANDARD PELVIC/THIGH **PAD SIZES**

Н	L
3	4
4	5
4	6
5	6



RECOMMENDED ADDUCTOR PAD SIZES

Н	L	Н
3	4	4
4	5	4
4	6	5
5	6	5

- Standard is D-forward—specify if D-up is preferred
- For full-length pads, specify height and length

LATERAL SUPPORT | MODIFICATIONS

THORACIC PAD FOAM MODIFICATIONS

MOD2550	Replace foam with ¼" Gel in Thoracic Support
---------	--

MOD2650	½" T-Foam in Thoracic Support
MOD2750	½" Pudgee in Thoracic Support
MOD2850	½" Sunmate in Thoracic Support
MOD2861	1" Sunmate in Thoracic Support
MOD3100	Waterfall Foam Thoracic Sppt

THORACIC PAD STRUCTURAL MODIFICATIONS

MOD2950 T	aper Front edge of	Lateral Support
-----------	--------------------	------------------------

Replace Multi-ply with 1/4" ABS MOD3000 MOD3001 Flat metal base, up to 5x6" MOD3011 Curved metal base, up to 5x6"

MOD3150 Contour-fit Pad MOD3175 Hinged Pad

MOD5510 **Curved Lateral Pad** MOD6803† Curved Back Bracket

† Applies to Modular Hardware only

Modifications also available on Li'l Kiddos Seating

PELVIC/THIGH PAD FOAM MODIFICATIONS

MOD3550	Replace foam with $\frac{1}{4}$ " Gel in Pelvic Support
MOD3650	½" T-foam in Pelvic Support
MOD3750	½" Pudgee in Pelvic Support
MOD3850	½" Sunmate in Pelvic Support
MOD3861	1" Sunmate in Pelvic Support
MOD3500	Waterfall Foam Pelvic Sppt

PELVIC/THIGH PAD STRUCTURAL MODS

MOD3025	D forward pad shape
---------	---------------------

MOD3050 Replace multi-ply with 1/4" ABS

MOD3055 Contoured Pelvic support for Anti-thrust seat

LATERAL SUPPORT | BRACKET COVERS



LI'L KIDDOS | LATERAL BRACKET COVERS





		COVERTENGIII	ACCOMINIODATES
SHORT	COV2650	4 1/2"	Flush brackets
LONG	COV2651	5 ½"	Offset or Height adapter


LATERAL SUPPORT | FIXED BRACKETS

	STANDARD PROFILE	MEDIUM PROFILE	HIGH PROFILE	REVERSE OFFSET
	1 1/3	5° 2'/4" 1	6 ³ / ₄ " 2 ³ / ₄ " 3 ³ / ₄ "	27
FLUSH	HDW2000	HDW2001	HDW2002	_
1" OFFSET	HDW2100	HDW2103	HDW2102	HDW2410
2" OFFSET	HDW2200	HDW2203	HDW2202	_
3" OFFSET	HDW2300	HDW2303	_	_

- Suitable for either thoracic or Suitable for either thoracic or Suitable for either thoracic pelvic applications
 - pelvic applications
- pelvic applications
- or pelvic applications
- Not available for Seats with X-brace cut outs

E-Z ACCESS swing-away

LATERAL SUPPORT | ADJUSTABLE BRACKET AND MODULAR BRACKETS



2" OFFSET	HDW2205	
	•	Suitable for either thoracic or pelvic applications

1" OFFSET

• Profile adjusts from 2 3/4" up to 4" high

	5½*	2" 51/2"	
FLUSH	HDW2540	HDW2500	
FLUSH: pediatric	HDW2540P	HDW2500P	
1" OFFSET	HDW2640	HDW2600	
1" OFFSET: pediatric	HDW2640P	HDW2600P	
2" OFFSET	HDW2740	HDW2700	
2" OFFSET: pediatric	HDW2740P	HDW2700P	

• Suitable for either thoracic or • Suitable for either thoracic or pelvic applications

FIXED MODULAR

- Pediatric (P) brackets are configured with both short pad mount plates and short mounting plates
- pelvic applications

LATERAL SUPPORT | MODULAR BRACKETS




HOW MODULAR HARDWARE WORKS

Each bracket is built of interlocking components that can be mixed and matched to achieve perfect positioning of secondary supports.

Modular hardware provides an economical and easy way to make changes to a client's system based upon their growth and evolving posture support needs.

For example: the bracket at right can be modified by replacing (D) PRT2502 Mounting Plate with HDW2501 for Rear Quick Adjust capabilities.

Replace (A) PRT Offset Pad Mount Plate and (B) Offset Adapter with PRT2625 Flush Pad Mount Plate for a flush lateral bracket.



	REAR QUICK ADJUST (RQA)	E-Z ACCESS swing away + RQA	E-Z ACCESS + FRONT Quick Adjust
	Up to 2" adjustment in ½" increments	Up to 2" adjustment in ½" increments 130° rotation	31/4"
FLUSH	HDW2513	HDW2511	_
FLUSH: short		HDW2511P	_
1" OFFSET	HDW2613	HDW2611	HDW2655
1" OFFSET: short		HDW2611P	_
2" OFFSET	HDW2713	HDW2711	HDW2755
2" OFFSET: short		HDW2711P	_

- Recommended for lateral thoracic applications only
- Adapter required for Curved Back application
- Recommended for lateral thoracic applications only
- Adapter required for Curved Back application
- Recommended for lateral thoracic applications only
- $\bullet~$ HDW 2655 adjusts to $\ensuremath{\textit{12}}\xspace"$ and 1"
- HDW2755 adjusts to 1/2", 1", and 11/2"
- Adapter required for Curved Back application


	20º ADJUSTABLE FIXED	E-Z ACCESS REMOVABLE	20º ADJ + REMOVABLE
	5 3/8" 21/4"	5½"	53/8" 27/4"
FLUSH	_	HDW2550	_
FLUSH: short	_	HDW2550P	_
1" OFFSET	HDW2659	HDW2551	HDW2650
1" OFFSET: short	HDW2659P	HDW2551P	HDW2650P
2" OFFSET	HDW2759	HDW2552	HDW2750
2" OFFSET: short	HDW2759P	HDW2552P	HDW2750P

- Suitable for either thoracic or pelvic applications
- Adapter required for Curved Back application
- Suitable for either thoracic or pelvic applications
- Adapter required for Curved Back application
- Suitable for either thoracic or pelvic applications
- Adapter required for Curved Back application

LATERAL SUPPORT | MODULAR BRACKET COMPONENTS

FLUSH OFFSET + 20° ADJUSTABLE 1/2" OFFSET ONLY **OFFSET** PAD MOUNT PLATES STD PRT2625 PRT2627 PRT2653 PRT2630 **SHORT** HDW2503 HDW2509 PRT2654


- Add height adapter for increased profile
- For offset bracket, use Offset Pad Mount Plate and Offset Adapter
- Thoracic or Pelvic application
- Requires offset adapter (shown below)
- Thoracic or Pelvic application
- Requires 20° adjustable offset adapter (shown below
- Thoracic or Pelvic application
- Use when only 1/2" offset is desirable
- Can not be configured with offset adapters or 20° adjustable offset adapters
- Thoracic or Pelvic application





- Use with modular brackets to accommodate thicker seats and backs
- Thoracic or Pelvic application
- Thoracic or Pelvic application

	BRACKET EXTENSION
SHORT	HDW2524
LONG	HDW2525
	• Connects had mount plate to

- Connects pad mount plate to mounting plate.
- Use with fixed modular bracket, RQA only bracket, 20° adjustable bracket







• For use only with other modular hardware components and packages — can not be used with standard fixed brackets



- Mounts to W/C upright canes, seat rail, or leg rest
- Fits 3/4, 7/8, and 1" tubing



LATERAL SUPPORT | MODULAR BRACKET COMPONENTS


REMOVABLE BRACKET BASE (female assembly)	REMOVABLE BRACKET, male assembly	E-Z ACCESS HINGE
11/4"	1" 2"	2" 2"/4"
PRT2601	PRT2650	PRT2500

• Hinge rotates 130 degrees

REAR	REAR QUICK ADJUST	
STD	HDW2501	
95 degree	HDW2535	
110 degree HDW2502 107 degree HDW2537		
		100 degree HDW2538

ADJUSTMENT RAIL				
	2"	HDW7102	10"	HDW7110
	3"	HDW7103	12"	HDW7112
16	4"	HDW7104	14"	HDW7114
160	6"	HDW7106	16"	HDW7116
	7"	HDW7107	18"	HDW7118
	8"	HDW7108		

- Can be rotated 180 degrees to move slot to medial
- Can be mounted in dual or staggered configurations





PRT0015

Modular Hardware Kit

Case filled with the most commonly used Modular Hardware components

1110 01011011	. real erroll of our porterior		
HDW2025	1/4" Bracket Spacer	HDW2528	3" Distal Extension Adapter
HDW2050	1/2" Bracket Spacer	PRT2502	Mounting Plate, standard
HDW2503	Pad Mount Plate, Short	PRT2601	E-Z Access Removable bracket base
HDW2504	Mounting Plate, short	PRT2625	Flush Pad Mount Plate
HDW2505	110 Degree Curved Back Adapter	PRT2626	1" Offset Adapter
HDW2508	100 Degree Curved Back Adapter	PRT2627	Offset Pad Mount Plate
HDW2509	Offset Pad Mounting Plate, Short	PRT2650	E-Z Access Removable bracket, male assembly
HDW2521	1" Height Adapter Plate	PRT2651	1" Offset Mount, 20 Degree
HDW2522	2" Height Adapter Plate	PRT2652	2" Offset Mount, 20 Degree
HDW2523	3" Height Adapter Plate	PRT2653	Offset Pad Mount Plate, 20°
HDW2524	1" Fixed Lateral Extension	PRT2654	Offset Pad Mount Plate, 20°, short
HDW2525	2" Fixed Lateral Extension	PRT2727	2" Offset Adapter
HDW2527	2" Distal Extension Adapter		

PRT0020

Fastener Kit

Assortment of t-nuts, screws, bolts, for mounting and adjusting Modular Hardware and other AES hardware

LI'L KIDDOS | LATERAL BRACKETS



	STANDARD FIXED	HIGH PROFILE FIXED	
FIXED BRACKETS	17/3	3%"	
FLUSH	HDW5100	HDW5102	
½" OFFSET	HDW5150	HDW5152	
1" OFFSET	HDW5110	HDW5112	

	MODULAR FIXED	REAR QUICK ADJUST	E-Z ACCESS
MODULAR BRACKETS		Up to 1½" adjustment range in ½" increments	1½" 5" 5" 130° rotation
FLUSH	HDW7540	HDW7513	HDW7500
½" OFFSET	HDW7545	HDW7553	HDW7505
1" OFFSET	HDW7640	HDW7613	HDW7600
1 ½" OFFSET	HDW7645	HDW7653	HDW7605
2" OFFSET	HDW7740	HDW7713	HDW7700

	RQA + E-Z ACCESS	REMOVABLE	REMOVABLE + 20° ADJUSTABLE
MODULAR BRACKETS	Up to 11/2" adjustment range in 1/2" increments 130° rotation (photo shown with hinge partially open)	5'	2 ½" (profile)
FLUSH	HDW7511	HDW3550	_
½" OFFSET	HDW7551	_	_
1" OFFSET	HDW7611	HDW3551	HDW3650
1 ½" OFFSET	HDW7651	_	_
2" OFFSET	HDW7711	HDW3552	HDW3750

LI'L KIDDOS HARDWARE COMPONENTS



PAD MOUNT PLATES

OFFSET (1" or greater) FLUSH OFFSET (½" only) PRT7625 PRT7627 PRT7653 PRT7624

- Add height adapter for increased profile
- For offset bracket, use Offset Pad Mount Plate and Offset Adapter
- Thoracic or Pelvic application
- Slot sized for 10-32
- Requires offset adapter (shown below)
- Thoracic or Pelvic application
- Slot sized for 10-32
- Requires 20° adjustable offset adapter (shown below)
- Use when wanting to achieve 1/2" offset *only*— cannot be combined with Offset Adapters shown in table below
- Thoracic or Pelvic application
- Slot sized for 10-32

	HEIGHT ADAPTER	OFFSET ADAPTER	20° ADJ OFFSET MOUNT
ADAPTERS	11/2"	11/2"	1½"
1/2"	_	[see PRT7624]	_
1"	PRT7521	PRT7626	PRT7651
11/2"	_	PRT7628	_
2"	PRT7522	PRT7629	PRT7652

- Use with modular brackets to accommodate thicker seats and backs
- Use with PRT7627 to achieve offset
- Thoracic or Pelvic application
- Requires PRT2653 (Offset, 20° adjustable pad mount plate)

BRACKET EXTENSION	BRACKET EXTENSION	
11/2"		
PRT7525		

- Connects pad mount plate to mounting plate.
- Use with fixed modular bracket, offset adapters, RQA only bracket, 20° adjustable bracket

MOUNTING PLATE
1½", 24"
PRT7502

E-Z ACCESS HINGE	REAR QUICK ADJUST	REMOVABLE BASE	REMOVABLE PAD MNT
11/2"	2%*	12.	11/5"
PRT7500	PRT7501	PRT7601	PRT7650

• Slots sized for 10-32


MEDIAL SUPPORT | PELVIC/THIGH SUPPORT PACKAGES

When ordering...

- **№** Specify pad shape and size; standard sizes & dimensions are listed
- ¥ ½" Medium SunMate foam is standard, please specify if other is desired
- **■** Specify upholstery fabric and color
- ¥ 14-20 T-nuts, unless otherwise noted
- **¥** Brackets made of aircraft-grade aluminum




- Includes pad and pull-away bracket. Specify pad shape and size.
- The pull-away bracket can be ordered with either $1 \frac{1}{2} 2^{\prime\prime}$ or $2 3 \frac{1}{2}^{\prime\prime}$ or $3 \frac{1}{2} 2^{\prime\prime}$ 5 ½" of height adjustment.
- To ensure bracket clearance, add $\frac{1}{2}$ to seat height for minimum adjustment
- L mount may be attached either proximally or distally.



- Includes pad and hide away bracket. Specify pad shape and size.
- To ensure bracket clearance, add ½" to seat height for minimum adjustment range.
- Hide-Away bracket may be ordered with either T mount or L mount. L'mount may be attached either proximally or distally.
- Includes pad and adjustable hide away bracket. Specify pad shape and size.
- The Adjustable Hide-Away bracket can be adjusted 2-3" or 4-5" depending on option ordered.
- To ensure bracket clearance, add 1/2" to seat height for minimum adjustment range.
- L mount may be attached either proximally or distally.

MEDIAL SUPPORT | PELVIC/THIGH SUPPORT PACKAGES



- Includes pad and swing-away bracket. Specify pad shape and size.
- To ensure bracket clearance, add ½" to seat height for minimum adjustment range
- Requires abductor height of 2" or more for pad attachment
- Pad attaches with wood screws
- Not recommended for use with 1/4" plywood seat base or thin pan base
- Includes pad and adjustable swing-away bracket. Specify pad shape and size.
- The Adjustable swing-away bracket offers either 2-3" or 3-4" of height
- To ensure bracket clearance, add $\frac{1}{2}$ " to seat height for minimum adjustment range
- ullet Not recommended for use with $\ensuremath{\ensuremath{\ensuremath{\mathcal{U}}}}$ plywood seat base or thin pan base
- Keystock requires mounting with wood screws

LI'L KIDDOS SUPPORT

MEDIAL THIGH



When ordering:

- **■** Specify pad size and shape; standard sizes are listed below
- № Medium Sunmate foam is standard; specify if other is desired
- **№** Specify upholstery fabric & color



ADJUSTMENT RANGE	11/4 - 23/4	23/4 - 41/4
STANDARD PAD PACKAGE	MTK31XX	MTK34XX
CUSTOM PAD PACKAGE	MTK33XX	MTK35XX


- To ensure bracket clearance, add 1/2" to seat height when ordering
- L Mount may be attached proximally or distally
- To ensure bracket clearance, add 1/2" to seat height when ordering
- L Mount may be attached proximally or distally




• To ensure bracket clearance, add 1/2" to seat height when ordering

MEDIAL SUPPORT | PADS

Medial Support pads provide a soft contact surface for the client and help keep client's lower extremities properly positioned to minimize adduction. Shape preference (wedge or oval) is determined by client's leg positioning needs and comfort.



WEDGE		А	В	С	D
S	ABD10XX	3	2	3	4
М	ABD11XX	4 1/2	2	3	5
L	ABD12XX	5	3	4	5
Custom	ABD19XX	Specify			



OVAL		А	В	С
XS	ABD37XX	2 1/2	3 1/2	3
S	ABD30XX	3	4	3
М	ABD31XX	4 1/2	5	3
L	ABD32XX	5	5	4
Custom	ABD39XX	Specify		,

MEDIAL SUPPORT | MODIFICATIONS

FOAM MODIFICATIONS

MOD3555 Replace foam ¼" Gel in Medial Support

MOD3655 ½" T-Foam in Support MOD3755 ½" Pudgee in Support

MOD3856 1" Sunmate in Support

MOD3857 Incontinence Film covering **Custom Foam Modification** MOD3859

MEDIAL SUPPORT | BRACKETS

	PULL-AWAY		
Adjustment			
11/2-2"	HDW4420		
2-3 ½"	HDW4435		
31/2-51/2"	HDW4455		

- To ensure bracket clearance, add 1/2" to seat height for minimum adjustment range
- L mount may be attached either proximally or distally

	HIDEAWAY W/ L MOUNT	HIDEAWAY W/ T MOUNT
Mounting Height		
1½"	HDW4215	HDW4216
21/2"	HDW4225	HDW4226
3"	HDW4230	HDW4231

- To ensure bracket clearance, add ½" to seat height for minimum adjustment range
- L mount may be attached either proximally or distally
- To ensure bracket clearance, add 1/2" to seat height for minimum adjustment range

	ADJUSTABLE HIDEAWAY	ADJUSTABLE SWING AWAY
Adjustment		
2-3"	HDW4245	HDW4115
3-4"	HDW4250	HDW4116

- To ensure bracket clearance, add ½" to seat height for minimum adjustment range
- To ensure bracket clearance, add ½" to seat height for minimum adjustment range
- Not recommended for use with 1/4" plywood seat base, or thin pan base
- Keystock requires mounting with wood screws

MEDIAL SUPPORT | BRACKETS



- To ensure bracket clearance, add $\frac{1}{2}$ " to seat To ensure bracket clearance, add $\frac{1}{2}$ " to seat height for minimum adjustment range
- Not recommended for use with ¼" plywood Not recommended for use with ¼" plywood seat base or thin pan base
- Keystock requires mounting with wood screws
- height for minimum adjustment range
- seat base or thin pan base
- Keystock requires mounting with wood screws

MEDIAL SUPPORT | BRACKET COMPONENTS



HIDEAWAY BASE
57/8"
PRT4200

	L MOUNT	T MOUNT
Mounting Height	3"	Center to center
1 ½"	PRT4215	PRT4216
2"	PRT4221	PRT4220
2 ½"	PRT4225	PRT4222
3"	PRT4230	PRT4231

HEAD SUPPORT | PACKAGES

Please reference Whitmyer Product Catalog (MKT3000-3) for full offering of Whitmyer head supports.



HSK2080

Flat Head Support Pad, Lift-Away Bracket

- Specify pad dimensions
- Specify upholstery fabric and color
- Available with Pediatric or Adult bracket
- 1" HR foam standard



HSK6832

Headrest Extension Pad w/ Adjustable Head Support bracket, 12"

- Specify pad dimensions
- Specify upholstery fabric and color
- 11/2" HR foam is standard—upgrading to 2" HR foam is recommended to allow for mounting offset




HSK6833

Curved Headrest Extension Pad w/ Adjustable Head Support bracket, 12"

- Specify pad dimensions
- Specify upholstery fabric and color
- Designed to match Curved Back shape
- 1" HR foam is standard—upgrading to 11/2" HR foam is recommended to allow for mounting offset


HEAD SUPPORT PADS



U STYLE HEADREST

HRU19XX Custom, Soft HRU29XX Custom, Reinforced

- Specify dimensions as referenced on drawing
- Specify upholstery fabric and color
- ½" Medium Sunmate is standard



HEADREST EXTENSION PAD

HRE12XX Flat, up to 18Wx10H HRE13XX Flat, up to 20W x15H

HRE17XX Curved, 18Wx10H

HRE18XX Curved, up to 20W x15H

HRE19XX Flat or Curved over 20Wx15H

- Requires additional adapter to interface with adjustable hardware
- Specify upholstery fabric and color
- 1" HR foam standard



LATERAL/POSTERIOR HEADREST **PADS**

HRP6831 Set of three pads HRP6839 Set of three pads, custom

- Specify pad size standard: 4 x 4, 5 x 5, 6 x 6, 6 ½ x 7 Specify dimensions if ordering custom
- Specify upholstery fabric and color
- Requires additional adapter to interface with adjustable hardware
- Hinge assembly not included
- 1/2" Medium Sunmate is standard

HEAD SUPPORT | MODIFICATIONS

MOD6100 Add Occipital Roll to Headrest MOD6150 Add Ear Relief to Headrest Pads

MOD6135 1/2" Pudgee in Headrest MOD6155 T-Foam in Headrest

HEAD SUPPORT | HARDWARE



HDW7880

Lift Away Head Support Bracket, Pediatric

• 91/2" height



HDW7883

Lift Away Head Support Bracket, Adult

• 12" height




HDW6830

Adjustable Posterior Head Support Bracket


• 12" height


HEAD SUPPORT | HARDWARE


COMPONENTS



	·
Standard	PRT6836
Pediatric	PRT6835
Hinged	PRT6834







HEADREST U BRACKET



Standard	PRT6833
Wide	PRT6838M
Pediatric	PRT6833P






• Center-to-center holes is 2"

TRUNK & PELVIC SUPPORT STRAPS

Please reference JAY Positioning Supports Catalog (MK-100093) for full offering of anterior trunk supports.



PELVIC STRAPS— **NEOPRENE WITH CAM BUCKLE**

		PAD W	PAD L	STRAP length	TOTAL Length
S	PSO40XX	1"	5 1/2"	10 ½"	29"
М	PSO41XX	1 3/4"	8 1/2"	10 ½"	31"
L	PSO42XX	2"	11 ½"	12"	36"

SOFT SUPPORT COMPONENTS & IARDWARE



BELT ANCHORS

	WEBBING SIZE
PRT6925	1″
PRT6926	1½"
PRT6927	2"




HDW6900 STRAP POSITIONER

• 1" webbing only

HARD300

Long Version



BELT PADDING


	Webbing	W	L
APB1801	1"	1 ½"	5 1/2"
APB1802	1"	11/2"	7 1/2"
APB1803	1½"	2"	7 1/2"
APB1804	1½"	21/4"	9 1/2"
APB1805	2"	23/4"	9 1/2"
APB1806	2"	23/4"	111/2"
APB1809	Custom: Specify dimensions		

Black Neoprene Standard.

MOD7001 Synthetic Sheepskin addition to Pads

LOWER EXTREMITY SUPPORT


Please reference JAY Positioning Supports Catalog (MK-100093) for full offering of foot boxes and shoe holders.



SINGLE LEG/FOOT BOX

			Footb	ox Dime	nsions
LEFT	RIGHT	W/C WIDTH	W	Н	D
TST7014L	TST7014R	14"	51/4"	9"	9"
TST7016L	TST7016R	16"	61/4"	10"	10"
TST7018L	TST7018R	18"	71/4"	10"	10"
TST7000		Custom; specify dimensions			

- OUTSIDE dimensions are referenced in table above
- Sold singly; order one each L (left) and R (right) for complete pair
- ½" HR foam is standard
- Hardware is NOT included
- Black vinyl upholstery is standard on contact surface and Veltex on non-contact surface
- Specify T-nut placement



FULL WIDTH LEG/FOOT BOX

		Footbox Dimensions		
	W/C WIDTH	W	Н	D
TST7114	14"	11"	9"	9"
TST7116	16"	13"	10"	10"
TST7118	18" 15" 10"			10"
TST7100	Custom; specify dimensions			

- OUTSIDE dimensions are referenced in table above
- ½" HR foam is standard
- Hardware is NOT included
- Black vinyl upholstery is standard on contact surface and Veltex on non-contact surface
- Specify T-nut placement



STABILIZER PACKAGES

KPG10XX	Stabilizer & Pull-away bracket Package, Standard
KPG19XX	Stabilizer & Pull-away bracke Package, Custom

- Provides pelvic stabilization through bilateral distal support of the femoral movement
- Design based on the clinical experience of Nancy Schuster, MS, OTR, ATP

STABILIZER BRACKET

HDW4840 Stabilizer Pull-away bracket

STABILIZER PADS

-	
KPB1014	Stabilizer Pad, 14" W
KPB1016	Stabilizer Pad, 16" W
KPB1018	Stabilizer Pad, 18" W
KPB19XX	Stabilizer Pad, Custom




KNEE ADDUCTOR PACKAGES

	LEFT	RIGHT
Short Arm	KAK21XXL	KAK21XXR
Long Arm	KAK31XXL	KAK31XXR
Custom	KAK39XXL	KAK39XXR

- Swingaway bracket
- Pad base is PlastiTech† material
- Pad foam is 1" Medium Sunmate
- Four-Way Stretch upholstery
- 3" Round Pad is Standard

†PlastiTech is a waterproof plastic base material that accepts wood screws and T-Nuts.




KNEE ADDUCTOR BRACKET


	LEFT	RIGHT
Short Arm	HDW4853	HDW4854
Long Arm	HDW4851	HDW4852

KNEE ADDUCTOR PADS

KAD10XX	Knee Adductor Pad, 3" Round
KAD11XX	Knee Adductor Pad, 4" Round
KAD19XX	Knee Adductor Pad Custom

LOWER EXTREMITY SUPPORT





PADDED CALF PANEL

		W	Н	W/C
S	PLS3010 14" 8"		16"	
М	PLS3020 16" 8"		18"	
L	PLS3030	18"	8"	20"
CUSTOM	PLS3090	Specify		fy

• Black Fabric coverd neoprene face, Black Vinyl on back side

TOE/ANKLE STRAPS WITH GROMMET ATTACHMENT

	PADDED	NON-PADDED	ADJUSTMENT RANGE
S	TST4000	TST1000	6 to 9"
М	TST4200	TST1500	8 to 12"
L	TST4500	TST2000	13 to 15"



LEG EXTENSION

	W	L
LRE10XX	14"	6
LRE11XX	16"	8
LRE12XX	18"	10
LRE19XX	Custom	

- Specify upholstery fabric
- 1½" HR foam standard
- T-nuts and hardware are NOT included



LEG EXTENSION MOUNTING BRACKET

HDW6828 Adjustable, Removable

COMPONENTS:


A PRT6828

B HDW6826L



LEG-REST KNEE PADS

TST0500 High Profile, 3" Diameter



LEG STRAP

PLS1000 Single, Specify W/C Size PLS2000 H-Style, Specify Size

UPPER EXTREMITY SUPPORT ARM & WRIST SUPPORT

Please reference JAY Positioning Supports Catalog (MK-100093) for full offering of upper extremity positioning supports.




ARM TROUGH PACKAGES

32440 Contoured Arm Trough with Attaching Hdw package

32450 Custom Contoured Arm Trough with Attaching Hdw package

- ABS base material
- Includes ¼ 20 T-nuts
- Specify Arm trough size when ordering; refer to the table at right for standard sizes



CONTOURED ARM TROUGH

		W	L
S	CAT5510	4"	13 ¼"
М	CAT5520	4 1/2"	16 1/4"
L	CAT5530	5"	19 1/4"
CUSTOM	JSTOM CAT5590 Specify		cify

- Outside dimensions are referenced
- Black neoprene is standard upholstery fabric
- Hardware not included
- ABS base material
- Two self-gripping (soft hook & loop) straps included



ARM TROUGH HARDWARE

32430 Arm Trough Attaching Hardware




PADDED WRIST STRAPS

		w	L	WRIST CIRCUMFERENCE
S	ATS7510	2 1/2"	10"	5½ to 7"
М	ATS7520	3"	12"	7½ to 9
L	ATS7530	3"	14	9½ to 11"
CUSTOM	ATS7590	Spe	cify	Specify

• Pads are constructed of vinyl-covered ComforFit material

MULTI-ADJUSTABLE ARM



MULTI-ADJUSTABLE ARM KIT

MAS2000

Standard Multi-adjustable Arm Kit


- Mounting opening adjusts from 3/4" round tube to $3^{3}/4$ " rectangular tube
- Includes set screws to lock in place
- Vertical extension range is 5" up to 20 1/2"
- Horizontal extension range is 4" up to 191/2"


MULTI-ADJUSTABLE ARM COMPONENTS


† 🚯 MAS2105	5" Long Rod Extension
† MAS2111	7.5" Long Rod Extension
† MAS2115	12" Long Rod Extension
† B MAS2120	16.5" Long Rod Extension
† () MAS2205	Short Adjusting Tube
† (1) MAS2210	Long Adjusting Tube
† 📵 MAS2300	Base Mount Assembly
† 🕞 MAS2305	90 Degree Tube Connector
† (MAS2310	Support Platform Mount Assembly
† 🚯 MAS2329	Support Platform (12W x 9 t/8"H), Standard
† () MAS2340	Quick Release Adapter Mount Block
† ① MAS2342	Detent Ring
† MAS2341	Ball Mount Tension Ring
MAS2320	Custom Support Platform Specify dimensions
MAS2124	21" Long Rod Extension
MAS2321	Mounting Plate for Dynavox, 5.0" x 2.5" Mounts in arm tube
MAS3000	Levo Application Kit Requires drilling arm tube

† Part is included in Kit (MAS2000)

PLATFORMS (TRAYS)







TRAY SPECIFICATIONS

	Α	В	С	D
XS	21	18	10	7
S	21	20	12	7
М	24	22	14	9
L	26	24	16	10
XL	28	26	18	10

WOOD PLATFORMS W/STRAPS


S	WCT0500
М	WCT0600
L	WCT0700
XL	WCT0800

- Hardware NOT included
- Hook & Loop straps included
- Does NOT include pencil edge

WOOD PLATFORMS W/ PENCIL **EDGE**

	NO SLOTS	SLOTS
XS	WCT1005	WCT1005S
S	WCT1000	WCT1000S
М	WCT2000	WCT2000S
L	WCT3000	WCT3000S
XL	WCT4000	WCT4000S


• Hardware NOT included



WOOD PLATFORMS WITH 2" EDGE

	NO SLOTS	SLOTS
XS	WCT1025	WCT1025S
S	WCT1020	WCT1020S
М	WCT2020	WCT2020S
L	WCT3020	WCT3020S
XL	WCT4020	WCT4020S


• Hardware NOT included



POLYCARBONATE PLATFORMS W/ PENCIL EDGE

	NO SLOTS	SLOTS
XS	WCT5005	WCT5005S
S	WCT5000	WCT5000S
М	WCT6000	WCT6000S
L	WCT7000	WCT7000S
XL	WCT8000	WCT8000S

• Hardware NOT included



POLYCARBONATE PLATFORMS WITH 2" EDGE

	NO SLOTS	SLOTS
XS	WCT5025	WCT5025S
S	WCT5020	WCT5020S
М	WCT6020	WCT6020S
L	WCT7020	WCT7020S
XL	WCT8020	WCT8020S

• Hardware NOT included

CUSTOM PLATFORM

WCT9000 Custom Upper Extremity Support Platform (no slots) WCT9000S Custom Upper Extremity Support Platformw/ Slots

- Specify dimensions & material
- Hardware NOT included

PLATFORM MOUNTING HARDWARE



THW1000 Toggle Hardware

- Slots required on tray for mounting
- Power Joystick option



THW2000 Slide Tube

- Mounts to armrest no slots required
- Tray attaches to tube arm with heavy-duty hook & loop strips



Z Slide Hardware

THW3600 11/8" Armpad Height THW3650 11/2" Armpad Height



THW3000 Wood Runners

- Runners include THW6000
- No slots required
- Profile 11/4"

THW3000T Wood Runners

- Runners include THW6000
- No slots required
- Profile 1½"



THW3100 S Style Runners

- Runners include THW6000
- Constructed of black ABS
- No slots required
- Profile 11/4"
- Runners include THW6000



THW3210 Slotted Channel Runner

- Runners include THW6000
- Constructed of black ABS
- No slots required
- Profile 11/4"
- Runners include THW6000



THW3500 Adjustable Lever

- No slots required
- Short adjustment rails are optional upgrade



THW3700 Top Drop Hardware


- No slots required
- Profile is adjustable to accommodate varying arm pad thicknesses



THW4500 Angle-Cam

- Slots required on tray for mounting
- Power Joystick option

PLATFORM ACCESSORIES



TAC10XX Protractor Wings

- Specify lateral or posterior mount
- Specify height and width
- Specify upholstery material, standard is black vinyl
- Hardware is not included



TAC11XX Removable Cover, Full

• Pencil Edge is omitted on full covers

TAC12XX Padding at Cutout Edge only

• Pencil Edge is included

TAC13XX Custom

TAC14XX Padded Cover, Anterior Half


• Pencil Edge is NOT included

TAC15XX Padded Cover, Posterior Half (Shown)

• Pencil Edge is included

TAC5600 Elbow/Forearm Pad

TAC5619 Elbow/Forearm Pad, Custom



TAC2000 Polycarbonate Easel

- Angle adjustable
- Removable



T-HANDLES

JTH1000 Short

• Handle height is 1"

JTH2000 Tall

• Handle height is 11/2"




TRAY STRAPS

THW5000 Hook & Loop D-Ring Release THW6000 Side Release Buckle

- Plastic 1" side release buckle
- Strap is 44" in total length

UPPER EXTREMITY SUPPORT

PLATFORM MODIFICATIONS




TAC4100L Joystick Cutout, Left TAC4100R Joystick Cutout, Right

STANDARD CUTOUT DIMENSIONS

	Α	В
XS	9"	4"
S	10"	4"
М	12"	4"
L	14"	4"
XL	16"	4"

TAC4150 Joystick Midline Cutout TAC4175 Joystick Cutout, Custom TAC4000 Window Insert

MISCELLANEOUS



Wheelchair Canopy - 12" - 15"

- Includes observation window
- Specify wheelchair width when orderingcanopy designed to fit chair widths from 12" to 15" wide
- Only available in vinyl, specify color

99131 - RT


Wheelchair Canopy - 16" - 18"

- Includes observation window
- Specify wheelchair width when ordering— canopy designed to fit chair widths from 16" to 18" wide
- Only available in vinyl, specify color




DATE	ACCT#	PO#
USER NAME	DEALER	CONTACT
PHONE	FAX	EMAIL
CLIENT MEASUREMENTS B CHARACTER STATE OF THE STATE OF T		A. SHOULDER WIDTH B. CHEST WIDTH C. HIP WIDTH D. WIDTH at KNEE E. SEAT to TOP of HEAD F. SEAT to TOP of SHOULDER G. SEAT to AXILLA H. SEAT to PSIS K. CHEST DEPTH
		L. BACK to ANTERIOR of ITs M. SEAT DEPTH N. SEAT to FOOT PLATE
INSERT Dimensions & Specif	ications	
	A. Seat Width B. Back Width C. Distance between Thoracic Laterals D. Distance between Pelvic Laterals E. Usable Seat Depth F. Total Seat Depth G. Back Height H. Seat to top of Lateral Thoracic Pad	INSERT INFORMATION Full Width Mount Back Seat Between Rail Mount Back Seat MONOGRAM: MOD5900 Name in BLOCK letters MOD5901 Name in Script letters MOD5902 Stock Design UPHOLSTERY FABRIC: Client Contact Surface Non-contact Surface: Black Vinyl (Standard) Other (Specify base material and color)
		· ·
MOBILITY BASE Dimensions	s & Specifications (Required)	
MAKE:	Seat· 🔲 🔲	INSTALLATION OPTIONS: ☐ Install System to Frame ☐ Install Transport Option Hardware MOD6800 Applies to Manual chairs only; requires frame
COLOR:		to be drop-shipped to AES
	PUSH HANDLE HEIGHT: W/C ARMREST TYPE: ☐ Flip-up ☐ Height Adjustable (mounting off seat)	

SEAT DIMENSION Specifications




SOLID SEAT INSERT

L	Α	В	
_			




I STYLE SEAT INSERT

Α	E	
В	F	
С	G	
D		



T STYLE SEAT INSERT


Α	С	
В	D	



OFFSET SEAT INSERT


□ LONG LEFT SIDE□ LONG RIGHT SIDE

Α	С	
В	D	




ANTI-THRUST SEAT INSERT

Α	D	
В	E	
С		




ADJUSTABLE SPLIT SEAT INSERT

Α	D	
В	E	
C	F	



BUILT-IN MEDIAL THIGH SUPPORT SEAT INSERT


Α	D	
В	E	
С		



CONTOUR SEAT INSERT


Α	G	
В	Н	
С	I	
D	J	
E	K	
F	L	

BACK DIMENSION Specifications




SOLID BACK INSERT

	Α	В	
ι			




I STYLE BACK INSERT

Α	D	
В	E	
С	F	




T STYLE BACK INSERT

Α	С	
В	D	




BI-ANGULAR BACK INSERT

Α	С	
В		




CURVED BACK INSERT

A B



ADJUSTABLE LUMBAR SUPPORT **BACK INSERT**

Α	С	
В	D	



CONTOUR BACK INSERT

Α	D	
В	E	
С		





