


# Rogue XP Aluminum Anodized Handrim



For further information see *Handrim Hardware Chart*. Quantity in *Handrim Hardware Kits* are for one wheel.<p>  
 Tab mounted handrims must be installed prior to the tire being installed.

Item Number	Part Number	Description	Quantity	Retail
8, 9, 10	505263	Handrim Hardware Kit, 3 Rivet Nut, EACH <i>Used with 18" Ki Spoke and 18" MAXX Spoke wheels</i>	1	\$19.31
7a	100206	Handrim Aluminum Anodized 3 Rivet Nut 18" <i>Used with 18" Ki Spoke and 18" MAXX Spoke wheels</i>	1	\$84.95
7a	502245	Handrim Aluminum Anodized 3 Rivet Nut 18" w/ Micro Grip <i>Used with 18" Ki Spoke and 18" MAXX Spoke wheels</i>	1	\$112.35
8, 9, 10	502168	Handrim Hardware Kit, 6 Rivet Nut, EACH <i>Used with Ki Spoke and MAXX Spoke wheels</i>	1	\$38.61
7b	200536	Handrim Aluminum Anodized 6 Rivet Nut 20" <i>Used with 20" Ki Spoke and 20" MAXX Spoke wheels</i>	1	\$84.95
7b	502246	Handrim Aluminum Anodized 6 Rivet Nut 20" w/ Micro Grip <i>Used with 20" Ki Spoke and 20" MAXX Spoke wheels</i>	1	\$133.11

Item Number	Part Number	Description	Quantity	Retail
7c	100560	Handrim Aluminum Anodized 6 Rivet Nut 22" <i>Used with 22" Ki Spoke and 22" MAXX Spoke wheels</i>	1	\$77.22
7c	502247	Handrim Aluminum Anodized 6 Rivet Nut 22" w/ Micro Grip <i>Used with 22" Ki Spoke and 22" MAXX Spoke wheels</i>	1	\$112.35
7d	200349	Handrim Hard Anodized 6 Rivet Nut 24" <i>Used with 24" Ki Spoke and 24" MAXX Spoke wheels</i>	1	\$82.74
7d	502248	Handrim Aluminum Hard Anodized 6 Rivet Nut 24" w/ Micro Grip <i>Used with 24" Ki Spoke and 24" MAXX Spoke wheels</i>	1	\$124.53
7e	200350	Handrim Hard Anodized 6 Rivet Nut 25" <i>Used with 25" Ki Spoke and 25" MAXX Spoke wheels</i>	1	\$82.74
7e	502249	Handrim Aluminum Hard Anodized 6 Rivet Nut 25" w/ Micro Grip <i>Used with 25" Ki Spoke and 25" MAXX Spoke wheels</i>	1	\$112.35
5, 6	502171	Handrim Hardware Kit, 5 Rivet Nut, EACH <i>Used with X Core Mag wheels</i>	1	\$33.10
4a	101897	Handrim Aluminum Anodized 5 Rivet Nut 20" <i>Used with 20" 5-Spoke X Core Mag wheels</i>	1	\$71.70
4a	502242	Handrim Aluminum Anodized 5 Rivet Nut 20" w/ Micro Grip <i>Used with 20" 5-Spoke X Core Mag wheels</i>	1	\$113.32
4b	101898	Handrim Aluminum Anodized 5 Rivet Nut 22" <i>Used with 22" 5-Spoke X Core Mag wheels</i>	1	\$71.70
4b	502243	Handrim Aluminum Anodized 5 Rivet Nut 22" w/ Micro Grip <i>Used with 22" 5-Spoke X Core Mag wheels</i>	1	\$113.32
4c	100975	Handrim Aluminum Anodized 5 Rivet Nut 24" <i>Used with 24" 5-Spoke X Core Mag wheels</i>	1	\$60.67
4c	502244	Handrim Aluminum Anodized 5 Rivet Nut 24" w/ Micro Grip <i>Used with 24" 5-Spoke X Core Mag wheels</i>	1	\$112.35
2a, 3	502174	Handrim Hardware Kit, 6 Tab, EACH <i>Used with Spinergy and Topolino wheels</i>	1	\$16.55
2b	502170	Handrim Hardware Kit, Superlight, EACH <i>Used with Superlight Spoke wheels</i>	1	\$16.55
1a	100754	Handrim Aluminum Anodized 6 Tab 24" <i>Used with 24" Superlight Spoke and Topolino wheels</i>	1	\$88.25
1a	502236	Handrim Aluminum Anodized 6 Tab 24" w/ Micro Grip <i>Used with 24" Superlight Spoke and Topolino wheels</i>	1	\$124.53
1b	502237	Handrim Aluminum Anodized 6 Tab 25" w/ Micro Grip <i>Used with 25" Superlight Spoke and Topolino wheels</i>	1	\$112.35

Item Number	Part Number	Description	Quantity	Retail
1c	103125	Handrim Aluminum Anodized Spinerger 6 Tab 20" <i>Used with 20" Spinerger Spox wheels</i>	1	\$93.77
1d	505584	Handrim Aluminum Anodized Spinerger 6 Tab 20" w/ Micro Grip <i>Used with 20" Spinerger Spox wheels</i>	1	\$126.53
1e	100827	Handrim Aluminum Anodized Spinerger 6 Tab 22" <i>Used with 22" Spinerger wheels</i>	1	\$93.77
1f	502238	Handrim Aluminum Anodized Spinerger 6 Tab 22" w/ Micro Grip <i>Used with 22" Spinerger wheels</i>	1	\$127.11
1g	100766	Handrim Aluminum Anodized Spinerger 6 Tab 24" <i>Used with 24" Spinerger wheels</i>	1	\$93.77
1h	502239	Handrim Aluminum Anodized Spinerger 6 Tab 24" w/ Micro Grip <i>Used with 24" Spinerger wheels</i>	1	\$127.09
1i	100767	Handrim Aluminum Anodized Spinerger 6 Tab 25" <i>Used with 25" Spinerger wheels</i>	1	\$93.77
1j	502240	Handrim Aluminum Anodized Spinerger 6 Tab 25" w/ Micro Grip <i>Used with 25" Spinerger wheels</i>	1	\$127.09
2a	100669	M5 x 10 BTN SHCS BLKZC w/Patch <i>Used with Spinerger and Topolino wheels</i>	6	\$1.66
2b	100536	M5 x 12 BTN SHCS ZC <i>Used with Superlight Spoke wheels</i>	6	\$2.21
3	100657	M5 Nylock Nut Blk Zn <i>Used with Spinerger and Topolino wheels</i>	6	\$0.56
5	100629	Handrim Spacer Alum Ano Riv <i>Used with X Core Mag wheels</i>	5	\$1.16
6	100724	M5 x 55 PH TRS HD M/S ZC <i>Used with X Core Mag wheels</i>	5	\$1.10
8	100653	Handrim Spacer Spoke <i>Used with Ki Spoke and MAXX Spoke wheels</i>	6	\$1.10
9	100698	Handrim Tab <i>Used with Ki Spoke and MAXX Spoke wheels. Tabs are included / installed with MAXX Spoke wheels.</i>	6	\$4.17
10	100654	M5 x 40 DIN7985 PH Pan HD M/S ZC <i>Used with Ki Spoke and MAXX Spoke wheels</i>	6	\$1.66

# Handrim Hardware Chart

Wheel		Handrim Connection				Handrim							
Wheel	Wheel Part No.	Connection Points	Tab	Spacer*	Screw*	Aluminum Anodized	Superlight	Plastic Coated	Projection	Natural Fit Standard	Natural Fit LT	Flex Rim	
18" Ki Spoke	200529	3	100698	100653	Aluminum, Plastic Ctd, Projection: 100654  Natural Fit: 100835	100206		101106					
20" Ki Spoke	200530					200536		200542					
22" Ki Spoke	200531					100560		100576	100569	200538	200201		
24" Ki Spoke	200532					200349	101161	100577	200547	100793	200202		
25" Ki Spoke	200533					200350		101870	200548	200539	200540		
26" Ki Spoke	200534					200351		100578	200549	100907	101454		
24" Superlight	101159	6			100536	100754	101161	100836		100830	100828		
25" Superlight	101160							101160	101091		101464	101460	
20" 5-Spoke X Core	101961	5		Aluminum, Projection, Natural Fit: 100629 Plastic Ctd: 101756	100724	101897		101963					
22" 5-Spoke X Core	101962						101898		101964				
24" 5-Spoke X Core	100960						100975		100976	200546			
24" 5-Spoke X Core	100960						101893				100768	100769	
20" Spinergy Spox	See Spinergy Spox Page	6			Screw: 100669 Nut: 100657	103125		103179					
22" Spinergy Spox						100827		100808		100889	100888		
24" Spinergy Spox						100766	101161	100615		100830	100828		
25" Spinergy Spox						100767	101160	100765		101464	101460		
26" Spinergy Spox						101477		101148		200200	100950		
22" Spinergy LX	See Spinergy LX Page	6			Screw: 100669 Nut: 100657	100827		100808		100889	100888		
24" Spinergy LX						100766	101161	100615		100830	100828	200213	
25" Spinergy LX						100767	101160	100765		101464	101460	200208	
26" Spinergy LX						101477		101148		200200	100950		
18" MAXX Spoke	109285	6	100698 (Not used w/ Superlight Handrim)	100653 (Not used w/ Superlight Handrim)	Aluminum, Plastic Ctd, Projection: 100654  Natural Fit: 100835  Superlight: Screw: 100669 Nut: 100657	100206		101106					
20" MAXX Spoke	108244					200536		200542					
22" MAXX Spoke	105135					100560		100576	100569	200538	200201		
24" MAXX Spoke	105136					200349	101161	100577	200547	100793	200202		
25" MAXX Spoke	107436					200350		101870	200548	200539	200540		
26" MAXX Spoke	107437					200351		100578	200549	100907	101454		

\* Spacer and screw part numbers listed in the chart are for standard handrim mounting. If using close mount handrim mounting on aluminum anodized, plastic coated or projection handrims, use part number 000240 for the spacer, part number 102616 for the screw on Ki Spoke and MAXX Spoke wheels or part number 100666 for the screw on 5-Spoke X Core wheels.