


Catalyst / Spark Seat Pan (Transit)



Part Number	Tape Size	Seat Pan Depth
004458	3" x 4"	10" - 11"
004458 & 004459	3" x 4" & 3" x 5"	12"
004459	3" x 5"	13"
004460	3" x 6"	14" - 22"

This is an Aluminum Seat Pan for Transit Use. Requires Seating System Cross Braces. The width and depth of the seat pan will measure approximately 0.5" shorter than the listed description.

Item Number	Part Number	Description	Quantity
1 - 6a, 6b, 7 - 10	120392	Transit Seat Pan Assembly, 10 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120393	Transit Seat Pan Assembly, 10 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120394	Transit Seat Pan Assembly, 10 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120395	Transit Seat Pan Assembly, 10 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120396	Transit Seat Pan Assembly, 10 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120397	Transit Seat Pan Assembly, 10 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120398	Transit Seat Pan Assembly, 10 x 18, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120399	Transit Seat Pan Assembly, 11 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120400	Transit Seat Pan Assembly, 11 x 13, Spark/Catalyst	1

Item Number	Part Number	Description	Quantity
1 - 5, 6c, 7 - 10	120401	Transit Seat Pan Assembly, 11 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120402	Transit Seat Pan Assembly, 11 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120403	Transit Seat Pan Assembly, 11 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120404	Transit Seat Pan Assembly, 11 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120405	Transit Seat Pan Assembly, 11 x 18, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120406	Transit Seat Pan Assembly, 12 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120407	Transit Seat Pan Assembly, 12 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120408	Transit Seat Pan Assembly, 12 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120409	Transit Seat Pan Assembly, 12 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120410	Transit Seat Pan Assembly, 12 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120411	Transit Seat Pan Assembly, 12 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120412	Transit Seat Pan Assembly, 12 x 18, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120413	Transit Seat Pan Assembly, 13 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120414	Transit Seat Pan Assembly, 13 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120415	Transit Seat Pan Assembly, 13 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120416	Transit Seat Pan Assembly, 13 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120417	Transit Seat Pan Assembly, 13 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120418	Transit Seat Pan Assembly, 13 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120419	Transit Seat Pan Assembly, 13 x 18, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120420	Transit Seat Pan Assembly, 14 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120421	Transit Seat Pan Assembly, 14 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120422	Transit Seat Pan Assembly, 14 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120423	Transit Seat Pan Assembly, 14 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120424	Transit Seat Pan Assembly, 14 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120425	Transit Seat Pan Assembly, 14 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120426	Transit Seat Pan Assembly, 14 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120427	Transit Seat Pan Assembly, 14 x 19, Spark/Catalyst	1

Item Number	Part Number	Description	Quantity
1 - 5, 6c, 7 - 10	120428	Transit Seat Pan Assembly, 14 x 20, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120429	Transit Seat Pan Assembly, 15 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120430	Transit Seat Pan Assembly, 15 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120431	Transit Seat Pan Assembly, 15 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120432	Transit Seat Pan Assembly, 15 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120433	Transit Seat Pan Assembly, 15 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120434	Transit Seat Pan Assembly, 15 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120435	Transit Seat Pan Assembly, 15 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120436	Transit Seat Pan Assembly, 15 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120437	Transit Seat Pan Assembly, 15 x 20, Spark/Catalyst	1
1 - 6a, 6b, 7 - 10	120438	Transit Seat Pan Assembly, 16 x 12, Spark/Catalyst	1
1 - 5, 6b, 7 - 10	120439	Transit Seat Pan Assembly, 16 x 13, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120440	Transit Seat Pan Assembly, 16 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120441	Transit Seat Pan Assembly, 16 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120442	Transit Seat Pan Assembly, 16 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120443	Transit Seat Pan Assembly, 16 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120444	Transit Seat Pan Assembly, 16 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120445	Transit Seat Pan Assembly, 16 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120446	Transit Seat Pan Assembly, 16 x 20, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120447	Transit Seat Pan Assembly, 17 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120448	Transit Seat Pan Assembly, 17 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120449	Transit Seat Pan Assembly, 17 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120450	Transit Seat Pan Assembly, 17 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120451	Transit Seat Pan Assembly, 17 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120452	Transit Seat Pan Assembly, 17 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120453	Transit Seat Pan Assembly, 17 x 20, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120454	Transit Seat Pan Assembly, 18 x 14, Spark/Catalyst	1

Item Number	Part Number	Description	Quantity
1 - 5, 6c, 7 - 10	120455	Transit Seat Pan Assembly, 18 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120456	Transit Seat Pan Assembly, 18 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120457	Transit Seat Pan Assembly, 18 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120458	Transit Seat Pan Assembly, 18 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120459	Transit Seat Pan Assembly, 18 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120460	Transit Seat Pan Assembly, 18 x 20, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120461	Transit Seat Pan Assembly, 19 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120462	Transit Seat Pan Assembly, 19 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120463	Transit Seat Pan Assembly, 19 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120464	Transit Seat Pan Assembly, 19 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120465	Transit Seat Pan Assembly, 19 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120466	Transit Seat Pan Assembly, 19 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120467	Transit Seat Pan Assembly, 19 x 20, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120468	Transit Seat Pan Assembly, 20 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120469	Transit Seat Pan Assembly, 20 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120470	Transit Seat Pan Assembly, 20 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120471	Transit Seat Pan Assembly, 20 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120472	Transit Seat Pan Assembly, 20 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120473	Transit Seat Pan Assembly, 20 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120474	Transit Seat Pan Assembly, 20 x 20, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120475	Transit Seat Pan Assembly, 21 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120476	Transit Seat Pan Assembly, 21 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120477	Transit Seat Pan Assembly, 21 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120478	Transit Seat Pan Assembly, 21 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120479	Transit Seat Pan Assembly, 21 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120480	Transit Seat Pan Assembly, 21 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120481	Transit Seat Pan Assembly, 21 x 20, Spark/Catalyst	1

Item Number	Part Number	Description	Quantity
1 - 5, 6c, 7 - 10	120482	Transit Seat Pan Assembly, 22 x 14, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120483	Transit Seat Pan Assembly, 22 x 15, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120484	Transit Seat Pan Assembly, 22 x 16, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120485	Transit Seat Pan Assembly, 22 x 17, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120486	Transit Seat Pan Assembly, 22 x 18, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120487	Transit Seat Pan Assembly, 22 x 19, Spark/Catalyst	1
1 - 5, 6c, 7 - 10	120488	Transit Seat Pan Assembly, 22 x 20, Spark/Catalyst	1
2 - 3, 4a, 5, 7 - 10	124403	Spark / Catalyst with Transit Aluminum Seat Pan Default Setup Hardware, 10"-14" Wide <i>Used with 10-14" wide chairs and without flip back armrests</i>	1
2 - 3, 4a, 4b, 5, 7 - 10	124404	Spark / Catalyst with Transit Aluminum Seat Pan Default Setup Hardware, 15"-18" Wide <i>Used with 15-18" wide chairs and without flip back armrests</i>	1
2 - 3, 4b, 5, 7 - 10	124405	Spark / Catalyst with Transit Aluminum Seat Pan Default Setup Hardware, 19"-22" Wide <i>Used with 19-22" wide chairs and without flip back armrests</i>	1
2 - 3, 4a, 5, 7 - 10	124406	Spark / Catalyst with Transit Aluminum Seat Pan Alternate Setup Hardware, 10"-14" Wide <i>Used with 10-14" wide chairs and with flip back armrests</i>	1
2 - 3, 4a, 4b, 5, 7 - 10	124407	Spark / Catalyst with Transit Aluminum Seat Pan Alternate Setup Hardware, 15"-18" Wide <i>Used with 15-18" wide chairs and with flip back armrests</i>	1
2 - 3, 4b, 5, 7 - 10	124408	Spark / Catalyst with Transit Aluminum Seat Pan Alternate Setup Hardware, 19"-22" Wide <i>Used with 19-22" wide chairs and with flip back armrests</i>	1
3	111408	Pin with Lanyard	4
4a	101708	M6 x 14 BTN SHCS BLKZC w/Patch	2
4b	102634	M6 X 1.0 X 16 BHCS BLK ZC	2
5	100658	M6 Nylock Nut Blk Zn	2
6a	004458	AB Tape, 3" x 4"	2
6b	004459	AB Tape, 3" x 5"	4
6c	004460	AB Tape, 3" x 6"	4
7	100746	M6 Flat Washer Blk Zn	2
8	000990	Seat Pan Clip CAT	2
9	102856	Seat Pan Stop Collar	2
10	102803	Washer Flat .250 X .500 X .027 Hx Blz	6