


# 1" Padded Auto Buckle Pelvic Positioning Belt



| Item Number | Part Number | Description                                                           | Quantity |
|-------------|-------------|-----------------------------------------------------------------------|----------|
| 1, 2a, 3, 4 | 111215      | Belt Pelvic Positioning 1" Auto Buckle Padded 53", Catalyst/Spark     | 1        |
| 1, 2b, 3, 4 | 111216      | Belt Pelvic Positioning 1" Auto Buckle Padded 53", Focus              | 1        |
| 1, 2c, 3, 4 | 111217      | Belt Pelvic Positioning 1" Auto Buckle Padded 53", RG ALX/Rogue/LW XP | 1        |
| 1, 2d, 3, 4 | 111218      | Belt Pelvic Positioning 1" Auto Buckle Padded 53", Clik/Ethos         | 1        |
| 1, 2e, 3, 4 | 111219      | Belt Pelvic Positioning 1" Auto Buckle Padded 53", Flip/Liberty       | 1        |
| 1           | 100613      | Belt Pelvic Positioning 1" Auto Buckle Padded 53"                     | 1        |
| 2a          | 100681      | M6 x 40 CLS8.8 DIN931 HHCS ZC<br><i>Used on Catalyst and Spark</i>    | 2        |
| 2b          | 101125      | M6 x 55 BTN SHCS BLKZC<br><i>Used on Focus</i>                        | 2        |

| Item Number | Part Number | Description                                                                        | Quantity |
|-------------|-------------|------------------------------------------------------------------------------------|----------|
| 2c          | 101781      | M6 x 50mm Partial Thread BHCS Blk ZC<br><i>Used on Tsunami, Rogue and Rogue XP</i> | 2        |
| 2d          | 100670      | M6 x 45 BTN SHCS BLK ZC<br><i>Used on Clik and Ethos</i>                           | 2        |
| 2e          | 101811      | M6 x 35 BHCS Blk Zn<br><i>Used on Flip and Liberty</i>                             | 2        |
| 3           | 100746      | M6 Flat Washer Blk Zn                                                              | 2        |
| 4           | 100658      | M6 Nylock Nut Blk Zn                                                               | 2        |