miniwalk



- Moves up and down during gait, allowing for a more normal heel strike, toe off and mid phase as the child moves forward
- Spring assist system to aide the child to "swing through" during gait
- Standard with castor swivel locks. Can be engaged for in line mobility
- Adjustable angle handlebar encourages more upright or forward trunk orientation
- AN EARLY START ... ages 1 1/2 to 7 years (refer to sizing table below)

Dynamic Gait Movement

The Miniwalk moves up and down with each step, allowing the user to come down into a heel strike and toe off. As the user steps forward, the spring lifts and assists with the swing through phase of gait

Wide range of accessories.....

- Anti tippers
 (req'd for no back and drag brakes)
- No back brakes
- Drag Brakes
- Midline Leg Separation Panel
- Lateral Hip Pads
- Castor Swivel Locks (standard)



Sizing	Mini 1	Mini 2	
Axilla Height Seat Height Max users weight	20-29 " 10.5-19" 90 lb	27-35" 14-25" 90 lb	

Finally... a Gait Trainer with Dynamic Movement to promote normal gait motion!

pacific rehab inc

Web www.pacificrehabinc.com Phone 888-222-9040 480-595-6591 Fax 480-575-7907 Email pacificrehabinc@aol.com Address: 31206 N 65th St Cave Creek, AZ 85331