


Discontinued Clik Caster Mount



This design last issued June 1st, 2016

Item Number	Part Number	Description	Quantity
2	002079	Caster Mount Sleeve LW	1
3	102667	M5 X 30 LHCS BK ZC FULL THREAD	1
6a	002561_ _ _	Caster Arm Short Right Klik <i>Used with 3" and 4" casters</i>	1
6b	002562_ _ _	Caster Arm Short Left Klik <i>Used with 3" and 4" casters</i>	1
6c	002069_ _ _	Caster Arm Medium Right Klik <i>Used with 5" casters</i>	1
6d	002070_ _ _	Caster Arm Medium Left Klik <i>Used with 5" casters</i>	1
6e	002071_ _ _	Caster Arm Long Right Klik <i>Used with 6" casters</i>	1
6f	002072_ _ _	Caster Arm Long Left Klik <i>Used with 6" casters</i>	1
7	102616	M5 X 25 BHCS BLK ZC	1
8	002462	M6 X 1.0 Dual Hole Nut	1
9	100746	M6 Flat Washer Blk Zn	1
10	102627	M6 X 30 SHCS BLK ZC	1
11	101180	M6 x 20 SHCS BLK ZC	1
12	101472	M6 x 45 SHCS BLK ZC	1