


Catalyst 5VX Tipping Lever (Step Tube)



New rear frames required with serial numbers prior to C5053765, C5HD003546, or CVX0091137 Not available with Heavy Duty Option

Item Number	Part Number	Description	Quantity
1 - 6	505420	Catalyst 5VX Tipping Lever (Step Tube) Assembly Kit	1
1, 2	003535	Tipping Lever (Step Tube) Assembly, CAT5VX	1
3	100860	M6 X 35 Hhcs P/Thd Zc	2
4	100658	M6 Nylock Nut Blk Zn	2
5	100746	M6 Flat Washer Blk Zn	2
6	003501	Tipping Lever (Step Tube)/Transit Wheel Sleeve, CAT5VX	1