

## Ethos Tires - Schwalbe Downtown HS 342



*Temporary replacement for Euro Trek Knobby 1 3/8" Tires. Please note tab mounted handrims must be installed prior to the tire being installed.*

Item Number	Part Number	Description	Quantity
1-3	111531	24" x 1 3/8" Schwalbe Downtown HS 342 Tire, Inner Tube and Rim Strip	1
1	111356	24" x 1 3/8" Schwalbe Downtown HS 342 (37-540) Tire, Grey/Black <i>Not available with airless insert</i>	1