

Model	Jazzy Passport ^{1,7}
Part Number	FJPT1700
Weight Capacity	250 lbs.
Maximum Speed ²	Up to 3.6 mph
Ground Clearance ⁴	3" at motor
Turning Radius ⁴	24.25"
Overall Length ⁴	40"
Overall Width ⁴	23.5"
Seat-to-Ground Height ⁴	19"
Seat-to-Foot Platform Height ⁴	14.5" (no adjustment)
Casters ⁴	8" solid
Drive Tires ⁴	11.25" perforated
Suspension	Limited
Range Per Charge ^{2,3}	Up to 7.9 miles (18 AH @ 250 lbs.) Up to 8.6 miles (18 AH @ 220 lbs.) Up to 4.6 miles (10.4 AH travel battery @ 250 lbs.) Up to 5 miles (10.4 AH travel battery @ 220 lbs.)
Total Weight Without Batteries ^{4,6}	59.8 lbs. (unit does not disassemble)
Heaviest Piece When Disassembled ⁴	N/A (unit does not disassemble)
Standard Seating ⁴	Material: Black, Foam Dimensions: 17.5" x 16.5"
Drivetrain	Two-motor, rear-wheel drive
Dual Braking System	Regenerative and electromechanical
Standard Electronics	50 amp, PG VSI
Battery Requirements ⁴	(1) 25.2V 18 AH Lithium ion (1) 25.2V 10.4 AH Lithium ion
Battery Weight ⁵	7.6 lbs. (18 AH) 5.2 lbs. 10.4 AH Lithium ion (air safe travel battery)
Battery Charger	2-amp, off-board
Waranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- ³ Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁴ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- ⁵ Battery weight may vary based on battery manufacturer.
- ⁶ Includes standard seat weight.
- ⁷ Power mobility device, not coded by dme pdac or does not meet criteria.

Jazzy[®] PASSPORT



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power chair only. Overall measurements will vary based on seating and accessory selections.



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® Passport and its lithium ion battery are permitted.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

JAZZYPASSPORT_031721

Jazzy[®] PASSPORT

Simply Fold and Travel

Airline
Approved
Travel Battery
Available



FDA Class II
Medical Device¹

AMERICA'S FAVORITE POWER CHAIR

gojazzy.com |    

MAKES TRAVELING A BREEZE



Folds in just a few simple steps
for effortless transport

Mesh cup holder



Under seat storage

Easy-access
free-wheel lever

FDA Class II
Medical Device¹

Experience the beauty of lightweight travel with the **Jazzy® Passport!** This compact power chair folds in just a few simple steps for effortless transport. With convenient storage the Jazzy Passport makes traveling a breeze.

- A compact and lightweight design
- Back seat pocket for extra storage while on the go
- Includes convenient, under seat mesh storage bag, mesh cup holder and 40" lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort
- Choice of 18 AH battery or 10.4 AH* travel battery

Colors



Black

Accessories

- XLR USB charger

Optional

- Stealth Simplicity G Cushion, 18"x18"



Military Patches



Folded Size
31" x 16" x 23.5"

*Airline safe.