

| | |
|--|---|
| Mark For: Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: () _____ Fax: () _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email | Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ <p style="text-align:center">ADDITIONAL SHIPPING INFORMATION</p> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: () _____ Fax: () _____ |
|--|---|

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|-----------------|----------|-------|------------|-------|---|----------|-------|------------|-------|
| 1. FRAME | | | | | | | | | |

| Base model | | | | |
|--|----------|-------|----------------------------------|------------|
| <i>E1236 coding still applies for EIF7 with sizes 14" or less. 300 lbs weight capacity</i> | | | | |
| <input type="checkbox"/> | 170WM112 | K0005 | Quickie QS5 X | \$2,550.00 |
| <input type="checkbox"/> | 170WM113 | K0005 | Quickie QS5 X - Transit included | \$2,924.00 |

| Frame Style | | | | |
|--|---------|--|------------------|--------------------|
| <i>See figure 1. Measure from front frame to floor</i> | | | | |
| <input type="checkbox"/> | 170FS31 | | Swing Away Frame | NC |
| | 170FS43 | | Fixed Front 75° | COMING SOON |
| | 170FS44 | | Fixed Front 85° | COMING SOON |

| Frame Inset (Fixed Front Only) | | | | |
|--|--------|--|----|--------------------|
| <i>See FIGURE 2 Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2. N/A for 14" and 15" Wide using composite angle adj footrest. 75° only available with 0" inset.</i> | | | | |
| | 170FF1 | | 0" | COMING SOON |
| | 170FF2 | | 1" | COMING SOON |

| Transit Option | | | | |
|-------------------------------------|---------|-------|----------------|----------|
| <i>ISO 7176-19 Transit Approved</i> | | | | |
| <input type="checkbox"/> | 170TO49 | K0108 | Transit System | \$374.00 |

| Frame Options | | | | |
|---|---------|--|--|----------|
| <i>Ergo Touch Points is a rubber cover on the knee bend area of the frame tubes</i> | | | | |
| <input type="checkbox"/> | 1709A11 | | Ergo Touch Points Swing Away hangers or Fixed front | \$150.00 |
| <input type="checkbox"/> | 1709A12 | | Frame Fold-N-Lock mechanism | NC |
| <input type="checkbox"/> | 1709A06 | | Twin LED lights | \$248.00 |



QS5 X - Swing Away



QS5 X - Fixed Front

Power Assist

To order an Empulse R90 Power Assist please see the R90 Order Form for Part #'s

| Option | HCPCS | DESCRIPTION | RETAIL |
|--------------------------|---------|--|------------|
| <input type="checkbox"/> | 209MW91 | Empulse R90 Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 2.5Ahr Removable Battery | \$7,300.00 |
| <input type="checkbox"/> | 209MW92 | Empulse R90 Range Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 5.0Ahr Removable Battery | \$7,400.00 |



[Learn more about the R90](#)

* For Empulse R90 wheelchair compatibility please review the R90 Compatatibility guide.

2. SIZE / SEATING

Frame Width

See **FIGURE 2** Measured from outside of seat tube to the outside of opposite seat tube.

| | | | | | |
|--------------------------|--------|-------|-----|--|----------|
| <input type="checkbox"/> | 170W4 | | 13" | | NC |
| <input type="checkbox"/> | 170W5 | | 14" | | NC |
| <input type="checkbox"/> | 170W6 | | 15" | | NC |
| <input type="checkbox"/> | 170W7 | | 16" | | NC |
| <input type="checkbox"/> | 170W8 | | 17" | | NC |
| <input type="checkbox"/> | 170W9 | | 18" | | NC |
| <input type="checkbox"/> | 170W10 | | 19" | | NC |
| <input type="checkbox"/> | 170W11 | E2201 | 20" | | \$583.00 |

Seat Sling

Standard seat sling includes utility bag sewn at the bottom of upholstery.

| | | | | | |
|--------------------------|--------|--|---------------------|--|----------|
| <input type="checkbox"/> | 170SL1 | | Standard Seat Sling | | NC |
| <input type="checkbox"/> | 170SL6 | | Omit Seat Sling | | -\$15.00 |

Seat Depth

See **Seat Squeeze Matrix** in Appendix section. Measured from back canes to front edge of upholstery. Frame Fold-N-Lock mechanism NA with 13"-14".


| | | | | | |
|--------------------------|---------|-------|-----|--|----------|
| <input type="checkbox"/> | 170SD6 | | 13" | | NC |
| <input type="checkbox"/> | 170SD7 | | 14" | | NC |
| <input type="checkbox"/> | 170SD8 | | 15" | | NC |
| <input type="checkbox"/> | 170SD9 | | 16" | | NC |
| <input type="checkbox"/> | 170SD10 | | 17" | | NC |
| <input type="checkbox"/> | 170SD11 | | 18" | | NC |
| <input type="checkbox"/> | 170SD12 | | 19" | | NC |
| <input type="checkbox"/> | 170SD13 | E2203 | 20" | | \$641.00 |

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

| | | | | | |
|--------------------------|---------|------------------------|--|--------|------------|
| <input type="checkbox"/> | 170JC20 | JAY Bundled Cushion | | | Price List |
| | | Cushion Name or Part # | | | |
| | | Width: | | Depth: | |



3. CASTERS / SEAT HEIGHT

Front Seat Height

| | | | | | | | |
|--------------------------|---------|-------|----|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 170SH41 | 14.0" | NC | <input type="checkbox"/> | 170SH18 | 18.0" | NC |
| <input type="checkbox"/> | 170SH42 | 14.5" | NC | <input type="checkbox"/> | 170SH20 | 18.5" | NC |
| <input type="checkbox"/> | 170SH12 | 15.0" | NC | <input type="checkbox"/> | 170SH21 | 19.0" | NC |
| <input type="checkbox"/> | 170SH13 | 15.5" | NC | <input type="checkbox"/> | 170SH23 | 19.5" | NC |
| <input type="checkbox"/> | 170SH14 | 16.0" | NC | <input type="checkbox"/> | 170SH24 | 20.0" | NC |
| <input type="checkbox"/> | 170SH15 | 16.5" | NC | <input type="checkbox"/> | 170SH26 | 20.5" | NC |
| <input type="checkbox"/> | 170SH16 | 17.0" | NC | <input type="checkbox"/> | 170SH27 | 21.0" | NC |
| <input type="checkbox"/> | 170SH17 | 17.5" | NC | <input type="checkbox"/> | 170SH29 | 21.5" | NC |
| | | | | <input type="checkbox"/> | 170SH36 | 22.0" | NC |

Caster Wheels (option A)

| | | | | | |
|--------------------------|---------|----------------|-----------------------------|--|----------|
| <input type="checkbox"/> | 170CW3 | | 4" Low Profile Polyurethane | | NC |
| <input type="checkbox"/> | 170CW44 | E2219 | 4" x 1.25" Semi Pneumatic | | \$63.00 |
| <input type="checkbox"/> | 170CW5 | | 5" Low Profile Polyurethane | | NC |
| <input type="checkbox"/> | 170CW7 | E2219 | 5" x 1.5" Semi-Pneumatic | | \$100.00 |
| <input type="checkbox"/> | 170CW34 | | 6" Polyurethane | | NC |
| <input type="checkbox"/> | 170CW33 | E2214 E2215 | 6" Pneumatic | | \$88.00 |
| <input type="checkbox"/> | 170CW35 | E2219 | 6" x 1.5" Semi Pneumatic | | \$100.00 |
| <input type="checkbox"/> | 170CW16 | | 7" Polyurethane | | NC |
| <input type="checkbox"/> | 170CW32 | | 8" Polyurethane | | NC |


Caster Fork

| | | | | | |
|--------------------------|---------|-------|---------------------------------------|--|----------|
| <input type="checkbox"/> | 170CF5 | | Multiple Position | | NC |
| <input type="checkbox"/> | 170CF69 | | Multiple Position Single Sided Fork | | \$197.00 |
| <input type="checkbox"/> | 170CF53 | E1015 | Frog Legs 0 weight to 100 lbs. | | \$507.00 |
| <input type="checkbox"/> | 170CF54 | E1015 | Frog Legs 101 lbs. weight to 190 lbs. | | \$507.00 |
| <input type="checkbox"/> | 170CF55 | E1015 | Frog Legs 191 lbs. weight to 265 lbs. | | \$507.00 |

Caster Wheels colors (option B)

Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

| | | | | | |
|--------------------------|---------|--------------------------|------------------------------|--|----------------------|
| <input type="checkbox"/> | 170CW64 | | 4" x 1.5" Aluminum Soft Roll | | \$216.00 |
| <input type="checkbox"/> | 170CW66 | | 5" x 1.5" Aluminum Soft Roll | | \$227.00 |
| <input type="checkbox"/> | 170CW68 | | 6" x 1.5" Aluminum Soft Roll | | \$239.00 |
| | Colors: | <input type="checkbox"/> | Match anodize color package | | Package price |
| | | <input type="checkbox"/> | Match frame color | | \$57.00 |
| | | <input type="checkbox"/> | Black anodize caster only | | \$58.00 |



Seat Squeeze Matrix

| Seat Depth | Max Seat Squeeze |
|------------|------------------|
| 13"-15" | 2.0" |
| 16"-19" | 2.5" |
| 20" | 3.0" |

4. REAR WHEELS / SEAT HEIGHT

| Rear Seat Height | | | | | | | |
|--------------------------|---------|-------|----|--------------------------|---------|-------|----|
| <input type="checkbox"/> | 170RH13 | 14.0" | NC | <input type="checkbox"/> | 170RH19 | 17.0" | NC |
| <input type="checkbox"/> | 170RH14 | 14.5" | NC | <input type="checkbox"/> | 170RH26 | 17.5" | NC |
| <input type="checkbox"/> | 170RH15 | 15.0" | NC | <input type="checkbox"/> | 170RH20 | 18.0" | NC |
| <input type="checkbox"/> | 170RH16 | 15.5" | NC | <input type="checkbox"/> | 170RH27 | 18.5" | NC |
| <input type="checkbox"/> | 170RH17 | 16.0" | NC | <input type="checkbox"/> | 170RH21 | 19.0" | NC |
| <input type="checkbox"/> | 170RH18 | 16.5" | NC | <input type="checkbox"/> | 170RH37 | 19.5" | NC |
| | | | | <input type="checkbox"/> | 170RH22 | 20.0" | NC |

| Axle | | | |
|--------------------------|---------|---------------------|----------|
| <input type="checkbox"/> | 170AX1 | Quick Release Steel | NC |
| <input type="checkbox"/> | 170AX11 | Quad Quick Release | \$144.00 |

Rear Wheel Type

Max weight capacity on Proton wheels is 220 lbs.

| | | | |
|--------------------------|----------|--------------------------------------|------------|
| <input type="checkbox"/> | 170RW103 | 20" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 170RW104 | 22" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 170RW102 | 24" Mag Wheel 5 Spoke | NC |
| <input type="checkbox"/> | 170RW90 | 20" Lite Spoke Wheel | \$234.00 |
| <input type="checkbox"/> | 170RW91 | 22" Lite Spoke Wheel | \$234.00 |
| <input type="checkbox"/> | 170RW92 | 24" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 170RW93 | 25" Lite Spoke Wheel | NC |
| <input type="checkbox"/> | 170RW94 | 26" Lite Spoke Wheel | \$234.00 |
| <input type="checkbox"/> | 170RW95 | 24" Spinergy SPOX Black (silver hub) | \$985.00 |
| <input type="checkbox"/> | 170RW96 | 25" Spinergy SPOX Black (silver hub) | \$985.00 |
| <input type="checkbox"/> | 170RW98 | 24" Spinergy LX Black (silver hub) | \$1,237.00 |
| <input type="checkbox"/> | 170RW99 | 25" Spinergy LX Black (silver hub) | \$1,237.00 |
| <input type="checkbox"/> | 170RW148 | 24" Proton Wheel | \$234.00 |
| <input type="checkbox"/> | 170RW149 | 25" Proton Wheel | \$234.00 |
| <input type="checkbox"/> | 170RW25 | Omit Rear Wheel and Axle * | -\$55.00 |

Spinergy Spoke Colors (Optional)

Please select one of the following Spinergy wheel options first. Includes black hub.

- Spinergy Spox spokes color:
- Black Blue Red White Yellow
- Spinergy LX spokes color:
- Black Blue Green Orange Pink
- Purple Red White Yellow

Rear Wheel Tire

| | | | | |
|--------------------------|---------|----------------|---|----------|
| <input type="checkbox"/> | 170RT1 | | Pneumatic | NC |
| <input type="checkbox"/> | 170RT24 | E2211 E2213 | Pneumatic Airless Insert - Infinergy® Lightweight - E-TPU material | \$114.00 |
| <input type="checkbox"/> | 170RT3 | | Full Poly | NC |
| <input type="checkbox"/> | 170RT38 | E2218 | SHOX Performance Poly | \$110.00 |
| <input type="checkbox"/> | 170RT25 | E2218 | Primo Orion Poly | \$126.00 |
| <input type="checkbox"/> | 170RT39 | E2211 E2212 | High-Pressure Clincher Black | \$110.00 |
| <input type="checkbox"/> | 170RT21 | E2211 E2212 | Schwalbe Marathon Plus Evolution | \$164.00 |

5. WHEEL LOCKS

Wheel Lock

No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth

| | | | |
|--------------------------|---------|-----------------------------|----------|
| <input type="checkbox"/> | 170WL1 | High Mount Push | NC |
| <input type="checkbox"/> | 170WL2 | High Mount Pull | NC |
| <input type="checkbox"/> | 170WL45 | Ergo Scissor | NC |
| <input type="checkbox"/> | 170WL88 | Fold Down Push (Knee lever) | \$86.00 |
| <input type="checkbox"/> | 170WL30 | Compact | \$164.00 |
| <input type="checkbox"/> | 170WL10 | Omit | NC |

| Wheel Dia. | RSH Range | |
|------------|-----------------------------|-----------------------------|
| | Min Squeeze (0° seat angle) | Max Squeeze (9° seat angle) |
| 20" | 14.0"-18.0" | 14.0"-17.5" |
| 22" | 14.0"-19.0" | 14.0"-18.5" |
| 24" | 15.0"-19.5" | 14.5"-19.0" |
| 25" | 15.5"-20.0" | 15.0"-19.5" |
| 26" | 16.0"-20.0" | 15.5"-20.0" |

Handrim

Natural Fits available in 22", 24", and 25". Natural Fit w/ 5 spoke MAG only available in 24". Not all handrims available on all wheels. See page 6 for matrices.

| | | | |
|--------------------------|---------|---|----------|
| <input type="checkbox"/> | 170HR1 | Aluminum Anodized | NC |
| <input type="checkbox"/> | 170HR29 | Aluminum Powder coated colored | \$82.00 |
| <input type="checkbox"/> | 170HR4 | Plastic Coated | \$129.00 |
| <input type="checkbox"/> | 170HR22 | Natural Fit - Standard Grip | \$549.00 |
| <input type="checkbox"/> | 170HR23 | Natural Fit - Super Grip | \$611.00 |
| <input type="checkbox"/> | 170HR25 | Natural Fit - No Thumb Grip | \$475.00 |
| <input type="checkbox"/> | 170HR26 | The Surge - Oval w/Gription Strip | \$493.00 |
| <input type="checkbox"/> | 170HR27 | The Surge LT - Oval w/Gription Strip | \$530.00 |
| <input type="checkbox"/> | 170HR28 | Q-Grip Neoprene Coated | \$408.00 |
| <input type="checkbox"/> | 170HR7 | Omit Handrim <i>Option may not be available on all wheels</i> | NC |

Camber

| | | | |
|--------------------------|---------|--------------|----|
| <input type="checkbox"/> | 170CI1 | 0° | NC |
| <input type="checkbox"/> | 170CI13 | 2° (default) | NC |
| <input type="checkbox"/> | 170CI3 | 4° | NC |

Center of Gravity Preset

| | | | |
|--------------------------|---------|----------------|----|
| <input type="checkbox"/> | 170CC20 | 1.0" (default) | NC |
| <input type="checkbox"/> | 170CC15 | 1.5" | NC |
| <input type="checkbox"/> | 170CC17 | 2.0" | NC |
| <input type="checkbox"/> | 170CC18 | 2.5" | NC |
| <input type="checkbox"/> | 170CC19 | 3.0" | NC |

Wheel Lock Option

| | | | |
|--------------------------|---------|--------------|----|
| <input type="checkbox"/> | 170WL12 | Do Not Mount | NC |
|--------------------------|---------|--------------|----|

Extension Handle

N/A w/ ergo scissor or compact wheel locks

| | | | | |
|--------------------------|--------|-------|-------------------------|---------|
| <input type="checkbox"/> | 170EH7 | E0961 | Extension Handle (Pair) | \$88.00 |
|--------------------------|--------|-------|-------------------------|---------|

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---------------------------------|----------|-------|---------------------------------|----------|--|----------|-------|------------|-------|
| 6. HANGER & FOOTREST | | | | | | | | | |
| Swing Away Hangers | | | | | Footrest Extension Tubes | | | | |
| <i>Swing Away frame only</i> | | | | | <i>See Extension matrix for leg length range</i> | | | | |
| <input type="checkbox"/> | 170H145 | | 60° Swing In/Out Hanger | NC | <input type="checkbox"/> | 170E1 | | Short | NC |
| <input type="checkbox"/> | 170H139 | | 70° Swing In/Out Hanger | NC | <input type="checkbox"/> | 170E2 | | Medium | NC |
| <input type="checkbox"/> | 170H141 | | 80° Swing In/Out Hanger | NC | <input type="checkbox"/> | 170E5 | | Long | NC |
| <input type="checkbox"/> | 170H160 | E0990 | Elevating Hanger - Ext Mount. | \$367.00 | | | | | |
| <input type="checkbox"/> | 170H159 | K0053 | Articulating Hanger - Ext Mount | \$431.00 | | | | | |
| <input type="checkbox"/> | 170H32 | | Omit Hanger and Footplate | -\$15.00 | | | | | |

| Footplates | | | | |
|--|---------|----------------|---|--------------------|
| <i>Heel loop or leg strap not included, See Footplate Options Section. Platform footplates flip up (sideways).</i> | | | | |
| <input type="checkbox"/> | 170F1 | | Composite | STD |
| <input type="checkbox"/> | 170F271 | K0040 E0951 | Composite Angle Adj-(16"-22"w) | \$211.00 |
| <input type="checkbox"/> | 170F5 | | Aluminum (ELR & ALR) | NC |
| <input type="checkbox"/> | 170F3 | K0040 | Alum Adult Angle Adjustable - Ext Mount | \$221.00 |
| <input type="checkbox"/> | 170F259 | K0040 | Multi Angle Adustable <input type="checkbox"/> Black NC <input type="checkbox"/> Std Frame color\$51 | \$221.00 |
| <input type="checkbox"/> | 170F50 | K0040 | Platform Footrest <input type="checkbox"/> Black NC <input type="checkbox"/> Std Frame color\$51 | COMING SOON |
| <input type="checkbox"/> | 107F335 | K0040 | Platform Ultra Lite Footrest | COMING SOON |
| <input type="checkbox"/> | 107F336 | | Platform Carbon Fiber footrest | COMING SOON |
| <input type="checkbox"/> | 170F7 | | Omit Footplate and Ext. Tube | -\$15.00 |

| Footrest Options | | | | |
|--------------------------|---------|-------|---|--------------------|
| <input type="checkbox"/> | 170HL1 | E0951 | Heel Loop <i>Swing away only</i> | \$47.00 |
| | 170HL5 | K0038 | Leg Strap | COMING SOON |
| | 107DI04 | | Extra Leg Strap | |
| <input type="checkbox"/> | 170IG2 | | Impact Guard Neoprene | \$105.00 |
| | 170DI06 | | Platform Position Plates | |
| | 170F35 | | Auto-Folding Platform <i>Carbon fiber only</i> | COMING SOON |

| Articulating / Elevating Hanger Cover | | | | |
|--|---------|--|----------------|----------|
| <input type="checkbox"/> | 170DI47 | | Neoprene Cover | \$397.00 |
| <input type="checkbox"/> | 170DI48 | | Gel Pad Cover | \$467.00 |

| Lower Leg Length (LLL) | | | | |
|---|-----|--------------------------|-----|--------------------------|
| <i>Lower Leg Length is measured from seat sling to the top rear of the footrest. See Appendix on Page 7. A minimum 2" ground clearance is required below the footrest. If the selected configuration results in less than 2" of ground clearance the footrest will be raised to 2" ground clearance, which may alter the lower leg length measurement. Available lower leg length is dependent on frame selection including adapter spacers. Reference LLL matrix in page 7</i> | | | | |
| <input type="checkbox"/> | 10" | <input type="checkbox"/> | 11" | <input type="checkbox"/> |
| <input type="checkbox"/> | 12" | <input type="checkbox"/> | 13" | <input type="checkbox"/> |
| <input type="checkbox"/> | 14" | <input type="checkbox"/> | 15" | <input type="checkbox"/> |
| <input type="checkbox"/> | 16" | <input type="checkbox"/> | 17" | <input type="checkbox"/> |
| <input type="checkbox"/> | 18" | <input type="checkbox"/> | 19" | |

| 7. BACKREST / SEATING | | | | |
|--|----------|-------|--|----------|
| Back Type | | | Backrest Upholstery | |
| <input type="checkbox"/> | 107BT123 | | Standard Fix (90°) | NC |
| <input type="checkbox"/> | 170BT47 | | Angle Adjustable Back (80° - 105°) | NC |
| Select desired setting <input type="checkbox"/> 80° <input type="checkbox"/> 85° <input type="checkbox"/> 90° <input type="checkbox"/> 100° | | | | |
| <input type="checkbox"/> | 107BT124 | | Fixed Half Fold down (to the front) | \$216.00 |
| <input type="checkbox"/> | 170BT43 | | Lumbar curve bend | \$87.00 |
| Select desired setting <input type="checkbox"/> 5° forward <input type="checkbox"/> 5° backward | | | | |
| <i>EXO standard and Tension adjustable upholstery binding is always black.</i> | | | | |
| <input type="checkbox"/> | 170BU95 | | Standard Upholstery EXO | NC |
| Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00 | | | | |
| <input type="checkbox"/> | 170BU70 | E2611 | Backrest Tension Adjust Ballistic <i>*Reference 186BU70 for billing</i> | \$367.00 |
| <input type="checkbox"/> | 170BU89 | E2611 | Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i> | \$387.00 |
| Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00 | | | | |
| <input type="checkbox"/> | 170BU90 | E2611 | Backrest Tension Adj EXO Pro 3DX <i>*Reference 186BU71 for billing</i> | \$387.00 |
| Color Binding: <input type="checkbox"/> Silver <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue \$26.00 | | | | |
| <input type="checkbox"/> | 170BU14 | | Omit Back Upholstery | -\$15.00 |

| Back Height Setting | | | | |
|---|---------|---------|---|----|
| <i>Half Fold Down N/A with the 10" - 13" back heights. Lumbar curve backrest with integral push handles available with 12" thru 20" back heights. Lumbar curve backrest with Fold Down or Omit push handles available with 14" thru 20" back heights.</i> | | | | |
| <input type="checkbox"/> | 170B161 | Setting | 10" - 12" Back Height <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" | NC |
| <input type="checkbox"/> | 170B162 | Setting | 12" - 14" Back Height <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" | NC |
| <input type="checkbox"/> | 170B163 | Setting | 14" - 16" Back Height <input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" | NC |
| <input type="checkbox"/> | 170B164 | Setting | 16" - 18" Back Height <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" | NC |
| <input type="checkbox"/> | 170B165 | Setting | 18" - 20" Back Height <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" | NC |

| Military Service Patches (Optional) | | | | |
|--|---------|--|-------------------------------------|----------|
| <input type="checkbox"/> | 170VU2 | | Army Patch | \$116.00 |
| <input type="checkbox"/> | 170VU4 | | Air Force Patch | \$116.00 |
| <input type="checkbox"/> | 170VU12 | | Coast Guard Patch | \$116.00 |
| <input type="checkbox"/> | 170VU1 | | Marines Patch | \$116.00 |
| <input type="checkbox"/> | 170VU3 | | Navy Patch | \$116.00 |
| <input type="checkbox"/> | 170VU5 | | Purple Heart Patch | \$116.00 |
| <input type="checkbox"/> | 170VU9 | | Purple Heart with Army Patch | \$181.00 |
| <input type="checkbox"/> | 170VU11 | | Purple Heart with Air Force Patch | \$181.00 |
| <input type="checkbox"/> | 170VU13 | | Purple Heart with Coast Guard Patch | \$181.00 |
| <input type="checkbox"/> | 170VU8 | | Purple Heart with Marines Patch | \$181.00 |
| <input type="checkbox"/> | 170VU10 | | Purple Heart with Navy Patch | \$181.00 |

| Push Handles | | | | |
|--|---------|--|---|----------|
| <i>Std EXO Uph N/A with Fold down. Height Adjustable push handles available with EXO / EXO Pro Tension A</i> | | | | |
| <input type="checkbox"/> | 170PH6 | | Integral Push Handles | NC |
| <input type="checkbox"/> | 170PH5 | | Fold Down Push Handles <i>Available w/12"-20" back height</i> | \$116.00 |
| <input type="checkbox"/> | 170PH21 | | Height adjustable push Handles <i>Available w/ 12"-20" back height. N/A w/Half fold down back</i> | \$175.00 |
| <input type="checkbox"/> | 170PH25 | | Height adjustable push Handles for JAY Backs <i>N/A w/Half fold down back</i> | \$175.00 |
| <input type="checkbox"/> | 170PH22 | | Omit Push Handles - <i>Straight Post</i> | NC |

| JAY Back (Optional) | | | | |
|---|---------|--|-----------------------------|------------|
| <i>To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s</i> | | | | |
| <input type="checkbox"/> | 170JB01 | | J3 Back (See J3 Order Form) | Price List |

| ✓ | Option # | HCPCS | Definition | Price | ✓ | Option # | HCPCS | Definition | Price |
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|
|---|----------|-------|------------|-------|---|----------|-------|------------|-------|

8. ARMREST / SIDE GUARDS

Padded Swing Away

Some limitations when used w/ angle adj back in the rear most setting. 10"-12" height ranges N/A w/Swing away anti tip. Please select a setting, will default to 11" if no selection is made.

| | | | | |
|--------------------------|----------|---------|--|----------|
| <input type="checkbox"/> | 170AA162 | | Padded Swing-Away Fixed height (13") | NC |
| <input type="checkbox"/> | 170AA18 | E0973 | Padded Swing-Away STD (9" to 14") | \$263.00 |
| | | Setting | <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14" | |

Armrest - Height Adjustable

| | | | | |
|--------------------------|---------|-------|---|----------|
| <input type="checkbox"/> | 170AA11 | E0973 | Single Post Hgt. Adj. Std (8.5"-13.5") | \$333.00 |
| <input type="checkbox"/> | 170AA29 | E0973 | Single Post Hgt. Adj. Low (7"-11") | \$333.00 |
| <input type="checkbox"/> | 170AA65 | E0973 | Flip Back Height Adj <i>LITE</i> (8.5"-13.5") | \$367.00 |

Armpad - For Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

| | | | | |
|--------------------------|---------|--|-----------------------------|----|
| <input type="checkbox"/> | 170AM23 | | Desk Length - Classic (10") | NC |
| <input type="checkbox"/> | 170AM24 | | Full Length - Classic (14") | NC |
| <input type="checkbox"/> | 170AM27 | | Ergo Pad (9.5"-13.5") | NC |

Omit Armrest

| | | | | |
|--------------------------|---------|--|--------------|----|
| <input type="checkbox"/> | 170AA14 | | Omit Armrest | NC |
|--------------------------|---------|--|--------------|----|

9. KEY OPTIONS & ACCESSORIES

Anti-Tip

| | | | | |
|--------------------------|--------|-------|--|----------|
| <input type="checkbox"/> | 170AT1 | E0971 | Anti Tips (Pair) | \$175.00 |
| <input type="checkbox"/> | 170AT | E0971 | Anti Tip swing away | |
| | | | <input type="checkbox"/> 170AT88 Right Mount | \$102.00 |
| | | | <input type="checkbox"/> 170AT89 Left Mount | \$102.00 |
| | | | <input type="checkbox"/> 170AT90 Pair | \$204.00 |

Additional Accessories

| | | | | |
|--------------------------|---------|--|------------------|---------|
| <input type="checkbox"/> | 170A1 | | Backpack - Adult | \$88.00 |
| <input type="checkbox"/> | 170A2 | | Backpack - Kids | \$70.00 |
| <input type="checkbox"/> | 1701A01 | | Seat Pouch Black | \$58.00 |

Additional Cover Accessories

| | | | | |
|--------------------------|---------|--|--------------------------|----------|
| <input type="checkbox"/> | 1708A02 | | Gel Cover Hanger/Frame | \$129.00 |
| <input type="checkbox"/> | 1708A01 | | Gel Cover Swing Away Arm | \$175.00 |

Positioning Belts

| | | | | |
|--------------------------|---------|-------|-------------------------------------|----------|
| <input type="checkbox"/> | 170ZZ15 | E0978 | JAY 2PT Belt 1.5" Push Btn 7.5" Pad | \$101.00 |
| <input type="checkbox"/> | 170ZZ26 | E0978 | JAY 2PT Belt 2" Push Btn 9.5" Pad | \$101.00 |
| <input type="checkbox"/> | 170ZZ14 | E0978 | JAY 2PT Belt 1.5" Push Btn No Pad | \$48.00 |
| <input type="checkbox"/> | 170A118 | E0978 | Belt 2" Aircraft Buckle Pad w/Clamp | \$139.00 |
| <input type="checkbox"/> | 170A120 | | Positioning Belt Slide Clamps | \$58.00 |

Side Guard

Side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist. Aluminum & carbon side guards some limitations with low RSH and large rear wheels

| | | | | |
|--------------------------|---------|-------|-------------------------------|----------|
| <input type="checkbox"/> | 170SG33 | K0108 | Plastic Standard | \$169.00 |
| <input type="checkbox"/> | 170SG29 | K0108 | Plastic Low | \$169.00 |
| <input type="checkbox"/> | 170SG35 | K0108 | Aluminum | \$286.00 |
| <input type="checkbox"/> | 170SG36 | K0108 | Aluminum w/ Fender | \$321.00 |
| <input type="checkbox"/> | 170SG12 | K0108 | Carbon Fiber | \$537.00 |
| <input type="checkbox"/> | 170SG41 | K0108 | Carbon Fiber Long w/LG Fender | \$607.00 |

Side Guard - Color Option

Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber

| | | | | |
|--------------------------|---------|--|------------------------------------|---------|
| <input type="checkbox"/> | 170GE01 | | Black | NC |
| <input type="checkbox"/> | 170GE02 | | Matching Frame- Standard Colors | NC |
| <input type="checkbox"/> | 170GE03 | | Matching Frame- Personalized Color | \$58.00 |

Spoke Guards

Only available w/ Lite Spoke wheels. Not available with 5 spoke MAG wheels

| | | | | |
|--------------------------|---------|-------|----------------------------|----------|
| <input type="checkbox"/> | 1704A01 | K0065 | Quickie Black Spoke Guards | \$134.00 |
| <input type="checkbox"/> | 1704A02 | K0065 | Quickie Clear Spoke Guards | \$134.00 |

Miscellaneous

| | | | | |
|--------------------------|----------|-------|--|----------|
| <input type="checkbox"/> | 1708A03 | | Caddy | \$58.00 |
| <input type="checkbox"/> | 170A2100 | E2207 | Crutch Holder | \$227.00 |
| <input type="checkbox"/> | 1708A21 | | Tool Kit | \$63.00 |
| <input type="checkbox"/> | 170A999 | | Quickie Water Bottle & Holder | \$63.00 |
| <input type="checkbox"/> | 1708A47 | | Mobile phone pocket | \$30.00 |
| <input type="checkbox"/> | 170FO4 | E0971 | Transit wheels (Pair) | \$122.00 |
| <input type="checkbox"/> | 170AT | E0971 | Tip Assist | |
| | | | <input type="checkbox"/> 170AT83 Right Mount | \$36.00 |
| | | | <input type="checkbox"/> 170AT84 Left Mount | \$36.00 |
| | | | <input type="checkbox"/> 170AT91 Pair | \$71.00 |

10. COLOR

Personalized Crossbrace Color

Option has an extended leadtime. Frame colors only, NA with Personalized frame colors. If no color pick crossbrace will be Gloss Black.

| | | | | | | | | | |
|--------------------------|---------|--------------------------|--|--|--|--|--|--|----------|
| <input type="checkbox"/> | 170GS15 | Crossbrace Color package | | | | | | | \$250.00 |
|--------------------------|---------|--------------------------|--|--|--|--|--|--|----------|

Standard Frame Colors

Handrim color option only available if Aluminum Coated colored option was selected.

| | | Frame | X-brace | Handrim |
|---------|----------------------|--------------------------|--------------------------|--------------------------|
| 170C140 | Aztec Gold | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C118 | Bead Blast | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C2 | Black | <input type="checkbox"/> | STD | <input type="checkbox"/> |
| 170C141 | Black Cherry | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C10 | Black Opal | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C38 | Blue Opal | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C36 | Candy Blue | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C14 | Candy Purple | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C35 | Candy Red | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C142 | Electric Blue | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C143 | Evergreen | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C44 | Glossy White | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C126 | Glow | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C144 | Green Apple | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C147 | Hot Sparkle Pink | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C148 | Mauve Pink | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C103 | Rootbeer | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C114 | Sparkle Silver | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C145 | Sunrise Orange | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C113 | Titanium Color Paint | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C8 | Yellow | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Matte Frame Colors

Handrim color option only available if Aluminum Coated colored options was selected.

| | | Frame | X-brace | Handrim |
|---------|---------------------|--------------------------|--------------------------|--------------------------|
| 170C116 | Matte Black | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C151 | Matte Black Cherry | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C150 | Matte Electric Blue | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C149 | Matte Evergreen | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| 170C152 | Matte Purple | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Personalized Frame Color

Only frame will have the personalized color; Crossbrace, handrim, footplate and Aluminum side guard not included (unless selected).

| | | | | | | | |
|--------------------------|---------|--|--|--|--|--|----------|
| <input type="checkbox"/> | 170C112 | Desert Camouflage | | | | | \$460.00 |
| <input type="checkbox"/> | 170C43 | Mossy Oak Camouflage | | | | | \$460.00 |
| <input type="checkbox"/> | 170C42 | Stars & Stripes | | | | | \$460.00 |
| <input type="checkbox"/> | 170C108 | Zebra (Choose above base color) Std frame colors with White unless above base paint is chosen N/A w/ Black, Black Opal & Matte Black | | | | | \$460.00 |

Anodized Color Package

Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog suspension Legs.

| | | | | | | | |
|--------------------------|---------|------------|--|--|--|--|----------|
| <input type="checkbox"/> | 170EB08 | Blue | | | | | \$187.00 |
| <input type="checkbox"/> | 170EB13 | Gold | | | | | \$187.00 |
| <input type="checkbox"/> | 170EB09 | Lime Green | | | | | \$187.00 |
| <input type="checkbox"/> | 170EB10 | Purple | | | | | \$187.00 |
| <input type="checkbox"/> | 170EB11 | Red | | | | | \$187.00 |
| <input type="checkbox"/> | 170EB12 | Ti Grey | | | | | \$187.00 |

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

| | | | | | | | |
|--------------------------|--------|--------|--|--|--|--|----|
| <input type="checkbox"/> | 170D16 | Black | | | | | NC |
| <input type="checkbox"/> | 170D17 | Blue | | | | | NC |
| <input type="checkbox"/> | 170D18 | Orange | | | | | NC |
| <input type="checkbox"/> | 170D19 | Red | | | | | NC |
| <input type="checkbox"/> | 170D20 | Silver | | | | | NC |

Built For Me

*To inquire about any unique personalized for Quikie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011*



11. APPENDIX

Wheel - Handrims Matrix

| Handrim | Rear Wheel | | | |
|----------------------------------|-------------|----------------|----------|--------|
| | 5 Spoke Mag | Lite Spoke | Spinergy | Proton |
| Aluminum Ano / Painted | 20,22,24 | 20,22,24,25,26 | 24,25 | 24,25 |
| Plastic Coat | 20,22,24 | 22,24 | 24,25 | - |
| Natural Fit - Standard Grip | 24 | 22,24,25,26 | 24,25 | - |
| Natural Fit - Super Grip | 24 | 22,24,25,26 | 24,25 | - |
| Natural Fit - No Thumb Grip | 24 | 22,24,25,26 | 24,25 | - |
| Surge - Oval w/Gription Strip | - | 22,24,25,26 | 24,25 | - |
| Surge LT - Oval w/Gription Strip | - | 22,24,25,26 | 24,25 | - |
| Q Grip Neoprene Coated | - | 24,25 | 24,25 | - |
| Omit Handrim | 20,22,24 | - | 24,25 | - |

Wheel - Tire Matrix

| Tire | Rear Wheel | | | |
|-----------------------------|-------------|-------------|----------|--------|
| | 5 Spoke Mag | Lite Spoke | Spinergy | Proton |
| Pneumatic | 20,22,24 | 20,22,24,26 | 24 | 24 |
| Pneumatic w/ Airless Insert | 20,22,24 | 20,22,24,26 | 24 | - |
| Full Poly | 20,22,24 | 20,22,24 | 24 | - |
| SHOX Performance Poly | 20,22,24 | 20,22,24,25 | 24,25 | - |
| Primo Orion Poly | 20,22,24 | 20,22,24 | 24 | - |
| High-Pressure Clincher Blk | 24 | 22,24,26 | 24 | 24 |
| Schwalbe Marathon | 22, 24 | 24,25,26 | 24,25 | 24,25 |

20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm

11. APPENDIX

Note 1: Fork Hole 1+ indicates that a 10mm spacer is used to achieve additional Front Seat Height (FSH).

Note 2: FSH values are rounded estimates and assume flat seat heights (0 degree seat angle). Increasing seat angle can affect final FSH.

| FSH for 4" Casters | | | | | | |
|--------------------|-----------|---------------|------|------|------|------|
| Caster Link | Fork Size | Hole Position | | | | |
| | | 4 | 3 | 2 | 1 | 1+ |
| Hemi | Multi 90 | - | - | 14.0 | 14.5 | 15.0 |
| | Multi 115 | - | 14.5 | 15.0 | 15.5 | 16.0 |
| | Multi 131 | - | 15.0 | 15.5 | 16.0 | 16.5 |
| | Multi 158 | - | 16.0 | 16.5 | 17.0 | 17.5 |
| | SS Low | - | - | 14.0 | 14.5 | 15.0 |
| | SS Std | - | 14.5 | 15.0 | 15.5 | 16.0 |
| | Frog Legs | - | 15.5 | 16.0 | 16.5 | - |
| Short | Multi 90 | - | - | 15.5 | 16.0 | 16.0 |
| | Multi 115 | - | 16.0 | 16.5 | 17.0 | 17.0 |
| | Multi 131 | - | 16.5 | 17.0 | 17.5 | 18.0 |
| | Multi 158 | - | 17.5 | 18.0 | 18.5 | 19.0 |
| | SS Low | - | - | 15.5 | 16.0 | 16.0 |
| | SS Std | - | 16.0 | 16.5 | 17.0 | 17.0 |
| | Frog Legs | - | 17.0 | 17.5 | 18.0 | - |
| Tail | Multi 90 | - | - | 16.5 | 17.0 | 17.5 |
| | Multi 115 | - | 17.0 | 17.5 | 18.0 | 18.5 |
| | Multi 131 | - | 17.5 | 18.0 | 18.5 | 19.0 |
| | Multi 158 | - | 18.5 | 19.0 | 19.5 | 20.0 |
| | SS Low | - | - | 16.5 | 17.0 | 17.5 |
| | SS Std | - | 17.0 | 17.5 | 18.0 | 18.5 |
| | Frog Legs | - | 18.0 | 18.5 | 19.0 | - |

| FSH for 5" Casters | | | | | | |
|--------------------|-----------|---------------|------|------|------|------|
| Caster Link | Fork Size | Hole Position | | | | |
| | | 4 | 3 | 2 | 1 | 1+ |
| Hemi | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | 15.5 | 16.0 | 16.0 |
| | Multi 131 | 15.0 | 15.5 | 16.0 | 16.5 | 17.0 |
| | Multi 158 | 16.0 | 16.5 | 17.0 | 17.5 | 18.0 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | 15.0 | 15.5 | 16.0 | 16.0 |
| | Frog Legs | - | - | 16.5 | 17.0 | - |
| Short | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | 16.5 | 17.0 | 17.5 |
| | Multi 131 | 16.5 | 17.0 | 17.5 | 18.0 | 18.0 |
| | Multi 158 | 17.5 | 18.0 | 18.5 | 19.0 | 19.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | 16.0 | 16.5 | 17.0 | 17.5 |
| | Frog Legs | - | - | 17.5 | 18.0 | - |
| Tail | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | 18.0 | 18.5 | 19.0 |
| | Multi 131 | 17.5 | 18.0 | 18.5 | 19.0 | 19.5 |
| | Multi 158 | 18.5 | 19.0 | 19.5 | 20.0 | 20.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | 17.5 | 18.0 | 18.5 | 19.0 |
| | Frog Legs | - | - | 19.0 | 19.5 | - |

| FSH for 6" Casters | | | | | | |
|--------------------|-----------|---------------|------|------|------|------|
| Caster Link | Fork Size | Hole Position | | | | |
| | | 4 | 3 | 2 | 1 | 1+ |
| Hemi | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | 16.0 | 16.5 |
| | Multi 131 | - | 16.0 | 16.5 | 17.0 | 17.0 |
| | Multi 158 | 16.5 | 17.0 | 17.5 | 18.0 | 18.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | 15.5 | 16.0 | 16.5 |
| | Frog Legs | - | - | - | 17.0 | - |
| Short | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | 17.5 | 18.0 |
| | Multi 131 | - | 17.5 | 18.0 | 18.5 | 18.5 |
| | Multi 158 | 18.0 | 18.5 | 19.0 | 19.5 | 19.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | 17.0 | 17.5 | 18.0 |
| | Frog Legs | - | - | - | 18.5 | - |
| Tail | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | 19.0 | 19.5 |
| | Multi 131 | - | 18.5 | 19.0 | 19.5 | 20.0 |
| | Multi 158 | 19.0 | 19.5 | 20.0 | 20.5 | 21.0 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | 18.5 | 19.0 | 19.5 |
| | Frog Legs | - | - | - | 20.0 | - |

| FSH for 7" Casters | | | | | | |
|--------------------|-----------|---------------|------|------|------|------|
| Caster Link | Fork Size | Hole Position | | | | |
| | | 4 | 3 | 2 | 1 | 1+ |
| Hemi | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 17.0 | 17.5 | 18.0 |
| | Multi 158 | 17.0 | 17.5 | 18.0 | 18.5 | 19.0 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |
| Short | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 18.5 | 19.0 | 19.0 |
| | Multi 158 | 18.5 | 19.0 | 19.5 | 20.0 | 20.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |
| Tail | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 19.5 | 20.0 | 20.5 |
| | Multi 158 | 19.5 | 20.0 | 20.5 | 21.0 | 21.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |

| FSH for 8" Casters | | | | | | |
|--------------------|-----------|---------------|------|------|------|------|
| Caster Link | Fork Size | Hole Position | | | | |
| | | 4 | 3 | 2 | 1 | 1+ |
| Hemi | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 17.5 | 18.0 | 18.5 |
| | Multi 158 | - | 18.0 | 18.5 | 19.0 | 19.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |
| Short | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 19.0 | 19.5 | 19.5 |
| | Multi 158 | - | 19.5 | 20.0 | 20.5 | 20.5 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |
| Tail | Multi 90 | - | - | - | - | - |
| | Multi 115 | - | - | - | - | - |
| | Multi 131 | - | - | 20.0 | 20.5 | 21.0 |
| | Multi 158 | - | 20.5 | 21.0 | 21.5 | 22.0 |
| | SS Low | - | - | - | - | - |
| | SS Std | - | - | - | - | - |
| | Frog Legs | - | - | - | - | - |

| Lower Leg Length Matrix | | | | |
|--|-----------|---|-----|--|
| <i>Note for Swing Away Hangers only: subtract 1" from the max achievable LLL for long extension tubes if the Multi Angle Adjustable footplate is selected.</i> | | | | |
| Hanger | Extension | LLL (inches) (measured from top of seat sling to top of footrest) | | Min Allowable FSH (due to 2" ground clearance requirement) |
| | | Min | Max | |
| 60° Swing In/Out | Short | 10 | 14 | No lower FSH limit |
| | Med | 11 | 16 | |
| | Long | 12 | 18 | |
| 70° Swing In/Out | Short | 10 | 16 | No lower FSH limit |
| | Med | 12 | 17 | |
| | Long | 13 | 19 | |
| 80° Swing In/Out | Short | 11 | 15 | No lower FSH limit |
| | Med | 12 | 17 | |
| | Long | 14 | 19 | |
| ELR | Short | 11 | 16 | No lower FSH limit |
| | Med | 11 | 18 | |
| | Long | 13 | 19 | |
| ALR | Short | 11 | 16 | No lower FSH limit |
| | Med | 13 | 18 | |
| | Long | 15 | 19 | |
| 75° Fixed Front | Short | 13 | 18 | |
| | Med | 13 | 20 | |
| | Long | 14 | 20 | |
| 85° Fixed Front 0" inset | Short | 13 | 18 | |
| | Med | 13 | 20 | |
| | Long | 15 | 20 | |
| 85° Fixed Front 1" inset | Short | 13 | 18 | |
| | Med | 15 | 20 | |
| | Long | 17 | 20 | |

COMING SOON



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122